**en** **Instructions for use - Infusion set**

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
V86  
L86  
S86  
B86

I86  
I99  
SMS Primary Set  
Needle-free SMS Primary Set

**(en) Exemplary product illustration(s) - deviations in product configuration are possible!** **(de) Exemplarische Produktabbildung(en) - Abweichungen in der Produktkonfiguration sind möglich!** **(cs) Příklad(y) zobrazení produktu - konfigurace produktu se může lišit!** **(da) Produktillustrationseksempel(er) - afvigelser i produktkonfiguration kan forekomme!** **(es) Ilustración(es) de ejemplo del producto - es posible que haya diferencias entre las ilustraciones de ejemplo y la configuración del producto!** **(fr) Illustration(s) du dispositif à titre d'exemple - des différences de configuration du dispositif sont possibles!** **(it) Illustrazione/! del prodotto a titolo esemplificativo - è possibile che vi siano differenze nella configurazione del prodotto!** **(nl) Voorbeeldillustratie(s) van product - de configuratie van uw product kan afwijken!** **(no) Produktillustrasjon(er) - forskjellige produktkonfigurasjoner kan forekomme!** **(pl) Przykładowa(-e) ilustracja(-e) produktu - zastrzega się możliwość zmian w konfiguracji produktu!** **(pt) Ilustração(ões) exemplificativa(s) do(s) produto(s) - são possíveis desvios na configuração do produto!** **(sv) Exempelbild(-er) av produkten - avvikelser i produktens konfiguration är möjliga!**

- 1 **en** Closure-piercing device with protective cap **de** Einsteckteil mit Schutzkappe **cs** Uzavírací-propichovací díl s ochrannou krytkou **da** Spike-enhed med beskyttelseshætte **es** Dispositivo de perforación de cierre con caperuza de protección **fr** Dispositif de perforation du bouchon avec capuchon de protection **it** Dispositivo di perforazione della chiusura con cappuccio protettivo **nl** Afsluitprikker met beschermkapje **no** Gjennomhullingsenhet med beskyttelseshette **pl** Łącznik kołcowy z kapturkiem ochronnym **pt** Dispositivo de perfuração com tampa de proteção **sv** Insticksdel med skyddshätta
- 2 **en** Closure clamp **de** Schließklemme **cs** Uzavírací svorka **da** Lukkeklemme **es** Abrazadera de cierre **fr** Clamp de fermeture **it** Morsetto di chiusura **nl** Afsluitklem **no** Lukkeklemme **pl** Zacisk zamykający **pt** Braçadeira de fecho **sv** Avstängningsklämma
- 3 **en** Connection port **de** Anschlussstück **cs** Připojovací port **da** Tilslutningsport **es** Puerto de conexión **fr** Port de connexion **it** Porta di connessione **nl** Aansluitpoort **no** Tilkoblingsport **pl** Port przyłączeniowy **pt** Porta de conexão **sv** Anslutningsport
- 4 **en** Level ring **de** Niveauring **cs** Hladinový kroužek **da** Niveauring **es** Anillo de nivel **fr** Repère circulaire de niveau **it** Anello di livello **nl** Niveauring **no** Nivåmerke **pl** Pierścień poziomy **pt** Anel de nível **sv** Nivåring
- 5 **en** Drip chamber **de** Tropfkammer **cs** Kapací komůrka **da** Dråbekammer **es** Cámara de goteo **fr** Chambre compte-gouttes **it** Camera di gocciolamento **nl** Druppelkamer **no** Dryppkammer **pl** Komora kroplowa **pt** Câmara de gotejamento **sv** Droppkammare
- 6 **en** Roller clamp **de** Rollenklemme **cs** Tlačka **da** Rulleklemme **es** Abrazadera de rodillos **fr** Clamp à roulette **it** Roller clamp **nl** Rolklem **no** Rulleklemme **pl** Zacisk rolkowy **pt** Grampo de rolos **sv** Rullklämma
- 7 **en** Tubing **de** Schlauch **cs** Hadička **da** Slang **es** Tubo flexible **fr** Tubulure **it** Tubazione **nl** Slangen **no** Slange **pl** Przewód **pt** Tubos **sv** Slangar
- 8 **en** Injection site **de** Zuspritzmöglichkeit **cs** Injekční místo **da** Injektionssted **es** Lugar de la inyección **fr** Site de perfusion **it** Sito di iniezione **nl** Injectieplaats **no** Injeksjonssted **pl** Miejsce wkłucia **pt** Local da injeção **sv** Injektionsställe
- 9 **en** Patient connection port with protective cap **de** Patientenanschluss mit Schutzkappe **cs** Připojovací port pacienta s ochrannou krytkou **da** Patienttilslutningsport med beskyttelseshætte **es** Puerto de conexión del paciente con caperuza de protección **fr** Port de connexion patient avec bouchon de protection **it** Porta di connessione del paziente con cappuccio protettivo **nl** Patiëntaansluitpoort met beschermkapje **no** Pasientforbindelsesport med beskyttelseshette **pl** Port przyłączeniowy pacjenta z kapturkiem ochronnym **pt** Porta de conexão do doente com tampa de proteção **sv** Patientanslutningsport med skyddshätta


**en** **Symbols** - Explanation of symbols available at: **de** **Symbole** - Symbolerklärung verfügbar unter: **cs** **Symboly** - Vysvětlení symbolů je k dispozici na stránce: **da** **Symboler** - Symbolforklaring tilgængelig på: **es** **Símbolos** - Explicación de los símbolos disponible en: **fr** **Symboles** - Explication des symboles disponible à l'adresse Internet suivante: **it** **Simboli** - Spiegazione dei simboli disponibile all'indirizzo: **nl** **Symbolen** - Uitleg van de symbolen beschikbaar op: **no** **Symboler** - Forklaring av symboler tilgjengelig på: **pl** **Symboly** - Objaśnienia symboli dostępne na stronie: **pt** **Símbolos** - Explicação dos símbolos disponível em: **sv** **Symboler** - Förklaring av symboler som finns på:

<https://www.codancompanies.com/download-centre/>

en Notes on product use (graphic brief instructions) de Hinweise zur Produkthanwendung (grafische Kurzanleitung) cs Poznámky k použití produktu (grafický stručný návod) da Bemærkninger om produktanvendelse (kort grafisk vejledning) es Notas sobre el uso del producto (breves instrucciones gráficas) fr Remarques relatives à l'utilisation du dispositif (guide rapide graphique) it Note sull'uso del prodotto (brevi istruzioni grafiche) nl Opmerkingen bij gebruik van het product (beknopte grafische instructies) no Merknader om bruk av produktet (grafisk hurtigveiledning) pl Wskazówki odnośnie zastosowania (skrócona instrukcja graficzna) pt Notas sobre a utilização do produto (breves instruções gráficas) sv Information om produktens användning (kortfattade anvisningar i bildform)


### V86/L86

- vented closure-piercing device


### B86/S86

- non-vented closure-piercing device


### I86/I99

- vented closure-piercing device

### I86/I99

- non-vented closure-piercing device

## Infusion set

suitable for use with compatible pumps as well as gravity feed

### Intended purpose

Transfer of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions by intravenous infusion.

### Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions.

### Contraindications

The product is not designed and approved for:

- Home use
- Use by unqualified, untrained personnel

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of:

- Blood and blood components
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions with regard to incompatibilities (interaction material/pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

### Patient target group

Infusion set (20 drops/ml):

Patients of at least 2 years of age and/or a weight of 12 kg.

Infusion set with micro drip orifice (60 drops/ml) and with identifier "-MT" in the product name:

Patients from birth/from birth weight (including premature babies)

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

### User

Qualified medical and pharmaceutical personnel only.

For use in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions.

### Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions carry implicit risks and should therefore only be carried out if they improve the patient's health.

### General safety instructions

• The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g., leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • In case of

consecutive administration of different drugs with the same infusion set (cf. exemplary product illustrations C and D), rinsing with a neutral compatible solution is recommended after each medication to avoid potential incompatibility of drug mixtures, maintain the integrity of the line and prevent microbiologic growth. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections of users.

### Notes on safe use of product

• Please observe any further instructions for use that are enclosed with this medical device or the instructions for use of any compatible medical devices used in parallel. • The product is designed for use with compatible pressure infusion pumps with generated pressures up to 200 kPa (2 bar). The product can also be used under gravity and with a pressure cuff with generated pressures up to 50 kPa (0.5 bar). Restrictions in the type of application are possible due to the product structure (e.g. products with integrated pressure-activated valve cannot be used by gravity feed). • Notes on use with pumps: the user shall ensure that the pump is approved for the administration of fluids, pharmaceuticals and parenteral nutrition solutions and calibrated with the appropriate infusion device before commissioning. Always observe the instructions for use from the pump manufacturer. Only open the roller clamp during pump operation! • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to use, existing closure plugs and disconnectable connections must be checked for tightness and retightened if necessary. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented) before using the medical device. • A fluid filter < 20 µm (according to EN ISO 8536-8) is required for infusion. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e. protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70% or isopropyl alcohol 70% in combination with chlorhexidine gluconate 2% can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Avoid excessive injection pressure during injections of solutions into the infusion device/line, especially when using small syringes. Handling errors can result in leakages.

### Use of product

Please also observe the exemplary product illustration(s) as well as the graphic brief instructions under notes on product use and all further information/images for product use on and in the product packaging.

General use (cf. exemplary product illustrations A/B/C/D):

1. Check if the sterile packaging is damaged or has been unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. In case of rigid/semi-rigid infusion containers/bottles, an infusion device with a vented closure-piercing device (1) (A/D) is required. In case of infusion bags, infusion devices with non-vented (1) (B/C) as well as vented closure-piercing devices (1) (A/D) can be used. Close the ventilation cap of vented closure-piercing devices (1) (A/D) first.
4. Close the roller clamp (6).
5. Remove the protective cap and insert the vented closure-piercing (1) (A/D) straight into the uprightly positioned rigid/semi-rigid infusion container/bottle or insert the non-vented closure-piercing device (1) (B/C) into the port of infusion bag with a rotating movement. Please ensure that the wall of the bag remains intact during piercing.

6. Fill the drip chamber (5) up to the level ring (4).
7. In the case of rigid/semi-rigid infusion containers/bottles, open the ventilation cap of the vented closure-piercing device (1) (A/D). For infusion bags, the ventilation cap remains closed.
8. Open the roller clamp (6) and prime the remaining infusion device up to the patient connection port (9) so that the product is applied without any air inclusions. Close the roller clamp (6).
9. In case of gravity application: proceed according to point 10. When using with pump: follow the pump manufacturer's instructions for use. Open the door of the infusion pump. Insert the tubing into the infusion pump. Roller clamp (6) has to be positioned on the patient connection port (9) side. Close the door of the infusion pump.
10. Remove the protective cap from the patient connection port (9) and connect the patient connection port (9) to the patient access.
11. For starting the infusion, open the roller clamp (6) and adjust the required medication and patient parameters.
12. Check the infusion at regular intervals.

Use of product configurations with connection port above/integrated into the drip chamber (cf. **exemplary product illustration C**, e.g. Drip SWAN®) for secondary line(s)/single-use syringe(s):

Use for administration with secondary line:

13. Stop the infusion and close the clamps (2) (6).
14. Disinfect the connection port (3) and connect the secondary line.
15. Open the closure clamp(s) of the secondary line and the roller clamp (9) of the infusion set.
16. For starting the infusion, check or adjust the medication and patient parameters.

For successive administration of infusion solutions through several secondary lines, the respective unused line shall be closed with the corresponding closure clamp. After each medication administration, flushing with a neutral, compatible solution is recommended.

Use for administration with pre-filled single-use syringe:

13. Stop the infusion and close the clamps (2) (6).
14. If applicable remove the product from the infusion pump.
15. Disinfect the connection port (3) and connect the single-use syringe.
16. Open the roller clamp (6) completely.
17. Slowly inject the solution from the single-use syringe into the drip chamber (5). The solution in the system is administered as a bolus. For flushing the connection port, (3) approx. 1 ml sterile air from the single-use syringe can be used.
18. Close the roller clamp (6) and open the closure clamp (2).
19. If applicable insert the product into the infusion pump.
20. For starting the infusion, check or adjust the medication and patient parameters.

Use of product configurations with multiple closure-piercing devices (cf. **exemplary product illustration D**):

The various closure-piercing devices (1) can be used for successive administration of infusion solutions. The respective unused line(s)

shall be closed with the corresponding closure clamp(s) (2). After each medication administration, flushing with a neutral compatible solution is recommended.

Regarding products equipped with Wing-Valve drip chamber (5) (D), the drip chamber can be filled or regulated in the liquid level by pressing the Wing-Valve wings together.

#### Product compatibility

• The closure-piercing device complies with EN ISO 8536-8 and can be connected with all compatible infusion containers. • All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/ EN ISO 80369-7. • The product is compatible with volumetric infusion pumps with a pressure alarm limit not exceeding 2.0 bar, that are designed for infusion sets with a nominal tubing diameter of 3.0/4.1 mm (inner/outer diameter) as flow element.

#### Residual risks

• Misconnections of the Luer connectors with medical devices of other application areas (e.g. enteral nutrition) are possible. • In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion. • The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

#### Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account.

The product is designed for a period of use of maximum 96 hours. A maximum period of use of 96 hours is recommended for non-lipid-containing solutions and 24 hours for lipid-containing solutions. Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

#### Product disposal

• The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • The closure-piercing device can be inserted into the notch on the roller clamp to avoid any injuries. • Protective medical gloves must be worn during disposal.

#### Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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