


CODAN Volumetric Infusion Pumps

A717V, A718V
A717V Plus, A718V Plus




Fact Sheet

Version Mar 2021
For Firmware version 5.50



en



Technical Data	Description	Remarks
Performance Data		
Product variants	A71XV, A71XV Plus	
Dosing units	g, mg, µg, ng, IU, mmol	Per hour, minute, kilogram, day
Flow rate	0.1 to 1,200 ml/h	Steps of 0.1 ml/h from 0.1 to 999.9 ml/h rates, steps of 1 ml/h from 1,000 to 1,200 ml/h rate.
Infusion volume	0.1 to 999.9 ml/1,000 to 9,999 ml	Steps of 0.1 ml from 0.1 to 999.9 ml volume/steps of 1 ml from 1,000 to 9,999 ml volume.
Infusion time	1 min to 99 h 59 min	Steps of 1 min
Programming	Rate in ml/h, volume in ml, duration in h and min	
Automatic rate calculation	Based on the volume and time	Rate in ml/h volume in ml, infusion duration in h and min
Bolus administration	Manually or automatically	With/without dose limit, limits via MED DB
Bolus rate/volume	0.1 to 1,200 ml/h/0.1 to 999 ml	Without infusion interruption
Purging of the infusion line	0.1 to 1,200 ml/h	
Accuracy of delivery rate	±5%	Under normal conditions according to IEC 60601-2-24
Keep Vein Open (KVO) rate	0.1 to 10.0 ml/h	Qualified administration sets (included in the pump menu)
Standby mode	1 to 60 min	Depending on infusion rate; deactivation possible
History and Interface		
History and event log	Up to 1,000 events	With timestamp
Nurse call	Connection to the central system	All A71XV variants. 24 V/0.2 A, configurable pulse, potential-free
RS-232 interface	Configuration interface	For use with ARGUSservice
Communication interface	Connection to docking stations/via docking stations	Connection to Ethernet is possible (i.e. PDMS)
Ethernet and Wi-Fi	Health Level 7 (HL7), IHE Profiles: DEC, CT, SN	A71XV Plus only, (Wi-Fi optional, cannot be retro-enabled)
Bluetooth	Secure 1 to 1 connection	For Patient ID scan and Barcode Assisted Infusion (BAI)
Operation		
Applications	Intensive care, anaesthesia, standard care, oncology, blood transfusion, parenteral nutrition, neonatal care	Dose calculation
Display type	Monochrome LCD, 98.4 x 26.2 mm, 240 x 64 Pixel	Illuminated, adjustable brightness
Night Mode	A feature to automatically dim the brightness of the device backlight and front panel LEDs.	
Display contents	Infusion rate in units and ml/h infused & target volume (VTBI), infused & remaining time, medication setup, bolus info, battery capacity, pressure level, and other menus.	
Menu overview	Intuitive and configurable for specific applications	Available in 18 languages, an overview of infusion data
Menu configuration	ARGUSservice software/temp settings by the operator	Client-specific configuration possible
Infusion lines	CODAN standard and cyto-sets and other validated sets	Tube materials: PVC, PVC-free, No DEHP/display of the set name
	A718V, A718V Plus: CODAN IVIP® dedicated set	With anti-free-flow clamp
Safety		
Pressure release	Automatic pressure release after occlusion	Configurable
Air bubble detection	Individual air bubbles	Configurable from 50 to 1,000 µl
	accumulated air bubbles	Configurable from 100 to 2,000 µl within 8 to 64 min
Drop detector	10 to 65 drops/ml	In transport mode, the alarm is less sensitive
Alarms, pre-alarms	Visual and audible signals for rapid detection	With text explanation and symbols, large alarm window
Alarm types	Occlusion (up and downstream), air bubbles, no drops and too many drops, end of infusion, standby, low battery, end of battery, maintenance due etc.	
Key lock	Activation/deactivation via numerical code	
Pressure limit for occlusion	100 to 1,200 mbar (also in mmHg, kPa, cmH ₂ O, Psi)	Configurable steps
DEERS: Medication database, and Barcode Assisted Infusion (BAI)	MED DB for A71XV and A71XV Plus	On request
	BAI with E-VTBI & Rinsing modes. A71XV Plus variants only.	Packaged bundle, optional (cannot be retro-enabled)
Power Supply		
Battery	1,800 mAh, NiMH	Maintenance-free
Battery operation	6 hours at 25 ml/h	
	5.5 hours at 25 ml/h for the Plus variant	
Battery charging time	< 6 hours in standby mode (at room temperature)	
Mains power supply	100 to 240 V, 25 VA, 50 - 60 Hz	
Product Specification		
Dimensions	190 x 160 x 130 mm (W x L x D)	Without combination clamp
Housing	ASA high-performance plastic	
Operating/medication temperature	5°C to 40°C/18°C to 30°C	
Temperature for storage/transport	0°C to 40°C/-20°C to 60°C	
Relative humidity	20% to 90% without condensation	
Atmospheric pressure	Operation:	750 to 1,060 mbar
	Storage and transport:	500 to 1,060 mbar
Weight with battery	1.78 kg	Including battery and power supply unit
Maintenance	Safety Standard Check (SSC): every 24 months	Or after 10,000 operating hours
Compliance		
IP Code	IP22	Protection against ingress of solid foreign objects with a diameter of ≥ 12.5 mm, and ingress of water in a position with 15° inclination.
Wi-Fi and Bluetooth Frequency	Single-band 2.4 GHz	
Protection class	II	Electric shock protection according to IEC 60601-1
Applied part type	CF	According to IEC 60601-1
MDD classification	IIb	
Electrical safety	IEC 60601-1, IEC 60601-2-24	
Electromagnetic compatibility	IEC 60601-1-2	
Usability	IEC 60601-1-6, IEC 62366	
Alarm systems/SW life cycle	IEC 60601-1-8, IEC 62304	
Accessories		
Docking stations	CODAN docking stations in 3 sizes for different combinations	Version M: PDMS link + power supply Version P: only power supply
Attachment system	Combination clamp	Suitable for pole and rail systems
Software	ARGUSservice	Program for simple maintenance and configuration



Manufacturer
CODAN ARGUS AG
 Oberneuhofstrasse 10 · 6340 Baar, Switzerland
 Tel. +41 (0) 41 785 09 44
 www.codancompanies.com
The decisive connection



SIMPLIFIED EU DECLARATION OF CONFORMITY

Hereby, CODAN ARGUS AG (Baar) declares that the radio equipment type A71XV Plus is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at www.codanargus.com