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de

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Xtrans® electrical equipment

en

Please read this carefully before use.

Indication for use

The Xtrans® electrical equipment is intended for invasive measurement of the blood pressure in different compartments of the circulatory system.

Directions for use

A short mechanical and electrical functional test must be carried out by a qualified member of staff before each use. If one or more components of the reusable equipment does not look perfect or function smoothly, this must be replaced. The equipment must be checked at regular intervals as follows:

- Visual check that the equipment is undamaged
 - Electrical and mechanical functional check of the reusable equipment
- CODAN pvb does not accept any liability if these directions for use are not observed.

Product description

Xtrans® pressure transducers can only be connected to the monitors in conjunction with re-usable contactboards and appropriate monitor connectors. They are made from high quality materials and are designed to give long service in hospitals if handled and cleaned correctly.

The electrical operation of contactboards and monitor connectors should be checked by hospital staff every time a new pressure monitoring system is installed (zero point balancing of the pressure transducer). If zero point balancing cannot be carried out, despite repeated attempts and correct connections, the contactboard and/or monitor connector and, if necessary, the pressure transducer, should be replaced.

Assembly/dismantling

The contactboards (which are made up of two housing sections) can be dismantled and subsequently reconnected to enable the easy cleaning of them. Before removing the top section, the plastic clips at the side have to be pressed and the top removed. The four fold contactboard also has a locking screw fitted to the back of the housing, which must be completely removed before dismantling.

After cleaning, the top housing section must be placed onto the bottom section again and locked firmly in place on both sides with an audible “click”. In the case of the four fold contactboard, the locking screw must then be firmly screwed into the bottom section of the board again (only tighten per hand, otherwise there is a danger of over-tightening).

If the top section of the housing cannot be properly cleaned or is damaged, it can be replaced by a new one. Top sections can be re-ordered separately.

Cleaning

The re-usable components should be wiped clean with a surface disinfectant and a soft cloth. Use a commercially available alcohol based disinfectant.

If contactboards become heavily soiled they may be placed in instrument disinfection solution (containing no aldehydes) in accordance with the manufacturer’s instructions. Afterwards, the contactboards must be rinsed in pure alcohol and completely dried. Always ensure that the cable connector of the board is not placed in the solution with the contactboard. If this does happen, the connector must be completely dried and checked by the medical technician before it is used again.

The contact pins of the contactboards can also be cleaned with a soft brush, (e.g. toothbrush) or sponge, after having removed the top section of the housing.

Cleaning restrictions

- Cleaning the contact pins with hard objects, such as a knife or a screw-driver, may cause damage, which can result in malfunction or failure of the entire measuring system.
- The electrical components are not to be cleaned by ultrasound under any circumstances.
- Hot cleaning procedures at temperatures above 65°C (e.g. in the washing machine) must not be used.
- The contactboard and monitor connector must not be electrically operated during cleaning.
- The following agents are not to be used for cleaning: spirits, benzenes, aldehydes, ethers, phenols, acids of any kind, caustic substances, phenol-based disinfectants and peroxide compounds.
- Contact boards and monitor connectors should only briefly come into contact with disinfectant for cleaning purposes. Avoid under all circumstances storing in a disinfectant bath.

Restrictions of use

- Contact board and monitor connector cables should not be kinked or fixed with cable ties or routed over sharp edges.
- Don't disconnect the plug connection between the contactboard and the monitor connector by pulling on the electric cable.
- Depending on the type and version, the plug connection between the contactboard and monitor connector can be equipped with a push-pull lock. This connection is released by **pushing** and then **pulling** the housing of the connector, and must not be disconnected by twisting.

Storage

If contactboards and monitor connectors are not used, they must be stored in a dry, dark place (cupboard) at room temperature.

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