IBPM DPT-6000 (DPT-6200, DPT-6400, DPT-6990) DPT-9000/Xtrans® (DPT-9400 & DPT-9990)

de	Gebrauchsanweisung – Deutsch	6
en	Instructions for use – English	24
fr	Notice d'utilisation – Français	42
it	Istruzioni per l'uso – Italiano	62
sv	Bruksanvisning – Svenska	82
no	Bruksanvisning – Norsk	100
es	Instrucciones de uso – Español	118
pt	Instruções de utilização – Português	136
nl	Gebruiksaanwijzing – Nederlands	156
da	Brugsanvisning – Dansk	174

DPT-6000 (DPT-6200, DPT-6400 & DPT-6990)
DPT-9000/Xtrans® (DPT-9400 & DPT-9990)

Spezifikationen | Specifications | Caractéristiques | Specifiche | Specifikationer | Spesifikasjoner | Características técnicas | Especificações | Specificaties |
 Spezifikationen

	2 V – 15 V (≤ 5 kHz)		0 °C – +40 °C
	-30 mmHg – +360 mmHg		15% – 90%
	4000 mmHg (max. 1 s)		600 hPa – 1060 hPa (-300 m – 4000 m)
	5 μ V/V/mmHg (± 1 %)		0,0001 kWh (max. 0,6 W)
	± 4 % / ± 4 mmHg		5 kV (<10 s)
			400 Ω – 10/+20 %
			356 Ω ± 10 %


EN 60601-1:2006 + Cor.:2010 + A1:2013 + AC:2014 + A1:2013/AC:2014 + A12:2014 + A2:2021
 EN 60601-1-2:2015
 EN 60601-1-6:2010 + A1:2015 + A2:2021
 EN 60601-1-12:2015 + A1:2020
 EN 60601-2-34:2014

EN ISO 8536-4:2020
 EN ISO 8536-9:2015
 EN ISO 8536-10:2015
 EN 1789:2020

EN 60529:1991 + A2:2013/AC:2019-02
 EN ISO 80369-1:2018
 EN ISO 80369-7:2021
 EN 13718-1:2014 + A1:2020

Medical Device Regulation (EU) 2017/745 (MDR)
 Directive 2011/65/EU & Directive (EU) 2015/863 (RoHS)
 Regulation (EC) 1907/2006 (REACH)
 EN 62366-1:2015 + AC:2015 + A1:2020
 ASTM F2503-20

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en **Instructions for use – English** **24–41**

1	Overview and foreword	24
2	Safety and protective measures	24
3	Intended use	25
4	Device description	26
5	Combination with medical devices from other companies	29
6	Provision	30
7	Preparation	31
8	Filling	31
9	Initial use	34
10	Zero calibration	36
11	Transmission quality check	36
12	Disinfection instructions	36
13	Duration of use	36
14	Technical specifications and operating conditions	37
15	Storage and transportation	39
16	Residual risks	39
17	Replacing the pressure monitoring set/sterile accessories	39
18	Ending use of a pressure monitoring set	39
19	Disposal	41
20	Customer services	41
21	Troubleshooting and fault rectification	41

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Instructions for use – English

1 Overview and foreword

Read this information carefully before use and keep it for future reference.

These Instructions for Use (IFU) contain important information about the intended use of the pressure monitoring set manufactured by CODAN pvb Critical Care with pressure transducers from the **DPT-6000** (DPT-6200, DPT-6400 & DPT-6990) or **DPT-9000/Xtrans®** (DPT-9400 & DPT-9990) series and sterile accessories. The configuration of this pressure monitoring set may differ from the device configurations on page 3 of these Instructions for Use, and therefore the precise guidelines and procedures must be determined by the responsible organisation. The responsible organisation (e.g. the ward manager, the responsible doctor or the clinical institution) designates the entity responsible for using the pressure monitoring set, the sterile accessories and/or the entire monitoring system in which the CODAN pvb Critical Care pressure monitoring set and/or sterile accessories are combined with other medical devices for invasive blood pressure monitoring.

Any serious incidents that occur during use of, or in connection with, CODAN pvb Critical Care devices and that have led, could have led or might lead to the death or temporary or permanent serious deterioration in a patient's, user's or other person's state of health, or a serious risk to public health must be reported to the manufacturer CODAN pvb Critical Care and to the competent federal authority.

Information labelled with one of the following symbols indicates a hazard,

If this hazard is not avoided, it can result in damage to property or the environment;

Hazard with a low risk level which, if not avoided, can result in minor or moderate injury;

Hazard with a medium risk level which, if not avoided, can result in death or severe injury;

Hazard with a high risk level which, if not avoided, can result in death or severe injury.

Further explanations of the symbols used on the packaging and in the Instructions for Use can be found in the tables on **Page 192 and subsequent pages**.

2 Safety and protective measures

Read and comply with all of the following safety instructions. If these instructions are not observed/complied with, there are significant safety risks for patients, users and third parties.

- 2.1  Comply with all instructions (including safety instructions and warnings) for the intended use of all **(third-party) devices** used in combination with the pressure monitoring set in order to guarantee the safe use of the pressure monitoring set and/or accessories. Prior to use, carefully read all the accompanying information provided by the manufacturers of the (third-party) devices.
- 2.2  **Before performing any therapeutic interventions on a patient** which are to be initiated on the basis of blood pressure values displayed on the control monitor, the plausibility of the measured values must be checked to rule out a device error or malfunction as the cause of measured values exceeding or falling below the set limit values. To do this, check the other patient displays (e.g. electrocardiography, oxygen saturation) and vital functions (e.g. pulse, respiration check). In addition, inspect the entire set-up of the pressure monitoring system at each change of shift (**see section 21, page 41 – Troubleshooting and fault rectification**).
- 2.3  **Exposed contacts** of electrical accessories must not be touched at the same time as the patient while the accessories are making an electrically conductive connection to a control monitor. **Risk of electric shock!** Wear medical protective gloves when preparing, using and disassembling the electrical accessories and ensure that an electrically conductive connection cannot be established between the patient and the exposed contacts.
- 2.4  **Spilling of liquids**, in whatever quantity, onto electrical contacts can lead to the loss of basic safety, including the essential performance characteristics. Make sure that no liquid can wet the electrical contacts of the accessories and/or pressure monitoring set during preparation, use or disassembly. In the event of liquid contact, replace the pressure monitoring set immediately, clean and dry the electrical accessories and check their function.
- 2.5  When **using a defibrillator**, any contact with the pressure monitoring system is prohibited. **Risk of electric shock!** Avoid any contact with the monitoring system when resuscitation is being performed using a defibrillator.
- 2.6  The pressure transducer types DPT-6200, DPT-6400, DPT-6990, DPT-9400 & DPT-9990 do not have a **flushing system** or equivalent component with which rapid flushing could be carried out and/or continuous flushing could be set. The responsible organisation must make sure that the pressure monitoring set is always free of **air pockets** and blood residue, and that the patency of the patient access port is maintained at all times. Air pockets in the pressure monitoring set negatively affect the signal transmission and can lead to an air embolism inside the patient's vascular system. **Clotted blood (coagulum)** inside the tubing system can cause the patient access port to become clogged and a loss of function of the pressure monitoring set. Moreover, infiltration of coagulum into the patient's circulatory system can cause vascular occlusion. The (incipient) **occlusion of the patient access port** can lead to deterioration or interruption of the signal transmission.

2.7

Before and during emergency use of the device within appropriate emergency vehicles (see 3.8 Intended working environment) during professional primary (emergency) and secondary (intensive care and ambulance) transport to or between clinical healthcare facilities, it is important to ensure that the Luer connectors and the sampling ports of the pressure monitoring set are free of dust, fluff or dirt, which can lead to leaks and/or contamination. To prevent soiling or ingress/leakage of liquids, air and/or pathogens, open-end Luer connectors must be closed with non-vented protective caps.

3 Intended use

3.1 Intended purpose

The pressure monitoring set is intended for use on a patient requiring continuous monitoring of their blood pressure via one or more vascular access ports (IBPM).

3.2 Indications

The use of CODAN pressure monitoring sets is always indicated whenever the medical condition requires, at the discretion of the treating physician, continuous monitoring of the blood pressure based on the patient presentation, specific diagnosis, diagnostic findings and/or pathophysiological considerations, among other things.

3.3 Clinical benefit

The pressure monitoring set allows for measurement of the blood pressure level in one or more compartments of a patient's circulatory system with subsequent signal transmission to a patient monitor to make possible the visualisation of the measured values on a screen (in the form of digits and/or curves).

Therefore, the clinical benefit of the pressure monitoring set is to enable the user to

- continuously monitor and document one or more blood pressure value(s) of a patient in order to be able to detect if the defined limit values have been exceeded and/or irregularities due to (pre-existing) illnesses or blood loss exist, and be able to react to them with appropriate treatment;
- make a diagnosis of (pre-existing) illnesses of the cardiovascular system and/or evaluate therapies with drugs that have a cardiovascular effect.

3.4 Contraindications

- Allergies to residues from an ethylene oxide [EO] sterilisation process
- Electrolyte imbalances (if necessary, the saline solution must be adjusted prior to use)
- All contraindications to an arterial or venous puncture (e.g. blood coagulation disorders or occlusive diseases)

3.5 Intended patient group

Pressure monitoring sets with **DPT-6200, DPT-6400, DPT-6990, DPT-9400 and DPT-9990** pressure transducers are used on critically ill and/or ventilated patients from birth, without any weight restriction, who require monitoring of their blood pressure levels in the intravascular system.

3.6 Intended user profile

The pressure monitoring set may only be used by healthcare professionals who have been trained in accordance with the instruction for Use for this medical device. The user must not have any mental or physical disabilities (e.g. colour blindness) that could prevent him/her from using the pressure monitoring set and/or accessories as intended, in accordance with the directions in these Instructions for Use.

3.7 Type and duration of intended physical contact

Depending on the configuration of the pressure monitoring set and the therapy administered for the specific medical condition/disease progression, the pressure monitoring set is connected to one or more vascular access line(s). This/these access line(s) is/are connected to one or more arterial or venous vessels in one or more area(s) of the patient's body, mainly the extremities and/or neck. The extended exposure time, i.e. the cumulative sum of individual, multiple or repeated contact times, may exceed 24 hours but not 30 days.


3.8 Intended working environment

The pressure monitoring set is used within clinical healthcare institutions for monitoring, diagnostics, anaesthesia and intensive care.

Additionally, the pressure monitoring sets with pressure transducers from the **DPT-6000** (DPT-6200, DPT-6400 & DPT-6990) or **DPT-9000/Xtrans®** (DPT-9400, DPT-9990) series are intended for use within a defined electromagnetic environment (see Section 14, page 37 – Technical Specifications and Operating Conditions).

3.9 Intended monitoring system (physical functionality)

The physical functionality of the monitoring system consists of conducting a pressure transmitted via intravascular access ❶ and pressure measurement lines ❷ through a fluid column ❸ to a pressure transducer ❹ and converting it into electrical signals via a Wheatstone full bridge. With the aid of reusable electrical accessories ❺ these signals are transmitted to a control monitor ❻.

Instructions for use – English

4 Device description

The medical devices from the manufacturer CODAN pvb Critical Care are developed, manufactured and placed on the market in accordance with the fundamental safety and performance requirements of the Regulation (EU) on Medical Devices 2017/745 (MDR), the Regulation 1907/2006 (REACH), the Directives 2011/65/EU and (EU) 2015/863 (RoHS) and the requirements of the European Standards EN 62366-1:2015 + AC:2015 + A1:2020, and they carry a CE marking. The electrical accessories and pressure monitoring sets with pressure transducers from the **DPT-6000** or **DPT-9000/Xtrans®** series meet the requirements of the European Standards EN 60529:1991 + A2:2013/AC:2019-02, EN 60601-1:2006 + Cor.:2010 + A1:2013 + AC:2014 + A1:2013/AC:2014 + A12:2014 + A2:2021, EN 60601-1-2:2015, EN 60601-1-6:2010 + A1:2015 + A2:2021 and EN 60601-2-34:2014. The Luer connectors of the sterile accessories and pressure monitoring sets are manufactured in accordance with the requirements of the Standards EN ISO 80369-1:2018 and EN ISO 80369-7:2021. CODAN pvb Critical Care medical devices are delivered in protective/transport packaging validated according to EN ISO 11607-1:2020 and EN ISO 11607-2:2020 to protect the devices from mechanical damage and general soiling during storage and transportation. Moreover, the pressure monitoring sets and sterile accessories are delivered in a single or double sterile barrier system (sterile packaging) validated according to the specified standards to protect against contamination with pathogens after the sterilisation process (in accordance with EN ISO 11135:2014 + A1:2019) and to maintain the sterile condition of the devices.

For a better understanding of the following information about the medical devices from the manufacturer CODAN pvb Critical Care, please refer to the configuration overviews on page 3 of these Instructions for Use.

4.1 Pressure monitoring set

The pressure monitoring set is a sterile, active and portable class IIb medical device. The essential performance characteristics of the pressure monitoring set in accordance with EN 60601-2-34:2014 include accuracy of the pressure measurements, defibrillation-proof and resistance to electrosurgery interference. The basic safety features of the pressure monitoring set include:

- Leak tightness;
- Pressure resistance;
- Tensile/electrical strength;
- Disinfectant resistance;
- Protection against mechanical stresses;
- Electrical safety.

A fully configured pressure monitoring set comprises

- Components for connecting the pressure monitoring set to a liquid medium for filling and irrigation **(I & III)**;
- One or more pressure transducers for converting one or more physical pressures (blood pressure) into one or more electrical signals **(II)**;
- Components for transmitting one or more physical pressures (blood pressure) to one or more pressure transducers **(IV)**.

The smallest possible configuration for a pressure monitoring set is a single pressure transducer **(II)** of the **DPT-6000** (DPT-6200, DPT-6400 & DPT-6990) or **DPT-9000/Xtrans®** series (DPT-9400 & DPT-9990). All other sections **(I, III & IV)** are optional and must be supplemented with the sterile accessories from the manufacturer CODAN pvb Critical Care (**see section 4.2, page 27 – Sterile accessories**). The following table lists the components of a pressure monitoring set and their main functional and operating features, based on the configuration overviews on page 3 of these Instructions for Use.

Pressure monitoring lines are required to connect different components and for pressure transmission. The pressure monitoring lines can be made from PVC (polyvinyl chloride), PE (polyethylene), or PUR (polyurethane) with different Luer connectors, lengths and (internal/external) diameters and they can be installed at various positions on the pressure monitoring set. Due to their versatility, the pressure monitoring lines are not explained in further detail in the following table. More information about the pressure monitoring lines from the manufacturer CODAN pvb Critical Care will be provided by the distributor and/or manufacturer upon request.

Overview (p. 3)		Component				
Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
-	-	Protective/transport packaging	-	<ul style="list-style-type: none"> • Protects the medical device from mechanical damage and general soiling during storage and transportation • Identifiability and traceability 	1	A sterile barrier system with external protective packaging is labelled with the symbol 20
-	-	Sterile packaging (sterile barrier system)	-	<ul style="list-style-type: none"> • Protection against the ingress of pathogens • Maintaining the sterile condition • Identifiability and traceability • Easy and residue-free opening 	1–2	<ul style="list-style-type: none"> • A single sterile barrier system is labelled with the symbol 21 • A double sterile barrier system is labelled with the symbol 22
A	III	Filling and irrigation system with 4-way stopcock	-	Connecting the pressure monitoring set to a saline bag for filling and rinsing the pressure monitoring set with saline solution	1	-
A	III	Drip chamber macro (vented or non-vented)	Filling and irrigation system with 4-way stopcock	<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) • Setting the liquid level • Flow control (20 drops per millilitre) • Particle filtration (15 µm) 	1	-
		Drip chamber micro (vented or non-vented)		<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) • Setting the liquid level • Flow control (60 drops per millilitre) • Particle filtration (15 µm) 		
		Spike (vented or non-vented)		<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) 		
A	III	Roller clamp	-	Flow control (opening and closing the flow channel)	1	-
A	III	4-way stopcock	-	Setting of four possible flow directions (360°)	1	-

Overview (p. 3)		Component				
Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
B	I	Filling and irrigation system	-	Connecting the pressure monitoring set to a saline bag for filling and rinsing the pressure monitoring set with saline solution	1	-
B	I	Drip chamber macro (vented or non-vented)	Filling and irrigation system	<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) • Setting the liquid level • Flow control (20 drops per millilitre) • Particle filtration (15 µm) 	1	-
		Drip chamber micro (vented or non-vented)		<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) • Setting the liquid level • Flow control (60 drops per millilitre) • Particle filtration (15 µm) 		
		Spike (vented or non-vented)		<ul style="list-style-type: none"> • Connection to a saline bag (vented or non-vented) 		
B	I	Roller clamp	-	• Flow control (opening and closing the flow channel)	1	-
A-B	II	Pressure transducer types DPT-6200, DPT-6400 or DPT-9400 (Xtrans®)	-	Conversion of a physical pressure (blood pressure) into an electrical signal	1	-
A-B	II	1-way stopcock	DPT-6200, DPT-6400, DPT-9400	Flow control (opening and closing the flow channel)	1	-
A-B	II	Electrical contacts		• Connecting the pressure transducer to an electrical accessory for signal transmission	1	-
A-B	II	3-way zeroing stopcock		<ul style="list-style-type: none"> • Zero calibration of the pressure transducer • Setting of three possible flow directions (180°) 	1	-
C	II	Pressure transducer types DPT-6990 or DPT-9990 (Xtrans®)	-	Conversion of a physical pressure (blood pressure) into an electrical signal	1	-
C	II	Electrical contacts	DPT-6990, DPT-9990	Connecting the pressure transducer to an electrical accessory for signal transmission	1	-
C	II	1-way stopcock (pressure transducer addition)	-	Flow control (opening and closing the flow channel)	1-2	<ul style="list-style-type: none"> • Only in combination with pressure transducers type DPT-6990 or DPT-9990 • The pressure transducer additions (1-way stopcock & 4-way stopcock) can be present in the ratios 0:2, 1:1 or 2:0
C	II	4-way stopcock (pressure transducer addition)	-	<ul style="list-style-type: none"> • Zero calibration of the pressure transducer • Setting of four possible flow directions (360°) 	1-2	<ul style="list-style-type: none"> • Only in combination with pressure transducers type DPT-6990 or DPT-9990 • The pressure transducer additions (1-way stopcock & 4-way stopcock) can be present in the ratios 0:2, 1:1 or 2:0
A, C	IV	4-way stopcock	-	Setting of four possible flow directions (360°)	1	-
B	IV	4-way stopcock combination	-	Setting of 16 possible flow directions (360°)	1	-

4.2 Sterile accessories

Sterile accessories (class IIa medical devices) are individual, sterile packaged components of a pressure monitoring set from the manufacturer CODAN pvb Critical Care, not including individual pressure transducers of the DPT-6000 or DPT-9000/Xtrans® series and consumables. The sterile accessories are used to supplement, extend or change the configuration and/or components of a pressure monitoring set. These include filling systems, blood sampling systems (for pressure transducer types DPT-6003, DPT-6009, DPT-9003 and DPT-9009), pressure monitoring lines and one-way or multi-way stopcocks with or without connected pressure monitoring lines.

Additions, extensions or changes to the configuration and/or components of a pressure monitoring set are always the responsibility of the responsible organisation and must only be implemented using the sterile accessories from the manufacturer CODAN pvb Critical Care. Additions, extensions or changes to the configuration and/or components using third-party devices are prohibited, because the basic safety and essential performance characteristics of the pressure monitoring set cannot then be guaranteed. Moreover, additions, extensions or changes to the configuration and/or components of the pressure monitoring set between the sampling port and patient access port may mean that sufficient volume can no longer be drawn using a reservoir (NBSS, NBSS^{XL}, ABSS & ABSS^{XL} configurations) for undiluted blood samples to be taken at the sampling port.

A list of all the available sterile accessories and other information will be provided by the distributor and/or manufacturer upon request.

4.3 Non-sterile accessories

Reusable mounting accessories, comprising holders, mounting plates and modular holder systems (mounting heads, claws and poles), are non-sterile class I medical devices used for assembling pressure monitoring sets and electrical accessories on round poles or rail systems.

Single-use mounting accessories, comprising velcro straps and fasteners, are non-sterile class I medical devices used for assembling pressure monitoring sets, reservoirs and electrical accessories on the patient's arm or for securing pressure monitoring lines.

Reusable electrical accessories are non-sterile class I medical devices that enable the transmission of signals between one or more pressure transducers and a control monitor. These include contact boards and monitor connectors for DPT-9000/Xtrans® and interface cables and cable management for DPT-6000.

A list of all the available non-sterile accessories and other information will be provided by the distributor and/or manufacturer upon request or can be downloaded from: <https://www.codancompanies.com/downloadcenter/>

Instructions for use – English

4.4 Consumables

Consumables are required for the intended use of a pressure monitoring set on a patient within a monitoring system and may have to be used up, disposed of or replaced before or during use. Illustrations in these Instructions for Use labelled with the following symbols inform the user about

➔ a used material that must be disposed of according to regulations (see section 19, page 41 – Disposal);

➔ a new consumable that is required for the intended use of a pressure monitoring set and must be provided.

Since the configuration of this pressure monitoring set may differ from the device configurations on page 3 of these Instructions for Use, the responsible organisation must identify all the consumables for the pressure monitoring set prior to use and provide a sufficient quantity of these consumables to the user (see section 5, page 29 – Combination with medical devices from other companies).

4.4.1 Protective caps

Non-vented protective caps are closures with one or more Luer connectors that are attached to the free Luer ports of a pressure monitoring set and/or sterile accessory to protect against the ingress/leakage of liquids, air and/or pathogens before and/or during use.

Please note that a pressure monitoring set and/or sterile accessory may be delivered with **vented protective caps** attached to the Luer ports. These vented caps must be replaced with **non-vented protective caps** prior to use to prevent the ingress/leakage of liquids, air and/or pathogens.

The list below shows all the protective caps that may be found in a CODAN pvb Critical Care pressure monitoring set when it is delivered.

Protective caps from CODAN pvb Critical Care		Colour coding	
Illustration	Luer type	Vented	Non-vented
	Male	Red, blue	Neutral
	Male	Red	Neutral
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Male & female	–	Neutral, blue, red

4.4.2 Consumables of the monitoring system

The list below shows all the consumables required for the intended use of the pressure monitoring set within a monitoring system:

- Cannulae and/or catheter(s) (vascular access);
- Bag with isotonic (0.9%) saline solution (for pressure monitoring sets with a filling and Irrigation system);
- Disposable syringes with isotonic (0.9%) saline solution (for pressure monitoring sets without a filling and irrigation system).

4.4.3 Other consumables

The list below shows all the other consumables required for the intended use of the pressure monitoring set:

- Disposable syringe(s) and/or equivalent sampling containers for volume shift;
- Disinfection wipes/swabs and/or disinfectant (see section 12, page 36 – Disinfection instructions);
- Medical protective gloves.

4.5 Colour coding

Take note of the colour coding of the pressure monitoring lines and one-way or multi-way stopcocks. Red stopcocks and pressure monitoring lines are intended for arterial, blue for venous, yellow for pulmonary-arterial and green for left atrial applications. Neutral stopcocks and pressure monitoring lines may be used individually for these applications. Components with neutral colour coding must be suitably labelled prior to use to prevent mix-ups or incorrect use, for example using coding pins or colour-coded stickers.

Colour coding	Application/catheter position	Examples of puncture sites
Red	Arterial	Radial artery, ulnar artery, brachial artery, femoral artery, dorsalis pedis artery
Blue	Venous	Internal jugular vein, external jugular vein, cephalic vein, basilic vein, subclavian vein
Yellow	Pulmonary-arterial	Internal jugular vein
Green	Left atrial	Right femoral vein

5 Combination with medical devices from other companies

Defined below are the minimum requirements for medical devices from other companies that must be met to ensure the safe and intended use of a pressure monitoring set with pressure transducers of the **DPT-6000** (DPT-6200, DPT-6400 & DPT-6990) or **DPT-9000/Xtrans®** (DPT-9400 & DPT-9990) series and/or CODAN pvb Critical Care accessories within a monitoring system. Moreover, it must be proven that the medical devices from other companies were manufactured and authorised in accordance with all applicable device-specific standards and that they enable the pressure monitoring set and/or accessories to be used in accordance with the operating conditions specified in these Instructions for Use (**see section 14, page 37 – Technical specifications and operating conditions**).

Any deviation from the listed requirements is always the responsibility of the responsible organisation. The responsible organisation must ensure that combining a pressure monitoring set and/or accessories from the manufacturer CODAN pvb Critical Care GmbH with one or more medical devices from another company, including all connections/interfaces, is safe and does not cause any reciprocal interference, which could lead to reduction or loss of the basic safety and/or essential performance characteristics of one or more devices within the monitoring system.

5.1 General requirements

Medical devices from other companies must be developed, manufactured and placed on the market in compliance with the fundamental safety and performance requirements of the Medical Devices Regulation (EU) 2017/745 (MDR) or the Medical Devices Directive 93/42/EEC (MDD), and they must also be CE marked.

5.2 Medical electrical equipment

Any user who combines medical electrical equipment from another company with a pressure monitoring set from CODAN pvb Critical Care, with or without the aid of accessories, is the System Configurator and thereby responsible for ensuring that the monitoring system verifiably meets the specified requirements.

5.2.1 Control monitors

Control monitors intended to be combined with CODAN pvb Critical Care electrical accessories and pressure monitoring sets within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standards EN 60601-1, EN 60601-1-2, EN 60601-1-6 and EN 60601-2-34.

5.3 Sterile and non-sterile accessories

Sterile and non-sterile accessories are used for **additions, extensions or changes to the configuration and/or components of the pressure monitoring set** and for signal transmission or attachment (**see section 4, page 26 – Device description**). Only CODAN pvb Critical Care accessories are authorised for use with a pressure monitoring set with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series. In addition, the accessories should only be used in combination with CODAN pvb-Critical Care pressure monitoring sets. Non-compliance can lead to reduction or loss of the basic safety and/or essential performance characteristics.

5.4 Consumables

Consumables are required for the intended use of a pressure monitoring set on a patient within a monitoring system and may have to be used up, disposed of and replaced before or during use (**see section 4.4, page 28 – Consumables**).

5.4.1 Consumables with Luer connectors

The Luer connectors of consumables (protective caps, catheters, cannulae or disposable syringes) must be verifiably manufactured and authorised in accordance with the current requirements of the standards EN ISO 80369-1 and EN ISO 80369-7 with a 6% (Luer) cone for intravascular or hypodermic applications.

Consumables with Luer connectors that are not manufactured and authorised in accordance with the above-mentioned standards may allow the ingress/leakage of liquids, air and/or pathogens when priming and/or using a pressure monitoring set. Before using the consumables, verify that they are manufactured and authorised in accordance with the specified standards.

5.4.1.1 Catheters and cannulae

Catheters and cannulae intended to be combined with CODAN pvb Critical Care pressure monitoring sets within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standards EN ISO 10555-1 (catheters) or EN ISO 7864 (cannulas).

5.4.1.2 Disposable syringes

Disposable syringes intended to be combined with CODAN pvb Critical Care pressure monitoring sets for volume shifts must verifiably meet the current requirements of the European standard EN ISO 7886-1.

5.4.2 Other consumables

5.4.2.1 Bags containing isotonic (0.9%) saline solution

Saline bags intended to be combined with CODAN pvb Critical Care pressure monitoring sets with a filling and Irrigation system within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard EN ISO 15747.

5.4.2.2 Medical protective gloves

Medical protective gloves intended to be worn before, during and after use of CODAN pvb Critical Care accessories and/or pressure monitoring sets must verifiably meet the current requirements of the European standard EN 455-1.

Instructions for use – English

5.5 Other (third-party) devices

5.5.1 Mounting poles & rails

Mounting rails intended to be combined with CODAN pvb Critical Care reusable mounting accessories within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard DIN EN ISO 19054 and have a maximum width of 25 mm and a maximum thickness of 10 mm. The tubing of the mounting poles or infusion stands must have a diameter of between 15 mm and 35 mm. The rectangular profile of the mounting poles or infusion stands must measure 30 x 30 mm.


6 Provision

To prepare correctly, follow the instructions below.

6.1 Unpacking

6.1.1  Pressure monitoring sets and sterile accessories with **damaged protective and/or sterile packaging (sterile barrier system)** must not be used. They must be disposed of because the basic safety and essential performance characteristics and sterility of the contents can no longer be guaranteed and contamination with pathogens is possible (**see section 19, page 41 – Disposal**). Before unpacking, make sure that the protective and sterile packaging is undamaged and inspect the integrity of the sealed seams.

6.1.2  Pressure monitoring sets and sterile accessories must not be used **after the expiry date**. They must be disposed of because the basic safety and performance characteristics and sterility of the contents can no longer be guaranteed and contamination with pathogens is possible (**see section 19, page 41 – Disposal**). Before unpacking, make sure that the expiry date stated on the protective and/or sterile packaging of the device has not passed.

6.1.4  Using the monitoring system and its components without observing the **hygiene guidelines** specified by the responsible organisation may lead to the spread/transmission of pathogens and thus pose a danger to patients, users and third parties. The entire monitoring system and its components should only be used in accordance with the specified hygiene guidelines and wearing **medical protective gloves**.


6.1.5  **After the sterile packaging is opened**, the pressure monitoring set and/or sterile accessories must immediately be prepared aseptically, filled and used in order to minimise the time that the device is exposed to potential pathogens. The responsible organisation must determine the appropriate length of time that may pass after the device is removed from the sterile packaging and before it is used.

If the sterile packaging is opened by mistake, the pressure monitoring set and/or sterile accessory must be aseptically processed, primed and used immediately in accordance with the specifications of the responsible organisation. Otherwise it should be disposed of because sterility cannot be guaranteed and contamination with pathogens is possible (**see section 19, page 41 – Disposal**).

6.1.6  Pressure monitoring sets and accessories with **visible defects** and/or **mechanical damage** must not be used. They must be disposed of because the essential performance characteristics and basic safety can no longer be guaranteed (**see section 19, page 41 – Disposal**). After unpacking, make sure that the devices are undamaged.

6.1.7  **Loose Luer connections** may allow the ingress/leakage of liquids, air and/or pathogens during filling and/or use. Before connecting the pressure monitoring set to the patient, check that all the Luer connections are tight and, if necessary, tighten them by hand (without tools). It is prohibited to use tools/instruments to tighten the Luer connections because they may break if excessive force is applied, resulting in an open flow channel.

6.2 Providing of the monitoring system


7 Preparation

To prepare correctly, follow the instructions below.

- 7.1.1  The pressure monitoring set is not approved for the **administration of medicinal products and/or other liquid media (e.g. contrast media)** because these may damage the components of the pressure monitoring set, meaning that the basic safety and essential performance characteristics are not guaranteed. The use of 0.9% saline solution (possibly heparinised) to fulfil the intended use is permitted.

7.2 Connecting a pressure monitoring set with filling system to a saline bag (DPT-6200, DPT-6400, DPT-9400)


- 7.2.1 Please note that the instructions in section 7.2 are only applicable for pressure monitoring sets with a filling and irrigation system and pressure transducers of the types DPT-6200, DPT-6400 & DPT-9400. Pressure monitoring sets without a filling and irrigation system should be primed using disposable syringes (see section 7.3, page 31 – Connecting a pressure monitoring set without a filling and irrigation system to a disposable syringe (DPT-6990 & DPT-9990)).

- 7.2.3 When filling the pressure monitoring set, make sure that no more than 50 mmHg pressure is applied to the connected saline bag because otherwise tiny air bubbles can form inside the pressure monitoring set, which may negatively affect the pressure transmission and cause an air embolism in the patient's vascular system. It is always recommended to gravity fill the pressure monitoring set to prevent the formation of tiny air bubbles.

7.3 Connecting a pressure monitoring set without filling system to a disposable syringe (DPT-6400 & DPT-9990)

- 7.3.1 Please note that the instructions in section 7.3 are only applicable for pressure monitoring sets without a filling and irrigation system and for pressure transducers of the types DPT-6400 & DPT-9990. Pressure monitoring sets with a filling and irrigation system should be primed using a saline bag (see section 7.2, page 31 – Connecting a pressure monitoring set with filling and Irrigation system to a saline bag (DPT-6200, DPT-6400, DPT-9400)).

8 Filling

Follow the instructions below to correctly prime the pressure monitoring set.

8.1 General information

- 8.1.1  **Air pockets** in the pressure monitoring set negatively affect the signal transmission and can lead to an air embolism inside the patient's vascular system. All the components of a pressure monitoring set must be completely primed with saline solution and free of air bubbles before the pressure monitoring set is connected to a patient access port. Any remaining air pockets in the pressure monitoring set must be loosened by carefully tapping the affected components and then flushed out of the pressure monitoring set. Moreover, throughout the entire duration of use of the pressure monitoring set, it must be ensured that all the components are free of air pockets. To do this, follow the illustrations and instructions in section 8.


- 8.1.2  → When using a pressure monitoring set with filling system and connected saline bag, this picture tells the user to prime/flush a component with saline using gravity or gentle pressure on the bag (50 mmHg max.).

- 8.1.3  → When using a pressure monitoring set without filling system, this picture tells the user to fill/flush a component with saline solution using the connected disposable syringe.


Instructions for use – English

ENGLISH




8.2 Filling a pressure monitoring set with priming system (Configuration A)




8.3 Filling a pressure monitoring set with priming system (Configuration B)

8.3.3




8.3.4




8.4 Filling a pressure monitoring set without priming system (Configuration C)




8.4.1

8.4.2

8.4.3



Instructions for use – English

9 Initial use



Follow the instructions below for the correct initial use of the pressure monitoring set.

9.1 Assembly

9.1.1 Assembly of a pressure monitoring set with DPT-6000 (DPT-6200, DPT-6400 & DPT-6990) pressure transducer

9.1.2 Assembly of a pressure monitoring set with DPT-9000/Xtrans® (DPT-9400 & DPT-9990) pressure transducer

9.2 Positioning

9.3 Patient connection

9.3.1 General


9.3.1.1  Before using the pressure monitoring set on a patient, make sure that all the vented protective caps have been disposed of and replaced with non-vented protective caps (**see section 4.4, page 28 – Consumables**). Vented protective caps may allow the ingress/leakage of liquids, air and/or pathogens.

9.3.1.2  Connecting a pressure monitoring set to an **unsecured patient's access port** can lead to uncontrolled bleeding, a severe drop in blood pressure and the need for a blood transfusion. Before connecting a pressure monitoring set, make sure that the patient's access port is closed to prevent uncontrolled bleeding and to guarantee safe initial use.


9.3.1.3  Before **transporting or repositioning the patient**, the roller/sliding clamps of the filling and irrigation system should be closed, since otherwise air can get into the patient's circulatory system via the tubing system and this can lead to an embolism. The user must ensure that the entire pressure monitoring set is free of air pockets for the duration of its use.

9.3.1.4 Depending on the configuration of the pressure monitoring set and the therapy of the patient's specific medical condition/disease, the pressure monitoring set can be connected to one or more vascular access ports in different regions of the body (**see section 3.7, page 25 – Type and duration of intended physical contact**).


9.3.2 Connecting a pressure monitoring set with filling system (Configuration A)

9.3.3 Connecting a pressure monitoring set with filling system (Configuration B)


9.3.4 Connecting a pressure monitoring set with filling system (Configuration C)

Instructions for use – English

10 Zero calibration


For correct zero calibration, follow the instructions below.

- 10.2  At the start of a shift, after every **change in the positioning height** of the pressure transducer or patient and after any interruption in signal transmission between pressure transducer and control monitor (e.g. due to a power failure or change of electrical accessory), another zero calibration must be performed to ensure accurate pressure monitoring.

11 Transmission quality check

To check the transmission quality, follow the instructions below.

- 11.2 After connecting the pressure monitoring set to a patient, if it is found that the signals displayed on the patient screen exhibit overshoot or are overdamped, the set-up of the pressure monitoring system must be checked. To do this, follow the instructions in **section 21, page 41 – Troubleshooting and fault rectification**.

12 Disinfection instructions

- 12.1  For disinfection, use **iodine-containing disinfectants** or standard **1-propanol, 2-propanol or ethanol-based disinfectants**. For reasons of user and patient safety, check with the manufacturer about the product-specific contact time of the disinfectant that must be complied with after application in order to adequately kill or inactivate pathogens. Disinfection without observing the required contact time may lead to the spread/transmission of pathogens and thus pose a danger to patients, users and third parties.

- 12.2  Following disinfection, make sure that the applied disinfectant has completely evaporated before using the components. If Luer interfaces are used without complete **evaporation of the disinfectant**, they may break, causing ingress/leakage of liquids, air and/or pathogens.

- 12.3 CODAN pvb Critical Care pressure monitoring sets are developed and manufactured in such a way that when the above-mentioned disinfectants are used and the product-specific contact and drying times and the maximum duration of use of the pressure monitoring set are complied with (**see section 13, page 36 – Duration of use**), any number of disinfections can be carried out without reducing the basic safety of the pressure monitoring set or damaging the affected component(s).

13 Duration of use

- 13.1  For reasons of hygiene, the **pressure monitoring set must be changed** at least every 96 hours and the cumulative duration of use of 30 days maximum should not be exceeded. Moreover, in the event of damage and/or any change in the performance of the individual components, e.g. due to environmental conditions, the pressure monitoring set must be replaced (**see section 17, page 39 – Replacing the pressure monitoring set/sterile accessories**).

- 13.2   **Reusing or reprocessing** the pressure monitoring set and sterile accessories is prohibited under all circumstances, since this can lead to lasting health damage due to infection or allergic/toxic reactions, accompanied by chronic long-term damage to tissues/organs. In addition, compliance with the recommended service life cannot be ensured if systems are reused. Moreover, material properties can be impaired by cleaning processes and as a result, the correct display of measured values, the leak-tightness of adhesively bonded and screw connections and the pressure resistance of the individual components, for example, can no longer be ensured. As a result, neither correct set-up nor proper use of the pressure monitoring set without air pockets and/or leakage of liquids is possible.

13.3 Resterilisation of the pressure monitoring set, accessories and/or individual components is prohibited. If resterilisation is carried out, the basic safety and essential performance characteristics of all the components cannot be ensured, which jeopardises patient safety.

14 Technical specifications and operating conditions

14.1 Electrical safety

14.1.1  The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series have an **IP37 enclosure protection rating**. They are therefore protected against accessing dangerous parts with a tool and against the harmful penetration of solid foreign bodies with a diameter ≥ 2.5 mm as well as water in the event of temporary submersion (30 minutes in water 1 m deep).

14.1.2  The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series have an **IK04 impact rating** and thus enable the measuring function and essential performance characteristics to be maintained after severe mechanical impact stress of up to 0.5 joule.

14.1.3  The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series are electrically isolated and have defibrillation-protected applied parts of the **CF (Cardiac Float) type**, and thus fulfil the strictest requirements for the insulation of accessible parts.

14.2 Conditions for continuous use

14.2.1 For the correct continuous use of the pressure monitoring set with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series, please take note of the technical specifications on page 2 of these Instructions for Use.

14.2.2  It is prohibited to use the pressure monitoring set and/or the accessories **outside the permitted conditions for continuous use** (for example, in the event of transient operating conditions with irregular or severely fluctuating environmental parameters) because it may lead to the loss of the essential performance characteristics and basic safety.

14.3 Electromagnetic compatibility (EMC)

The pressure monitoring sets with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series meet the requirements of the European standards EN 60601-1-2:2015 and EN 60601-2-34:2014 for proof of electromagnetic compatibility (EMC) as well as the additional requirements in accordance with EN 60601-1-12:2015 + A1:2020 for the use of medical electrical equipment in an emergency setting. They are intended for use within a defined and controlled electromagnetic environment. This environment and the tests required as proof of electromagnetic compatibility in accordance with EN 60601-1-2:2015, EN 1789:2020, EN 13718-1:2014+A1:2020, EUROCAE ED-14G:2015/RTCA/DO-160G:2010 + Cor.1:2014 (Section 8.5 Category S, Section 8.8 Category U, Section 20, Section 21 Category M), CISPR 11:2009+A1:2010 and EN 60601-2-34:2014 were successfully performed and are presented in the manufacturer's EMC guidance. The EMC guidance will be provided by the distributor and/or manufacturer upon request.

Temporary impairments due to electromagnetic interference are permissible, but normal operation must be ensured again after a maximum of 10 seconds. A reduction or loss of the essential performance characteristics due to **use outside the defined environmental conditions** and significantly greater electromagnetic interference can lead to an unacceptable risk for patients, users or third parties, for example:

- the loss of essential performance characteristics during use may lead to the patient experiencing circulatory crisis if there is a lack of monitoring;
- an electric shock may lead to the patient experiencing ventricular fibrillation and/or cardiovascular failure;
- the malfunction of other equipment may lead to cardiovascular failure or death of the patient due to the failure of vital equipment;
- if incorrect values are displayed (falsely high/falsely low), this may lead to severe systemic harm to the patient due to the incorrect or omitted administration of medication.

The responsible organisation must ensure that the pressure monitoring set with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series are only used within the environment described in the EMC guidance from the manufacturer CODAN pvb Critical Care.

Using the pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series **together with and/or immediately next to other stacked equipment** must be avoided because this may cause a reduction in the essential performance characteristics. However, if such an application is clinically necessary and unavoidable, all the equipment must be observed during use to make sure that it is functioning as intended.

The **use of electrical accessories from other companies** may result in increased electromagnetic interference emissions and/or decreased electromagnetic immunity of the pressure monitoring set, and lead to a reduction in the essential performance characteristics. Only use original CODAN pvb Critical Care accessories together with pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series (see section 4.3, page 27 – Non-sterile accessories).

Instructions for use – English

- 14.3.4  **Portable radio equipment and other RF communications equipment** (for example, mobile telephones including their antenna cables, external antennas and/or other accessories) must not be used within 30 cm of the electrical accessories and pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series. Failure to observe this distance may lead to a reduction in the performance characteristics of the pressure monitoring set.
- 14.3.5  The pressure monitoring sets can be used during the **use of an electrosurgical unit** (300 W monopolar cutting, 100 W monopolar soft coagulation, 337 kHz). However this may lead to brief malfunctions and a reduction in the essential performance characteristics of the pressure monitoring set. All equipment must be observed during use to ensure that it is functioning as intended.
- 14.3.6  Using a pressure monitoring set **without wearing medical protective gloves and/or outside the environment described in the EMC guidance** from the manufacturer CODAN pvb Critical Care may lead to a reduction or loss of the essential performance characteristics due to harmful electrostatic discharge. The pressure monitoring set should only be used with medical protective gloves and within the environment described in the EMC guidance from the manufacturer CODAN pvb Critical Care in order to prevent harmful electrostatic discharge.

14.4 Compatibility with electrosurgical units

The pressure monitoring sets can be used during use of an **electrosurgical unit**. However this may lead to brief malfunctions and a reduction in the essential performance characteristics of the pressure monitoring set (**see section 14.3, page 37 – Electromagnetic compatibility**).

14.5 Compatibility with magnetic fields (magnetic resonance imaging)

Non-clinical investigations have shown that the pressure monitoring sets with **pressure transducers from series DPT-6000 and DPT-9000 (Xtrans®)** are **MR conditional** according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). A patient can be scanned together with a pressure monitoring set within the MR environment, under the following conditions:

- Static magnetic field of 1.5 Tesla, 3 Tesla and 7 Tesla with
- Maximum spatial gradient magnetic field of 24,600 G/cm (246 T/m)
- Maximum gradient field product of 1,033,000,000 G²/cm (1,033 T²/m)

14.5.1

It was possible to maintain the measuring accuracy of pressure transducers ($\pm 4\%$ / ± 4 mmHg), as required by the European standard EN 60601-2-34, during MR-induced malfunction tests in which the devices were exposed to the static field B₀ at clinical level (200 mT), the time-varying magnetic field gradient (dB/dt) and the RF field (B₁, E) of 1.5 Tesla, 3 Tesla and 7 Tesla MR systems. The devices have passed all of the malfunction tests according to the proprietary test procedures at CODAN pvb Critical Care GmbH.

All components made of plastic or other non-metallic materials in the pressure monitoring sets (e.g. pressure monitoring lines, one-way or multi-way stopcocks and reservoir systems) are MR safe and can be safely placed in the MRI tube. The MR image quality may be compromised by an incorrect proton signal if the plastic components are placed within or around the imaging area.

The DPT-6000 and DPT-9000 (Xtrans®) series pressure transducers and micro drip chamber filling systems must not be placed inside the MRI tube.

Non-clinical investigations have shown that the **systems for filling using a micro drip chamber** are **MR conditional** according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). A patient can be scanned together with a filling system with micro drip chamber within the MR environment, under the following conditions:

14.5.2

- Static magnetic field of 1.5 Tesla, 3 Tesla and 7 Tesla with
- Maximum spatial gradient magnetic field of 24,600 G/cm (246 T/m)
- Maximum gradient field product of 1,033,000,000 G²/cm (1,033 T²/m)

The systems for filling using a micro drip chamber must not be placed inside the MRI tube.

Non-clinical investigations have shown that the following **sterile accessories made of plastic or other non-metallic materials** are **MR safe** according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20):

14.5.3

Systems for filling with a macro drip chamber, systems for filling with a spike, Heidelberger extensions, spiral lines, lines for filling, blood sampling adapters, pressure monitoring lines, arterial blood sampling systems (ABSS), arterial blood sampling systems XL (ABSS^{XL}), needle-free blood sampling systems (NBSS), needle-free blood sampling systems XL (NBSS^{XL}), reservoir systems, reservoir systems XL, sampling extensions, pressure monitoring set additions, CVP measuring limbs, one-way stopcocks, multi-way stopcocks, and sampling ports.

The listed product groups can be safely placed in the tube of the MR system. The MR image quality may be compromised by an incorrect proton signal if the plastic components are placed within or around the imaging area.

14.5.4

CODAN pvb Critical Care devices labelled with the symbol shown opposite are considered **MR unsafe** according to the established terminology of the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). No medical device labelled with this symbol may be used within the MR environment.

15 Storage and transportation

- 15.1  For the correct **storage and transportation** of the pressure monitoring sets with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series, please observe the storage and transportation symbols on page 2 of these Instructions for Use and on the protective or sterile packaging of the device. Storage or transportation of the pressure monitoring set and/or accessories outside of the permitted environmental conditions is prohibited because it may lead to reduced useful life and loss of the essential performance characteristics and basic safety.
- 15.2  Pressure monitoring sets and sterile accessories that are not immediately aseptically prepared, primed and used must not be transported or stored **if the sterile packaging has already been opened**. They must be disposed of because it is no longer possible to guarantee that the sterile condition has been maintained and contamination with pathogens is possible (**see section 19, page 41 – Disposal**).

16 Residual risks

Based on the completed risk assessment in accordance with the EN ISO 14971:2019 + A11:2021 standard, the following residual risks arise with regard to the use of pressure monitoring sets and sterile accessories from the manufacturer CODAN pvb Critical Care: ventricular fibrillation, cardiovascular failure or crisis, drop in blood pressure, drop in cardiac output, drop in oxygen saturation, blood loss, severe systemic damage, chronic long-term damage to tissues and organs, tissue necrosis, thrombosis, (air) embolisms, spread of pathogens, local or systemic infections, vasospasms, sepsis, transient fever, local and allergic reactions, hypersensitivity reactions, pulmonary oedema, malignant tumors, genetic damage, hypervolaemia, malpuncture, vascular injuries, damage to visual organs, burns, or death of the patient.

All residual risks as well as the overall risk fall within the acceptable range after risk mitigation and cannot be reduced any further. The medical/clinical benefit of the pressure monitoring sets and sterile accessories far outweighs the remaining residual risks or the remaining overall risk. It is the responsibility of the attending physician to inform the patient about possible residual risks.

17 Replacing the pressure monitoring set/sterile accessories

During replacement, the safety of the patient must be ensured. Before a pressure monitoring set or accessories can be replaced, the new one(s) must be primed such that they are free of air bubbles. After the patient access port has been closed, the pressure monitoring set or accessories can be disconnected by hand (without tools) and the new one(s) connected. Please follow the instructions in sections 6–13 concerning safe provision, preparation, filling, initial operation and use. The replaced pressure monitoring set and/or accessories must be disposed of (**see section 19, page 41 – Disposal**).

18 Ending use of a pressure monitoring set


18.1 General

- 18.1.1  To end use, follow the hygiene guidelines specified by the responsible organisation and wear **medical protective gloves** to prevent the spread/transmission of pathogens.
- 18.1.2  When use of a pressure monitoring set has to be terminated, for example because there is no longer a medical indication, or in the case of damage or a change in performance, or if the duration of use is exceeded, and the set has come into contact with **potentially infectious substances/materials** (e.g. blood, saliva and/or bodily excretions), it poses a risk and must be disposed of according to the regulations for contaminated/infectious medical waste (biohazard) to prevent the spread/transmission of pathogens (**see section 19, page 41 – Disposal**).

18.2 Uncoupling a pressure monitoring set from the patient access port

- 18.2.1  Disconnecting a pressure monitoring set from an **unsecured patient access port** can lead to uncontrolled bleeding, a severe drop in blood pressure and the need for a blood transfusion. Before uncoupling a pressure monitoring set, make sure that the patient access port is closed to prevent uncontrolled bleeding and to guarantee that it can be disconnected safely.

18.2.2 Uncoupling a pressure monitoring set (Configuration A)

19 Disposal

- 19.1  For disposal, follow the **hygiene guidelines** specified by the responsible organisation and wear **medical protective gloves** to prevent the spread/transmission of pathogens.
- 19.2 Pressure monitoring sets, sterile accessories, consumables and packaging/transportation materials (e.g. protective & sterile packaging, mounting brackets, paper sleeves/labels, rubber bands and transportation strips) that must be disposed of **prior to use** due to damage or because they have passed their expiry date should be disposed of according to the regulations of the responsible organisation.
- 19.3  Pressure monitoring sets, sterile accessories, consumables and packaging/transportation materials (e.g. protective & sterile packaging, mounting brackets, paper sleeves/labels, rubber bands and transportation strips) that have to be disposed of **during or after use** because there is no longer a medical indication, or in the case of damage or a change in performance, or if the duration of use is exceeded, and they have come into contact with **potentially infectious substances/materials** (e.g. blood, saliva and/or bodily excretions), they pose a risk and must be disposed of according to the regulations for contaminated/infectious medical waste (biohazard) to prevent the spread/transmission of pathogens.

20 Customer service

Contact the manufacturer or your distributor for more information.

 CODAN pvb Critical Care GmbH
Römerstraße 18 · 85661 Forstinning, Deutschland
Tel. +49 (0) 81 21 – 98 020
codan@codanpvb.de · www.codancompanies.com

21 Troubleshooting and fault rectification

Phase	Possible complications	Possible causes	Fault rectification Ⓐ / prevention Ⓑ	IFU
Filling	Filling of the pressure monitoring set not possible or only to a limited extent	Insufficient flushing	Ⓐ Check the pressurisation of the saline bag (Configurations A & B) and correct if required Ⓐ Identify twisted or kinked tubes and correct if required Ⓐ Check the stopcock positions and correct if required Ⓐ Check the patency of the filling and Irrigation system and correct if required	9 - 8 7
	Uncontrolled ingress/leakage of saline solution and/or air at Luer interfaces	Loose Luer interfaces	Ⓐ Check the Luer interfaces and tighten by hand if required	6
Pressure monitoring	Signal transmission or zero calibration not possible	Flow channel not patent	Ⓐ Identify twisted or kinked tubes and correct if required Ⓐ Check the stopcock positions and correct if required	- 8
		Incorrect position or incipient occlusion of the patient access port	Ⓐ Check the patency of the patient access port and replace if required Ⓐ Check the pressurisation of the saline bag (Configurations A & B) and correct if required	9 9
		Pressure transducer incorrectly connected to the electrical accessories	Ⓐ Check the connection between the pressure transducer and electrical accessories and correct if required	9
		Electrical accessories damaged or defective	Ⓑ Always handle, clean and store electrical accessories with the utmost care	-
	Pressure curve displayed on the control monitor without pressure monitoring values	Incorrect assembly and/or zero calibration of the pressure monitoring set	Ⓐ Check the positioning height of the pressure transducer and correct if required and carry out another zero calibration	9, 10
		Falsely high/low measured values (overshoot)	Flow channel not patent	Ⓐ Identify twisted or kinked tubes and correct if required Ⓐ Check the stopcock positions and correct if required
	Incorrect position or incipient occlusion of the patient access port		Ⓐ Check the patency of the patient access port and replace if required Ⓐ Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required	9 9
	Air bubbles in the tubing system		Ⓐ Inspect the pressure monitoring set and, if required, remove air bubbles by gently tapping the affected components	8
	Incorrect assembly and/or zero calibration of the pressure monitoring set		Ⓐ Check the positioning height of the pressure transducer and correct if required and carry out another zero calibration	9, 10
	Blood forced into the tubing system	Inadequate liquid column	Ⓐ Check the pressurisation of the saline bag (Configurations A & B) and correct if required Ⓐ Check the stopcock positions and correct if required	9 8
			Lockere Luer-Schnittstellen	Ⓐ Check the Luer interfaces and tighten by hand if required
		Use of unapproved sterile accessories to supplement, extend or change the configuration and/or components of the pressure monitoring set	Ⓑ Only use approved CODAN pvb Critical Care sterile accessories and contact the manufacturer or distributor for more information if required	4

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CODAN Medical AG