

# **IBPM** DPT-6000 (DPT-6003 & DPT-6009) DPT-9000/Xtrans® (DPT-9003 & DPT-9009)

<b>de</b>	<b>Gebrauchsanweisung – Deutsch</b>	<b>12</b>
<b>en</b>	<b>Instructions for use – English</b>	<b>41</b>
<b>bg</b>	<b>Инструкции за употреба – Български</b>	<b>70</b>
<b>cs</b>	<b>Návod k použití – Čeština</b>	<b>100</b>
<b>hr</b>	<b>Upute za upotrebu – Hrvatski</b>	<b>129</b>
<b>hu</b>	<b>Használati utasítás – Magyar</b>	<b>158</b>
<b>lv</b>	<b>Lietošanas instrukcija – Latviešu</b>	<b>187</b>
<b>pl</b>	<b>Instrukcja używania – Polski</b>	<b>216</b>
<b>ro</b>	<b>BInstrucțiuni de utilizare – Română</b>	<b>245</b>
<b>sl</b>	<b>Navodila za uporabo – Slovenščina</b>	<b>275</b>

**DPT-6000** (DPT-6003 & DPT-6009)  
**DPT-9000/Xtrans®** (DPT-9003 & DPT-9009)

Spezifikationen | Specifications | Спецификации | Specifikace | Specifikacije | Műszaki adatok | Specifikācijas | Specyfikacje | Specificații | Specifikacije

	2 V – 15 V ( $\leq 5$ kHz)		0 °C – +40 °C
	-30 mmHg – +360 mmHg		15% – 90%
	4000 mmHg (max. 1 s)		600 hPa – 1060 hPa (-300 m – 4000 m)
	5 $\mu$ V/V/mmHg ( $\pm 1$ %)		0,0001 kWh (max. 0,6 W)
	$\pm 4$ % / $\pm 4$ mmHg		5 kV (<10 s)
			400 $\Omega$ – 10/+20 %
			356 $\Omega$ $\pm 10$ %

EN 60601-1:2006 + Cor.:2010 + A1:2013 + AC:2014 + A1:2013/AC:2014 + A12:2014 + A2:2021  
 EN 60601-1-2:2015  
 EN 60601-1-6:2010 + A1:2015 + A2:2021  
 EN 60601-1-12:2015 + A1:2020  
 EN 60601-2-34:2014

EN ISO 8536-4:2020  
 EN ISO 8536-9:2015  
 EN ISO 8536-10:2015  
 EN 1789:2020

EN 60529:1991 + A2:2013/AC:2019-02  
 EN ISO 80369-1:2018  
 EN ISO 80369-7:2021  
 EN 13718-1:2014 + A1:2020

Medical Device Regulation (EU) 2017/745 (MDR)  
 Directive 2011/65/EU & Directive (EU) 2015/863 (RoHS)  
 Regulation (EC) 1907/2006 (REACH)  
 EN 62366-1:2015 + AC:2015 + A1:2020  
 ASTM F2503-20

**DPT-6000 (DPT-6003 & DPT-6009)**  
**DPT-9000/Xtrans® (DPT-9003 & DPT-9009)**

	A	B	C	D	E	F	X	Y	
I									
II									
III									
IV									
V									
VI									
VII									
VIII									
	DPT-6003 & DPT-9003					DPT-6003/6009 DPT-9003/9009		DPT-6003 & DPT-9003	



































DPT-6000 (DPT-6003 & DPT-6009)  
 DPT-9000/Xtrans® (DPT-9003 & DPT-9009)

# CLC – Classic Configuration

	A	B	C	D	E	F	X	Y	
I									
II									
III									
IV									
V									
VI									
VII									
VIII									
			DPT-6003 & DPT-9003			DPT-6003/6009 DPT-9003/9009		DPT-6003 & DPT-9003	

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









































## 1 Overview and foreword

Read this information carefully before use and keep it for future reference.

These Instructions for Use (IFU) contain important information about the intended use of the pressure monitoring set manufactured by CODAN pvb Critical Care with pressure transducers from the **DPT-6000** (DPT-6003 & DPT-6009), **DPT-9000/Xtrans®** (DPT-9003 & DPT-9009) series and the sterile accessories. The configuration of this pressure monitoring set may differ from the device configurations on pages 3–9 of these Instructions for Use, and therefore the precise guidelines and procedures must be determined by the responsible organisation. The responsible organisation (e.g. the ward manager, responsible doctor or clinical institution) designates the entity responsible for using the pressure monitoring set, the sterile accessories and/or the entire monitoring system in which the CODAN pvb Critical Care pressure monitoring set and/or sterile accessories are combined with other medical devices for invasive blood pressure monitoring.

Any serious incidents that occur during use of, or in connection with, CODAN pvb Critical Care devices and that have led, could have led or might lead to the death or temporary or permanent serious deterioration in a patient's, user's or other person's state of health, or a serious risk to public health must be reported to the manufacturer CODAN pvb Critical Care and to the competent federal authority.

Information labelled with one of the following symbols indicates a hazard,

If this hazard is not avoided, it can result in damage to property or the environment;

Hazard with a low risk level which, if not avoided, can result in minor or moderate injury;

Hazard with a medium risk level which, if not avoided, can result in death or severe injury;

Hazard with a high risk level which, if not avoided, can result in death or severe injury.

Further explanations of the symbols used on the packaging and in the Instructions for Use can be found in the tables on **Page 304**.

## 2 Safety and protective measures

Read and comply with all of the following safety instructions. If these instructions are not observed/complied with, there are significant safety risks for patients, users and third parties.

- 2.1  Comply with all instructions (including safety instructions and warnings) for the intended use of all **(third-party) devices** used in combination with the pressure monitoring set in order to guarantee the safe use of the pressure monitoring set and/or accessories. Prior to use, carefully read all the accompanying information provided by the manufacturers of the (third-party) devices.
- 2.2  **Before performing any therapeutic interventions on a patient** which are to be initiated on the basis of blood pressure values displayed on the control monitor, the plausibility of the measured values must be checked to rule out a device error or malfunction as the cause of measured values exceeding or falling below the set limit values. To do this, check the other patient displays (e.g. electrocardiography, oxygen saturation) and vital functions (e.g. pulse, respiration check). In addition, inspect the entire set-up of the pressure monitoring system at each change of shift (**see section 23, page 69 – Troubleshooting and fault rectification**).
- 2.3  **Exposed** contacts of the electrical accessories must not be touched at the same time as the patient while the accessories are making an electrically conductive connection to a control monitor. **Risk of electric shock!** Wear medical protective gloves when preparing, using and disassembling the electrical accessories and ensure that an electrically conductive connection cannot be established between the patient and the exposed contacts.
- 2.4  **Spilling of liquids**, whatever the quantity, onto electrical contacts can lead to the loss of basic safety, including the essential performance characteristics. Make sure that no liquids can wet the electrical contacts of the accessories and/or pressure monitoring set during preparation, use or disassembly. In the event of liquid contact, replace the pressure monitoring set immediately, clean and dry the electrical accessories and check their function.
- 2.5  When **using a defibrillator**, any contact with the pressure monitoring system is prohibited. **Risk of electric shock!** Avoid any contact with the monitoring system when resuscitation is being performed using a defibrillator.
- 2.6  The use of a pressure monitoring set together with a saline bag and pressure cuff described in the illustrations of these Instructions for Use only applies to pressure monitoring sets with filling system and pressure transducers of types DPT-6003 or DPT-9003 (**see section 7.2, page 50 – Connecting a pressure monitoring set with filling system to a saline bag (DPT-6003 & DPT-9003)**). Pressure monitoring sets with pressure transducers of types DPT-6009 or DPT-9009 should only be used with syringe pumps, as using them with a saline bag and pressure cuff may lead to a saline solution overdose (**see section 7.3, page 51 – Connecting a pressure monitoring set with filling line to an syringe (DPT-6003/6009 & DPT-9003/9009)**). To avoid administering an overdose of saline solution to patients with a strict infusion limit (e.g. neonates and children), the flow rate of the syringe pump must be adjusted in accordance with clinical standards. The male Luer connectors of the pressure transducer types DPT-6009 and DPT-9009 have yellow union nuts to make them clearly identifiable.

# Instructions for Use – English

2.7  In order to maintain the flow rate that is needed, depending on the patient population, to keep the catheter open when using pressure transducer types DPT-6003/6009 & DPT-9003/9009, the occlusion pressure of a syringe pump must be adjusted according to the desired output and clinical standards.

2.8  Before and during emergency use of the device within appropriate emergency vehicles (see 3.8 Intended working environment) during professional primary (emergency) and secondary (intensive care and ambulance) transport to or between clinical healthcare facilities, it is important to ensure that the Luer connectors and the sampling ports of the pressure monitoring set are free of dust, fluff and dirt which can lead to leaks and/or contamination. To prevent soiling or the ingress/leakage of liquids, air and/or pathogens, open-end Luer connectors must be closed with non-vented protective caps.

## 3 Intended use

### 3.1 Intended purpose

The pressure monitoring set is intended for use on a patient requiring continuous monitoring of their blood pressure via one or more vascular access ports (IBPM) and/or requiring blood withdrawals.

### 3.2 Indication

The use of CODAN pressure monitoring sets is always indicated whenever the medical condition requires, at the discretion of the treating physician, continuous monitoring of the blood pressure based on the patient presentation, specific diagnosis, diagnostic findings and/or pathophysiological considerations, among other things.

### 3.3 Clinical benefit

The pressure monitoring set allows for measurement of the blood pressure level in one or more compartments of a patient's circulatory system with subsequent signal transmission to a patient monitor to make possible the visualisation of the measured values on a screen (in the form of digits and/or curves). Furthermore, integrated sampling components or sampling systems (open or closed sampling systems) allow for blood sampling to take place (see section 13, page 57 – Blood sampling).

Therefore, the clinical benefit of the pressure monitoring set is to enable the user to

- continuously monitor and document one or more blood pressure values of a patient in order to establish whether they are above or below defined limit values and/or to detect any irregularities due to (pre-existing) illnesses or blood loss, and to react to them with appropriate treatment;
- make a diagnosis of (pre-existing) illnesses of the cardiovascular system and/or evaluate therapies with drugs that have a cardiovascular effect;
- collect one or more blood sample(s) for laboratory diagnosis, or for evaluation of a patient's metabolic, respiratory and/or haemodynamic situation.

### 3.4 Contraindication

- Allergies to residues from an ethylene oxide [EO] sterilisation process
- Electrolyte imbalances (if necessary, the saline solution must be adjusted prior to use)
- All contraindications to an arterial or venous puncture (e.g. blood coagulation disorders or occlusive diseases)

### 3.5 Intended patient group

Pressure monitoring sets with **DPT-6003** or **DPT-9003** pressure transducers in combination with a filling and irrigation system and a pressure cuff are used on critically ill and/or ventilated patients from 2 years of age (12 kg) who require monitoring of their blood pressure levels in the intravascular system and/or the withdrawal of blood samples.

Pressure monitoring sets with **DPT-6003/6009** or **DPT-9003/9009** pressure transducers in combination with a filling line and a syringe pump are used on critically ill and/or ventilated patients from birth, without any weight restriction, who require monitoring of their blood pressure levels in the intravascular system and/or the withdrawal of blood samples.

### 3.6 Intended user profile

The pressure monitoring set may only be used by healthcare professionals and emergency medical personnel with the appropriate specialised training, who have been trained in accordance with the Instructions for Use for this medical device. The user must not have any mental or physical disability (e.g. colour blindness) that could prevent him/her from using the pressure monitoring set and/or accessories as intended, in accordance with the directions in these Instructions for Use.

### 3.7 Type and duration of intended physical contact

Depending on the configuration of the pressure monitoring set and the therapy for the specific medical condition/disease progression, the pressure monitoring set is connected to one or more vascular access line(s). This/these vascular accessline(s) is/are connected to one or more arterial or venous vessels in one or more area(s) of the patient's body, mainly the extremities and/or neck. The extended exposure time, i.e. the cumulative sum of individual, multiple or repeated contact times, may exceed 24 hours but not 30 days.

### 3.8 Intended working environment

The pressure monitoring set is used within clinical healthcare institutions for monitoring, diagnostics, anaesthesia and intensive care. Additionally, it may be used outside of clinical healthcare facilities within appropriate emergency vehicles during professional primary (emergency) and secondary (intensive care and ambulance) transport, when going to or between clinical healthcare facilities. The envisaged primary and secondary emergency vehicles include


- Ambulances, rescue vehicles and emergency doctor's vehicles;
- Rescue helicopters (piston engine and jet propulsion);
- Rescue aircraft (fixed-wing turbojet or turboprop);
- Rescue ships, lifeboats, water ambulances.

Additionally, the pressure monitoring sets with pressure transducers from the **DPT-6000** (DPT-6003 & DPT-6009) or **DPT-9000/Xtrans®** (DPT-9003 & DPT-9009) series are intended for use within a defined electromagnetic environment (see Section 16, page 64 – Technical Specifications and Operating Conditions).

### 3.9 Intended monitoring system (physical functionality)

The physical functionality of the monitoring system consists of conducting a pressure transmitted via intravascular access ❶ and pressure measurement lines ❷ through a fluid column ❸ to a pressure transducer ❹ and converting it into electrical signals via a Wheatstone full bridge. With the aid of reusable electrical accessories ❺, these signals are transmitted to a control monitor ❻.

## 4 Device description

The medical devices from the manufacturer CODAN pvb Critical Care are developed, manufactured and placed on the market in accordance with the fundamental safety and performance requirements of the Regulation (EU) on Medical Devices 2017/745 (MDR), the Regulation 1907/2006 (REACH), the Directives 2011/65/EU and (EU) 2015/863 (RoHS) and the requirements of the European Standards EN 62366-1:2015 + AC:2015 + A1:2020, and they carry CE marking. The electrical accessories and pressure monitoring sets with pressure transducers from the **DPT-6000** or **DPT-9000/Xtrans®** series meet the requirements of the European Standards EN 60529:1991 + A2:2013/AC:2019-02, EN 60601-1:2006 + Cor.:2010 + A1:2013 + AC:2014 + A1:2013/AC:2014 + A12:2014 + A2:2021, EN 60601-1-2:2015, EN 60601-1-6:2010 + A1:2015 + A2:2021, EN 60601-1-12:2015 + A1:2020, EN 1789:2020, EN 13718-1:2014+A1:2020, EUROCAE ED-14G:2015/RTCA/DO-160G:2010 + Cor.1:2014 (Section 8.5 Category S, Section 8.8 Category U, Section 20, Section 21 Category M) and EN 60601-2-34:2014. The Luer connectors of the sterile accessories and pressure monitoring sets are manufactured in accordance with the requirements of the Standards EN ISO 80369-1:2018 and EN ISO 80369-7:2021. CODAN pvb Critical Care medical devices are delivered in protective/transport packaging validated according to EN ISO 11607-1:2020 and EN ISO 11607-2:2020 to protect the devices from mechanical damage and general soiling during storage and transportation. Moreover, the pressure monitoring sets and sterile accessories are delivered in a single or double sterile barrier system (sterile packaging) validated according to the specified standards to protect against contamination with pathogens after the sterilisation process (in accordance with EN ISO 11135:2014 + A1:2019) and to maintain the sterile condition of the devices.

For a better understanding of the following information about the medical devices of the manufacturer CODAN pvb Critical Care, refer to the configuration overviews on pages 3–9 of these Instructions for Use.

### 4.1 Pressure monitoring set

The pressure monitoring set is a sterile, active and portable class IIb medical device. The essential performance characteristics of the pressure monitoring set in accordance with EN 60601-2-34:2014 include accuracy of the pressure measurements, defibrillation-proof and resistance to electrosurgery interference. The basic safety features of the pressure monitoring set include:

- Leak tightness;
- Pressure resistance;
- Tensile/electrical strength;
- Disinfectant resistance;
- Protection against mechanical stresses;
- Electrical safety.

A fully configured pressure monitoring set comprises

- Components for connecting the pressure monitoring set to a liquid medium for filling and irrigation (**I–III**);
- One or more pressure transducers for converting one or more physical pressures (blood pressure) into one or more electrical signals (**IV**);
- Additional components for pressure transmission (**V**);
- Components for transmitting one or more physical pressures (blood pressure) to one or more pressure transducers and, if applicable, for the collection of blood samples (**VI–VIII**).

The smallest possible configuration for a pressure monitoring set is a single pressure transducer (**IV**) of the **DPT-6000** (DPT-6003 & DPT-6009) or **DPT-9000/Xtrans®** series (DPT-9003 & DPT-9009). All other sections (**I–III** & **V–VIII**) are optional and must be supplemented with the sterile accessories from the manufacturer CODAN pvb Critical Care (see section 4.2, page 46 – **Sterile accessories**). The following table lists the components of a pressure monitoring set and their main functional and operating features, based on the configuration overviews on pages 3–9 of these Instructions for Use.

Pressure monitoring lines are required to connect different components and for pressure transmission. The pressure monitoring lines can be made from PVC (polyvinyl chloride), PE (polyethylene), or PUR (polyurethane) with different Luer connectors, lengths and (internal/external) diameters and they can be installed at various positions on the pressure monitoring set. Due to their versatility, the pressure monitoring lines are not explained in further detail in the following table. More information about the pressure monitoring lines from the manufacturer CODAN pvb Critical Care will be provided by the distributor and/or manufacturer upon request.

Overview (p. 3–9)			Component				
Blood sampling	Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
All	-	-	Protective/transport packaging	-	<ul style="list-style-type: none"> <li>• Protects the medical device from mechanical damage and general soiling during storage and transportation</li> <li>• Identifiability and traceability</li> </ul>	1	A sterile barrier system with external protective packaging is labelled with the symbol 20

# Instructions for Use – English

Overview (p. 3–9)			Component				
Blood sampling	Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
All	-	-	Sterile packaging (sterile barrier system)	-	<ul style="list-style-type: none"> <li>Protection against the ingress of pathogens</li> <li>Maintaining the sterile condition</li> <li>Identifiability and traceability</li> <li>Easy and residue-free opening</li> </ul>	1–2	<ul style="list-style-type: none"> <li>A single sterile barrier system is labelled with the symbol 21</li> <li>A double sterile barrier system is labelled with the symbol 22</li> </ul>
All	A	I–III	AeroGuard® (filling and irrigation system with filter system)	-	<ul style="list-style-type: none"> <li>Connecting the pressure monitoring set to a saline bag for filling and rinsing the pressure monitoring set with saline solution</li> <li>Particle filtration and air separation</li> </ul>	1	Only in combination with pressure transducer types DPT-6003 or DPT-9003
All	A	I	Drip chamber macro (vented or non-vented)	AeroGuard®	<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (20 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>	1	-
			Drip chamber micro (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (60 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>		
			Spike (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> </ul>		
		II	Roller clamp		<ul style="list-style-type: none"> <li>Flow control (opening and closing the flow channel)</li> </ul>	1	-
		II	Filter system		<ul style="list-style-type: none"> <li>Particle filtration (1.2 µm)</li> <li>Elimination of air pockets</li> </ul>	1	-
All	B–D	I–III	Filling and irrigation system	-	Connecting the pressure monitoring set to a saline bag for filling and rinsing the pressure monitoring set with saline solution	1	Only in combination with pressure transducer types DPT-6003 or DPT-9003
All	B–D	I	Drip chamber macro (vented or non-vented)	Filling and irrigation system	<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (20 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>	1	-
			Drip chamber micro (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (60 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>		
			Spike (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> </ul>		
	B–D	1	Roller clamp		<ul style="list-style-type: none"> <li>Flow control (opening and closing the flow channel)</li> </ul>	1	-
	D	III	Sliding clamp		<ul style="list-style-type: none"> <li>Flow control (opening and closing the flow channel)</li> </ul>	1–4	Only in combination with pressure monitoring sets with two or more filling cannulas
			Free filling outflow with a male Luer connector		<ul style="list-style-type: none"> <li>Addition of a pressure transducer without a filling system or filling line (X/IV – CLC &amp; SWAN)</li> </ul>	1–3	-
All	E	I–III	ADD filling and irrigation system (basic)	-	<ul style="list-style-type: none"> <li>Connecting the pressure monitoring set to a saline bag for filling and rinsing the pressure monitoring set with saline solution</li> <li>Connecting the pressure monitoring set to an ADD filling system (extension) (Y/II – CLC &amp; SWAN)</li> </ul>	1	Only in combination with pressure transducer types DPT-6003 or DPT-9003
All	E	I	Drip chamber macro (vented or non-vented)	ADD filling and irrigation system (basic)	<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (20 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>	1	-
			Drip chamber micro (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> <li>Setting the liquid level</li> <li>Flow control (60 drops per millilitre)</li> <li>Particle filtration (15 µm)</li> </ul>		
			Spike (vented or non-vented)		<ul style="list-style-type: none"> <li>Connection to a saline bag (vented or non-vented)</li> </ul>		
		I	Roller clamp		<ul style="list-style-type: none"> <li>Flow control (opening and closing the flow channel)</li> </ul>	1	-
		III	4-way stopcock		<ul style="list-style-type: none"> <li>Connection to an ADD addition line (Y/II – CLC &amp; SWAN)</li> <li>Setting of four possible flow directions (360°)</li> </ul>	1	-
All	F	II–III	Filling line with female Luer connector	-	Connecting the pressure monitoring set to a disposable syringe for filling and rinsing the pressure monitoring set with saline solution	1–2	Only in combination with pressure transducer types DPT-6009 or DPT-9009
-, CLC, SWAN	Y	II–III	ADD filling and irrigation system (addition)	-	Connecting the pressure monitoring set to an ADD filling system (basic) (E/III) or an ADD filling system (addition) (Y/II – CLC & SWAN) for filling and rinsing the pressure monitoring set with saline solution	1	Only in combination with pressure transducer types DPT-6003 or DPT-9003
-, CLC, SWAN	Y	II	ADD addition line with a male Luer connector	ADD filling and irrigation system (addition)	Connection to a 4-way stopcock of an ADD filling system (basic) (E/III – CLC & SWAN) or an ADD filling system (addition) (Y/III – CLC & SWAN)	1	-
		III	4-way stopcock		<ul style="list-style-type: none"> <li>Connection to an ADD addition line of an ADD filling system (addition) (Y/III – CLC &amp; SWAN)</li> <li>Setting of four possible flow directions (360°)</li> </ul>	1	-

Overview (p. 3–9)			Component				
Blood sampling	Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
-, CLC, SWAN	A–F, X,Y	IV	Pressure transducer types DPT-6003 or DPT-9003 (Xtrans®)	-	<ul style="list-style-type: none"> <li>• Conversion of a physical pressure (blood pressure) into an electrical signal</li> <li>• Rapid flushing and continuous flushing</li> </ul>	1–4	-
NBSS, NBSS <sup>XL</sup> , ABSS, ABSS <sup>XL</sup>	A–F						
-, CLC, SWAN	A–F, X,Y	IV	Flush system	-	<ul style="list-style-type: none"> <li>• Rapid flushing for filling and the removal of air pockets and/or blood residue</li> <li>• Continuous flushing (3 ml/h at 300 mmHg) to prevent the patient access port from becoming clogged</li> </ul>	1	-
NBSS, NBSS <sup>XL</sup> , ABSS, ABSS <sup>XL</sup>	A–F						
-, CLC, SWAN	A–F, X,Y	IV	Electrical contacts	DPT-6003, DPT-9003	<ul style="list-style-type: none"> <li>• Connecting the pressure transducer to an electrical accessory for signal transmission</li> </ul>	1	-
NBSS, NBSS <sup>XL</sup> , ABSS, ABSS <sup>XL</sup>	A–F						
-, CLC, SWAN	A–F, X,Y	IV	3-way zeroing stopcock	-	<ul style="list-style-type: none"> <li>• Zero calibration of the pressure transducer</li> <li>• Setting of three possible flow directions (180°)</li> </ul>	1	-
NBSS, NBSS <sup>XL</sup> , ABSS, ABSS <sup>XL</sup>	A–F						
All	F	IV	Pressure transducer types DPT-6009 or DPT-9009 (Xtrans®)	-	<ul style="list-style-type: none"> <li>• Conversion of a physical pressure (blood pressure) into an electrical signal</li> <li>• Rapid flushing and continuous flushing</li> </ul>	1–2	-
All	F	IV	Flush system	DPT-6009, DPT-9009	<ul style="list-style-type: none"> <li>• Rapid flushing for filling and the removal of air pockets and/or blood residue</li> <li>• Continuous flushing (30 ml/h at 300 mmHg) to prevent the patient access port from becoming clogged</li> </ul>	1	-
			Electrical contacts		<ul style="list-style-type: none"> <li>• Connecting the pressure transducer to an electrical accessory for signal transmission</li> </ul>	1	-
			3-way zeroing stopcock		<ul style="list-style-type: none"> <li>• Zero calibration of the pressure transducer</li> <li>• Setting of three possible flow directions (180°)</li> </ul>	1	-
All	B	V	Temporary pressure monitoring	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Alternating between two pressure channels</li> </ul>	1	Only in combination with pressure monitoring sets with a single pressure transducer
All	B	V	2-way stopcock	Temporary pressure monitoring	Setting of two possible flow directions (90°)	1	<ul style="list-style-type: none"> <li>• 2-way stopcock with or without a connected or enclosed pressure monitoring line</li> <li>• 3-way or 4-way stopcock with a connected or enclosed pressure monitoring line</li> </ul>
			3-way stopcock		Setting of three possible flow directions (180°)		
			4-way stopcock		Setting of four possible flow directions (360°)		
All	C	V	S.E.R.O. (Selectable Eliminator of Resonance Overshoot) damping stopcock	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Switching an attenuator on and off to correct falsely high/low measured values (overshoot)</li> <li>• Setting of four possible flow directions (360°)</li> </ul>	1	-
CLC	A–Y	VI–VII	CLC (Classic Configuration)	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Blood sampling</li> </ul>	1	-
CLC	A–Y	VII	3-way stopcock	CLC	<ul style="list-style-type: none"> <li>• Blood sampling</li> <li>• Setting of three possible flow directions (180°)</li> </ul>	1–3	ArterioGuard® only in combination with a (pulmonary)-arterial set-up of a pressure monitoring set
			4-way stopcock		<ul style="list-style-type: none"> <li>• Blood sampling</li> <li>• Setting of four possible flow directions (360°)</li> </ul>		
			ArterioGuard® (4-way stopcock)		<ul style="list-style-type: none"> <li>• Blood sampling</li> <li>• Check valve protects against intra-arterial injections initiated by mistake</li> <li>• Setting of four possible flow directions (360°)</li> </ul>		
SWAN	A–Y	VI–VII	SWAN Switch (Swabable Adaptor Needle-free)	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Hygienically closed blood sampling</li> <li>• Setting of four possible flow directions (360°)</li> </ul>	1–3	-
NBSS	A–F	VI–VII	NBSS (Needle-free Blood Sampling System)	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Hygienically closed blood sampling</li> <li>• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood</li> </ul>	1	Only in combination with a (pulmonary)-arterial set-up of a pressure monitoring set
NBSS	A–F	VI	Reservoir (2 ml)	NBSS	<ul style="list-style-type: none"> <li>• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood (2 ml)</li> <li>• Protection against the ingress of pathogens</li> </ul>	1	-
		VII	SWAN Switch (Swabable Adaptor Needle-free)		<ul style="list-style-type: none"> <li>• Hygienically closed blood sampling</li> <li>• Setting of four possible flow directions (360°)</li> </ul>	1–3	-
NBSS <sup>XL</sup>	A–F	VI–VII	NBSS <sup>XL</sup> (Needle-free Blood Sampling System)	-	<ul style="list-style-type: none"> <li>• Transmission of the physical pressure (blood pressure) to a pressure transducer</li> <li>• Hygienically closed blood sampling</li> <li>• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood</li> </ul>	1	Only in combination with a (pulmonary)-arterial set-up of a pressure monitoring set

# Instructions for Use – English

Overview (p. 3–9)			Component				
Blood sampling	Column	Row	Designation	Part of	Main functional & operating features	Quantity	Note
NBSS <sup>XL</sup>	A–F	VI	1-way stopcock	NBSS <sup>XL</sup>	• Flow control (opening and closing the flow channel)	1	-
			Reservoir (6 ml)		• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood (6 ml) • Protection against the ingress of pathogens		
		VII	SWAN Switch (Swabable Adaptor Needle-free)		• Hygienically closed blood sampling • Setting of four possible flow directions (360°)	1–3	-
ABSS	A–F	VI–VII	ABSS (Arterial Blood Sampling System)	-	• Transmission of the physical pressure (blood pressure) to a pressure transducer • Hygienically closed blood sampling • Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood	1	Only in combination with a (pulmonary)-arterial set-up of a pressure monitoring set
ABSS	A–F	VI	Reservoir (2 ml)	ABSS	• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood (2 ml) • Protection against the ingress of pathogens	1	-
			1-way stopcock		• Flow control (opening and closing the flow channel)	1–2	-
		VII	Sampling port		• Hygienically closed blood sampling	1–2	The blood sampling adapter HG-01 is required for blood sampling
ABSS <sup>XL</sup>	A–F	VI–VII	ABSS <sup>XL</sup> (Arterial Blood Sampling System)	-	• Transmission of the physical pressure (blood pressure) to a pressure transducer • Hygienically closed blood sampling • Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood	1	Only in combination with a (pulmonary)-arterial set-up of a pressure monitoring set
ABSS <sup>XL</sup>	A–F	VI	1-way stopcock	ABSS <sup>XL</sup>	• Flow control (opening and closing the flow channel)	1	-
			Reservoir (6 ml)		• Drawing, short-term storage (duration of a blood sampling process) and return of mixed blood (6 ml) • Protection against the ingress of pathogens	1	-
		VII	Sampling port		• Hygienically closed blood sampling	1–2	The blood sampling adapter HG-01 is required for blood sampling
			1-way stopcock		• Flow control (opening and closing the flow channel)	1–2	-
All	A–Y	VIII	Patient connection with a male Luer connector	-	Connecting the pressure monitoring set to a patient access port	1–4	-

## 4.2 Sterile accessories

Sterile accessories (class IIIa medical devices) are individual, sterile packaged components of a pressure monitoring set from the manufacturer CODAN pvb Critical Care, not including individual pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series and consumables. The sterile accessories are used to supplement, extend or change the configuration and/or components of a pressure monitoring set. These include filling systems, blood sampling systems (only in combination with pressure transducer types DPT-6003, DPT-6009, DPT-9003 and DPT-9009), pressure monitoring lines and one-way or multi-way stopcocks with or without connected pressure monitoring lines.

**Additions, extensions or changes to the configuration and/or components of a pressure monitoring set** are always the responsibility of the responsible organisation and must only be implemented using the sterile accessories from the manufacturer CODAN pvb Critical Care. Additions, extensions or changes to the configuration and/or components using third-party devices are prohibited, because the basic safety and essential performance characteristics of the pressure monitoring set cannot then be guaranteed. It is recommended that the components between the pressure transducer and patient access port should not exceed 240 cm in length, as this can possibly affect the quality of the signal transmission. Moreover, additions, extensions or changes to the configuration and/or components of the pressure monitoring set between the sampling port and patient access port may mean that sufficient volume can no longer be drawn using a reservoir (NBSS, NBSS<sup>XL</sup>, ABSS & ABSS<sup>XL</sup> configurations) for undiluted blood samples to be taken at the sampling port.

A list of all the available sterile accessories and other information will be provided by the distributor and/or manufacturer upon request.

## 4.3 Non-sterile accessories

**Reusable mounting accessories**, comprising holders, mounting plates and modular holder systems (mounting heads, claws and poles), are non-sterile class I medical devices and they are used for assembling pressure monitoring sets and electrical accessories on round poles or rail systems.

**Single-use mounting accessories**, comprising velcro straps and fasteners, are non-sterile class I medical devices used for assembling pressure monitoring sets, reservoirs and electrical accessories on the patient's arm or for securing pressure monitoring lines.

**Reusable electrical accessories** are non-sterile class I medical devices that enable the transmission of signals between one or more pressure transducers and a control monitor. These include contact boards and monitor connectors for **DPT-9000/Xtrans®** and interface cables and cable management for **DPT-6000**.

A list of all the available non-sterile accessories and other information will be provided by the distributor and/or manufacturer upon request or can be downloaded from: <https://www.codancompanies.com/downloadcenter/>

## 4.4 Consumables

Consumables are required for the intended use of a pressure monitoring set on a patient within a monitoring system and may have to be used up, disposed of or replaced before or during use. Illustrations in these Instructions for Use labelled with the following symbols inform the user about

➔ a used material that must be disposed of according to regulations (see section 21, page 68 – Disposal);

→ a new consumable that is required for the intended use of a pressure monitoring set and must be provided.

Since the configuration of this pressure monitoring set may differ from the device configurations on pages 3–9 of these Instructions for Use, the responsible organisation must identify all the consumables for the pressure monitoring set prior to use and provide a sufficient quantity of these consumables to the user (see section 5, page 48 – **Combination with medical devices from other companies**).









#### 4.4.1 Protective caps

Non-vented protective caps are closures with one or more Luer connectors. They are attached to the free Luer ports of a pressure monitoring set and/or sterile accessory to protect against the ingress/leakage of liquids, air and/or pathogens before and/or during use.

Please note that upon delivery, **vented protective caps** may be attached to the Luer ports of a pressure monitoring set and/or sterile accessory. These caps must be replaced with **non-vented protective caps** prior to use to prevent the ingress/leakage of liquids, air and/or pathogens.

The list below show all the protective caps that may be found in a CODAN pvb Critical Care pressure monitoring set when it is delivered.

Protective caps from CODAN pvb Critical Care		Colour coding	
Illustration	Luer type	Vented	Non-vented
	Male	Red, blue	Neutral
	Male	Red	Neutral
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–
	Female	Neutral	–

#### 4.4.2 HG-01 for ABSS and ABSS<sup>XL</sup> configurations

The individually and sterile packaged blood sampling adapter HG-01 is a class IIa medical device, which enables blood sampling via the sampling port of the ABSS and ABSS<sup>XL</sup> (see section 13, page 57 – **Blood sampling**).

The **sampling port** of the ABSS and ABSS<sup>XL</sup> must not be punctured with needles or cannulas, as this may cause damage to the plastic membrane and lead to the ingress/leakage of liquids and/or air.  
To puncture the membrane, only use the blood sampling adapter HG-01 (REF: 74.6934) together with disposable syringes and/or equivalent sampling containers with a Luer connector.

#### 4.4.3 Consumables of the monitoring system

The list below shows all the consumables required for the intended use of the pressure monitoring set within a monitoring system:

- Cannula(s) and/or catheter(s) (vascular access);
- Bag with isotonic (0.9%) saline solution (for pressure monitoring sets with filling and irrigation system for use with a pressure cuff);
- Disposable syringe(s) for pressure monitoring sets with filling line for use with a syringe pump).

#### 4.4.4 Other consumables

The list below shows all the other consumables required for the intended use of the pressure monitoring set:

- Disposable syringe(s) and/or equivalent sampling containers for volume shift (see section 13, page 57 – **Blood sampling**);
- Disinfection wipes/swabs and/or disinfectant (see section 14, page 63 – **Disinfection instructions**);
- Medical protective gloves.

#### 4.5 Colour coding

Take note of the colour coding of the pressure monitoring lines and one-way or multi-way stopcocks. Red stopcocks and pressure monitoring lines are intended for arterial, blue for venous, yellow for pulmonary-arterial and green for left atrial applications. Neutral stopcocks and pressure monitoring lines may be used individually for these applications. Components with neutral colour coding must be suitably labelled prior to use to prevent mix-ups or incorrect use, for example using coding pins or colour-coded stickers.

Colour coding	Application/catheter position	Examples of puncture sites
Red	Arterial	Radial artery, ulnar artery, brachial artery, femoral artery, dorsalis pedis artery
Blue	Venous	Internal jugular vein, external jugular vein, cephalic vein, basilic vein, subclavian vein
Yellow	Pulmonary-arterial	Internal jugular vein
Green	Left atrial	Right femoral vein

# Instructions for Use – English

## 5 Combination with medical devices from other companies

Defined below are the minimum requirements for medical devices from other companies that must be met to ensure the safe and intended use of a pressure monitoring set with pressure transducers of the **DPT-6000** (DPT-6003 & DPT-6009) or **DPT-9000/Xtrans®** (DPT-9003 & DPT-9009) series and/or CODAN pvb Critical Care accessories within a monitoring system. Moreover, it must be proven that the medical devices from other companies were manufactured and authorised in accordance with all applicable device-specific standards and that they enable the pressure monitoring set and/or accessories to be used in accordance with the operating conditions specified in these Instructions for Use (**see section 16, page 64 – Technical specifications and operating conditions**).

Any deviation from the listed requirements is always the responsibility of the responsible organisation. The responsible organisation must ensure that combining a pressure monitoring set and/or accessories from the manufacturer CODAN pvb Critical Care GmbH with one or more medical devices from another company, including all connections/interfaces, is safe and does not cause any reciprocal interference, which could lead to reduction or loss of the basic safety and/or essential performance characteristics of one or more devices within the monitoring system.

### 5.1 General requirements

Medical devices from other companies must be developed, manufactured and placed on the market in compliance with the fundamental safety and performance requirements of the Medical Devices Regulation (EU) 2017/745 (MDR) or the Medical Devices Directive 93/42/EEC (MDD), and they must also be CE marked.

### 5.2 Medical electrical equipment

Any user who combines medical electrical equipment from another company with a pressure monitoring set from CODAN pvb Critical Care, with or without the aid of accessories, is the System Configurator and thereby responsible for ensuring that the monitoring system verifiably meets the specified requirements. Medical electrical equipment that is used in an emergency environment and combined with a pressure monitoring set, with or without the aid of accessories, must also verifiably meet the current requirements of the European standard EN 60601-1-12.

#### 5.2.1 Control monitors

Control monitors intended to be combined with CODAN pvb Critical Care electrical accessories and pressure monitoring sets within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standards EN 60601-1, EN 60601-1-2, EN 60601-1-6, EN 60601-1-12 (applicable for applications in an emergency environment) and EN 60601-2-34.

#### 5.2.2 Syringe/infusion pumps

Syringe/infusion pumps intended to be combined with CODAN pvb Critical Care pressure monitoring sets with filling line and a disposable syringe within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standards EN 60601-1, EN 60601-1-2, EN 60601-1-6, EN 60601-1-12 (applicable for applications in an emergency environment) and EN 60601-2-24.

### 5.3 Sterile and non-sterile accessories

Sterile and non-sterile accessories are used for **additions, extensions or changes to the configuration and/or components of the pressure monitoring set** and for signal transmission or attachment (**see section 4, page 43 – Device description**). Only CODAN pvb Critical Care accessories are authorised for use with a pressure monitoring set with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series. In addition, the accessories should only be used in combination with the CODAN pvb Critical Care pressure monitoring set. Non-compliance can lead to reduction or loss of the basic safety and/or essential performance characteristics.

### 5.4 Consumables

Consumables are required for the intended use of a pressure monitoring set on a patient within a monitoring system and may have to be used up, disposed of and replaced before or during use (**see section 4.4, page 46 – Consumables**).

#### 5.4.1 Consumables with Luer connectors

The Luer connectors of consumables (protective caps, catheters, cannulas, disposable syringes, blood sampling adapters or equivalent blood sampling containers) must be verifiably manufactured and authorised in accordance with the current requirements of the standards EN ISO 80369-1 and EN ISO 80369-7 with a 6% (Luer) taper for intravascular or hypodermic applications.

**Consumables with Luer connectors** that are not manufactured and authorised in accordance with the above-mentioned standards may allow the ingress/leakage of liquids, air and/or pathogens when filling and/or using a pressure monitoring set. Before using the consumables, verify that they are manufactured and authorised in accordance with the specified standards.

##### 5.4.1.1 Catheters and cannulas

Catheters and cannulas intended to be combined with CODAN pvb Critical Care pressure monitoring sets within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standards EN ISO 10555-1 (catheters) or EN ISO 7864 (cannulas).

##### 5.4.1.2 Disposable syringes

Disposable syringes intended to be combined with CODAN pvb Critical Care pressure monitoring sets for volume shifts must verifiably meet the current requirements of the European standard EN ISO 7886-1. Disposable syringes intended to be combined with CODAN pvb Critical Care pressure monitoring sets with filling line for use on a syringe pump for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard EN ISO 7886-2.

#### 5.4.2 Other consumables

##### 5.4.2.1 Bags containing isotonic (0.9%) saline solution

Saline bags intended to be combined with CODAN pvb Critical Care pressure monitoring sets with filling and irrigation system within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard EN ISO 15747.

**5.4.2.2 Medical protective gloves**

Medical protective gloves intended to be worn before, during and after use of CODAN pvb Critical Care accessories and/or pressure monitoring sets must verifiably meet the current requirements of the European standard EN 455-1.

**5.5 Other (third party) devices**

**5.5.1 Pressure cuffs**

The manometers of pressure cuffs intended to be combined with CODAN pvb Critical Care pressure monitoring sets with filling and irrigation system and a saline bag within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard EN ISO 81060- 1. Moreover, the manometers must have a measurement range of 0 mmHg–300 mmHg (or greater) and a measurement accuracy of at least  $\pm 3$  mmHg/ $\pm 2\%$ .

**5.5.2 Mounting poles & rails**

Mounting rails intended to be combined with CODAN pvb Critical Care reusable mounting accessories within a monitoring system for invasive blood pressure monitoring must verifiably meet the current requirements of the European standard DIN EN ISO 19054 and have a maximum width of 25 mm and a maximum thickness of 10 mm. The tubing of the mounting poles or infusion stands must have a diameter of between 15 mm and 35 mm. The rectangular profile of the mounting poles or infusion stands must measure 30 mm x 30 mm.


**6 Provision**

To prepare correctly, follow the instructions below.

**6.1 Unpacking**

6.1.1  Pressure monitoring sets and sterile accessories with **damaged protective and/or sterile packaging (sterile barrier system)** must not be used. They must be disposed of because the basic safety and essential performance characteristics and sterility of the contents can no longer be guaranteed and contamination with pathogens is possible (**see section 21, page 68 – Disposal**). Before unpacking, make sure that the protective and sterile packaging is undamaged and inspect the integrity of the sealed seams.

6.1.2  Pressure monitoring sets and sterile accessories must not be used after the **expiry date has passed**. They must be disposed of because the basic safety and performance characteristics and sterility of the contents can no longer be guaranteed and contamination with pathogens is possible (**see section 21, page 68 – Disposal**). Before unpacking, make sure that the expiry date stated on the protective and/or sterile packaging of the device has not passed.

6.1.4  Using the monitoring system and its components without observing the **hygiene guidelines** specified by the responsible organisation may lead to the spread/transmission of pathogens and thus pose a danger to patients, users and third parties. The entire monitoring system and its components should only be used in accordance with the specified hygiene guidelines and wearing **medical protective gloves**.

6.1.5  **After the sterile packaging is opened**, the pressure monitoring set and/or sterile accessories must immediately be prepared aseptically, filled and used in order to minimise the time that the device is exposed to potential pathogens. The responsible organisation must determine the appropriate length of time that may pass after the device is removed from the sterile packaging and before it is used.


6.1.6  If the packaging is opened by mistake, the pressure monitoring sets and/or sterile accessories must be aseptically processed, primed and used immediately in accordance with the specifications of the responsible organisation. Otherwise they should be disposed of because the sterility of the devices cannot be guaranteed and contamination with pathogens is possible (**see section 21, page 68 – Disposal**).

6.1.7  Pressure monitoring sets and accessories with **visible defects and/or mechanical damage** must not be used. They must be disposed of because the essential performance characteristics and basic safety can no longer be guaranteed (**see section 21, page 68 – Disposal**). After unpacking, make sure that the devices are undamaged.


6.1.8  **Loose Luer connections** may allow the ingress/leakage of liquids, air and/or pathogens during filling and/or use. Before connecting the pressure monitoring set to the patient, check that all the Luer connections are tight and, if necessary, tighten them by hand (without tools). It is prohibited to use tools/instruments to tighten the Luer connections because they may break if excessive force is applied, resulting in an open flow channel.

# Instructions for Use – English


## 6.2 Providing of the monitoring system

## 6.3 Providing of addition sets

## 6.4 Providing of ADD sets

## 7 Preparation


To prepare correctly, follow the instructions below.

7.1.1  The pressure monitoring set is not approved for the **administration of medicinal products and/or other liquid media (e.g. contrast media)** because these may damage the components of the pressure monitoring set, meaning that the basic safety and essential performance characteristics are not guaranteed. The use of 0.9% saline solution (possibly heparinised) to fulfil the intended use is permitted.

7.1.2  If the **reservoir plunger** is in an open position upon delivery, the pressure monitoring set must not be used and must be disposed of because the sterility of the device cannot be guaranteed (**see section 21, page 68 – Disposal**). Before filling the pressure monitoring set, check that the reservoir is in a closed position.

## 7.2 Connecting a pressure monitoring set with filling system to a saline bag (DPT-6003 & DPT-9003)


7.2.1 Please note that the instructions in section 7.2 are only applicable for pressure monitoring sets with filling and irrigation system and pressure transducer types DPT-6003 or DPT-9003. Pressure monitoring sets with pressure transducer types DPT-6009 or DPT-9009 should only be used with syringe pumps (**see section 7.3, page 51 – Connecting a pressure monitoring set with filling line to a syringe (DPT-6003/6009 & DPT-9003/9009)**).

7.2.3 When filling the pressure monitoring set, make sure that no more than 50 mmHg pressure is applied to the connected saline bag because otherwise tiny air bubbles can form inside the pressure monitoring set, which may negatively affect the pressure transmission and cause an air embolism in the patient's vascular system. It is always recommended to gravity fill the pressure monitoring set to prevent the formation of tiny air bubbles.

**7.3 Connecting a pressure monitoring set with filling line to a syringe (DPT-6003/6009 & DPT-9003/9009)**

7.3.1 Please note that the instructions in section 7.3 are only applicable for pressure monitoring sets with filling line and for pressure transducer of the types DPT-6003/6009 & DPT-9003/9009.


**8 Filling**

Follow the instructions below to correctly prime the pressure monitoring set.


**Air pockets** in the pressure monitoring set negatively affect the signal transmission and can lead to an air embolism inside the patient's vascular system. All the components of a pressure monitoring set must be completely primed with isotonic saline solution and free of air bubbles before the pressure monitoring set is connected to a patient access port. Any remaining air pockets in the pressure monitoring set must be loosened using the rapid flushing function of the pressure monitoring set and/or by carefully tapping the affected components and then flushing the air pockets out of the pressure monitoring set. Moreover, throughout the entire duration of use of the pressure monitoring set, it must be ensured that all the components are free of air pockets. To do this, follow the illustrations and instructions in section 8.

**8.1 Filling the outlets of an addition set**


**8.2 Filling a pressure transducer (DPT-6003/6009 & DPT-9003/9009)**


# Instructions for Use – English


## 8.3 Filling a 2-way stopcock


## 8.4 Filling an S.E.R.O. (Selectable Eliminator of Resonance Overshoot) damping stopcock


## 8.5 Filling a reservoir (2 ml)


## 8.6 Filling an XL reservoir (6 ml)

## 8.7 Filling a SWAN switch

### 8.7.1


### 8.7.2


Do not screw any **protective caps** onto the side port of the **SWAN switch 4-way stopcock** because this can damage the valve and may therefore lead to the ingress/leakage of liquids, air and/or pathogens.


**8.8 Filling the sampling port of an ABSS**


**8.9 Filling the sampling port of an ABSS<sup>XL</sup>**

**8.10 Filling an ArterioGuard®**

**8.11 Filling a 4-way stopcock**










# Instructions for Use – English

## 9 Initial use







Follow the instructions below for the correct initial use of the pressure monitoring set.

### 9.1 Assembly

#### 9.1.1 Assembly of a pressure monitoring set with DPT-6000 (DPT-6003 & DPT-6009) pressure transducer


#### 9.1.2 Assembly of a pressure monitoring set with DPT-9000/Xtrans® (DPT-9003 & DPT-9009) pressure transducer


9.1.2.3

**9.2 Positioning**

9.2.1

9.2.2

Make sure that there is always enough liquid in the **connected saline bag** because otherwise continuous flushing cannot be guaranteed and this may lead to the patient access port becoming clogged. If a new saline bag is needed, the roller/sliding clamps of the pressure monitoring set should be closed before removing the old bag. To connect a new saline bag, follow the illustrations in **section 7.2, page 50 – Connecting a pressure monitoring set with filling and irrigation system to a saline bag (DPT-6003 & DPT-9003)**.

9.2.3

Every hour and after each rapid flushing, check that a pressure of **300 mmHg** is applied to the connected saline bag, as otherwise no continuous flushing is ensured and this can cause the patient access port to become clogged.

**9.3 Patient connection**

9.3.1


Before using the pressure monitoring set on a patient, make sure that all the **vented protective caps** have been disposed of and replaced with non-vented protective caps (**see section 4.4, page 46 – Consumables**). Vented protective caps may allow the ingress/leakage of liquids, air and/or pathogens.

9.3.2

Connecting a pressure monitoring set to an **unsecured patient access port** can lead to uncontrolled bleeding, a severe drop in blood pressure and the need for a blood transfusion. Before connecting a pressure monitoring set, make sure that the patient access port is closed to prevent uncontrolled bleeding and to guarantee safe initial use.

9.3.3

9.3.4

Before **transporting or repositioning the patient**, the roller/sliding clamps of the filling and flush system should be closed, since otherwise air can get into the patient's circulatory system via the tubing system and this can lead to an embolism. The user must ensure that the entire pressure monitoring set is free of air pockets for the duration of its use.


9.3.5

Depending on the configuration of the pressure monitoring set and the therapy of the patient's individual medical condition/disease, the pressure monitoring set can be connected to one or more vascular access ports in different regions of the body (**see section 3.7, page 42 – Type and duration of intended physical contact**).

# Instructions for Use – English

## 10 Zero calibration

For correct zero calibration, follow the instructions below.

At the start of a shift, after every **change in the positioning height** of the pressure transducer or patient, and after any interruption in signal transmission between pressure transducer and control monitor (e.g. due to a power failure or change of electrical accessory), another zero calibration must be performed to ensure accurate pressure monitoring.


If, during zero calibration, saline solution leaks **from the lateral duct** of the 3-way zeroing stopcock **without the flush system being operated**, the pressure monitoring set must be inspected and replaced with a new one if required because this is indicative of uncontrolled flushing due to a defective flush system (see section 19, page 66 – Replacing the pressure monitoring set/sterile accessories).

## 11 Transmission quality check


To check the transmission quality, follow the instructions below.

### 11.1 Square wave test

If the completed square wave test reveals that the generated signals exhibit overshoot or are overdamped, the set-up of the pressure monitoring system must be checked. To do this, follow the instructions in section 23, page 69 – Troubleshooting and fault rectification.


### 11.2 Pressure signals

If, after connecting the pressure monitoring set to a patient, it is found that the signals displayed on the patient screen exhibit overshoot or are overdamped, the set-up of the pressure monitoring system must be inspected. To do this, follow the instructions in section 23, page 69 – Troubleshooting and fault rectification.

## 12 S.E.R.O. (Selectable Eliminator of Resonance Overshoot)

Take note of the following illustrations for the correct use of the S.E.R.O. damping stopcock

## 13 Blood sampling

### 13.1 General

- 13.1.1  Due to (incipient) **occlusion of the patient access port**, drawing or returning liquids using a reservoir (NBSS, NBSS<sup>XL</sup>, ABSS & ABSS<sup>XL</sup> configurations) or a blood sampling container can generate too great a vacuum or excess pressure in the tubing system of the pressure monitoring set. This may cause blood cells to be damaged (haemolysis) and blood test results could thus be distorted. Moreover, it can lead to the ingress/leakage of liquids, air and/or pathogens. **Prior to every blood sampling**, check that the patient access port is patent and continuous.
- 13.1.2  The **surface of the sampling port** must be disinfected before and after each use and must be free of any blood residue (spray and/or wipe disinfection) to prevent the ingress of pathogens and blood clots into the pressure monitoring set (**see section 14, page 63 – Disinfection instructions**). Prior to use, also ensure that the used disinfectant has completely evaporated. Disinfectant that has not evaporated can damage the components and lead to distortion of the blood analyses due to haemolysis of the blood samples.
- 13.1.3  Too **rapid drawing or returning of liquids** using a reservoir (NBSS, NBSS<sup>XL</sup>, ABSS & ABSS<sup>XL</sup> configurations) or a blood sampling container can generate too great a vacuum or excess pressure in the tubing system. This may cause blood cells to be damaged (haemolysis) and blood test results could thus be distorted. Moreover, it can lead to the ingress/leakage of liquids, air and/or pathogens. To carry out a complication-free volume shift, the reservoirs and blood sampling containers must always be operated carefully and without any noticeable resistance.
- 13.1.4  **Clotted blood (coagulum)** inside the tubing system can cause clogging of the patient access port and a loss of function of the pressure monitoring set. Moreover, infiltration of coagulum into the patient's circulatory system can cause vascular occlusion. To prevent coagulum from forming inside the pressure monitoring set, all the components involved in blood sampling should be checked for blood residue after the blood sampling process by means of a visual inspection, and all residue must be removed using the rapid flushing function. For pressure monitoring sets with an NBSS, NBSS<sup>XL</sup>, ABSS or ABSS<sup>XL</sup>, remove the reservoir from its designated fixture in order to better inspect the bottom of the plunger stopper for blood residue. If residue is detected, flush the component with pure saline solution while gently opening and closing the control knob.
- 13.1.5  Do not sample blood via the **S.E.R.O. damping stopcock** because this can lead to coagulum forming inside the component. **Clotted blood (coagulum)** can lead to the patient access port becoming clogged and to a functional impairment or loss of function of the component and/or the entire pressure monitoring set. Moreover, infiltration of coagulum into the patient's circulatory system can cause vascular occlusion.
- 13.1.6  No valid pressure transmission is possible throughout the entire blood sampling process due to the various manipulations of the components between the patient connection and pressure transducer with accompanying movement artefacts.


### 13.2 Blood sampling with CLC

Follow the instructions below for blood sampling with a 3-way stopcock, 4-way stopcock or ArterioGuard®.

- 13.2.1  Before each blood sampling of blood, a sufficient quantity of liquid must be removed from the pressure monitoring set using a suitable blood sampling container in order to ensure that undiluted blood can be withdrawn for blood analyses. A blood sample diluted with saline solution leads to invalid analysis results. The minimum volume to be withdrawn can be calculated by multiplying the configuration-dependent liquid volume between the patient vessel and the sampling port (**dead-space volume**) by a factor defined by the relevant healthcare institution. To work out the configuration-dependent dead-space volume, draw liquid into a blood sampling container until blood can be seen at the sampling port. The withdrawn volume is then read off using the scale on the blood sampling container, and multiplied by the defined factor. The determined result specifies the volume that must be removed prior to a withdrawal in order to ensure that undiluted blood samples can be taken. For blood gas analyses, it is recommended that twice the dead-space volume be removed prior to the blood sampling process. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set.

# Instructions for Use – English

13.2.2

## 13.3 Blood sampling with a SWAN switch

Follow the instructions below for blood sampling with a SWAN switch.

13.3.1

Before each blood sampling, a sufficient quantity of liquid must be removed from the pressure monitoring set using a suitable blood sampling container in order to ensure that undiluted blood can be withdrawn for blood analyses. A blood sample diluted with saline solution leads to invalid analysis results. The minimum volume to be withdrawn can be calculated by multiplying the configuration-dependent liquid volume between the patient vessel and the sampling port (dead-space volume) by a factor defined by the relevant healthcare institution. To work out the configuration-dependent dead-space volume, draw liquid into a blood sampling container until blood can be seen at the sampling port. The withdrawn volume is then read off using the scale on the blood sampling container, and multiplied by the defined factor. The determined result specifies the volume that must be removed prior to a withdrawal in order to ensure that undiluted blood samples can be taken. For blood gas analyses, it is recommended that twice the dead-space volume be removed prior to the blood sampling process. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set.

13.3.2

The SWAN switch 4-way stopcock must not be punctured with **needles or cannulas**, as this may damage the valve and cause ingress/leakage of liquids, air and/or pathogens. Only use blood sampling containers or adapters with a Luer interface manufactured in accordance with the applicable standards and which is CE marked.

13.3.3



To ensure correct use of the SWAN Switch 4-way stopcock, only blood sampling **containers and adapters** manufactured in accordance with the applicable standards and which are CE marked may be connected at the side port. Lockable Luer lock systems are recommended for this application. When using Luer slip systems, it must be ensured that they are connected to the side port very carefully, twisting by 90°. Inserting Luer slip systems into the valve on the side port without twisting may damage the valve.

13.3.4

A **45° position of the SWAN switch 4-way stopcock** is not permitted because the flow channel of the sampling point is not closed and this can lead to the ingress/leakage of liquids, air and/or pathogens. For complication-free blood sampling, follow the application steps in 13.3.5.

13.2.5

**13.4 Blood sampling with an NBSS or NBSS<sup>XL</sup>**

Follow the instructions below for blood sampling with an NBSS or NBSS<sup>XL</sup>

13.4.1

For blood gas analyses, it is recommended that twice the **dead-space volume** (configuration-dependent liquid volume between the patient vessel and the sampling port) be removed prior to the blood sampling process. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set. The blood sampling systems (NBSS, NBSS<sup>XL</sup>, ABSS & ABSS<sup>XL</sup>) within the pressure monitoring set of the manufacturer CODAN pvb Critical Care are configured in such a way that with a completely filled reservoir and when using a catheter with an external diameter of ≤ 1.3 mm (18 G/4 Ch) and a length of ≤ 20 cm, at least twice the dead-space volume can be removed. If, however, it is clinically necessary and unavoidable to use a catheter with an external diameter greater than 1.3 mm and/or a length greater than 20 cm, the dead-space volume of the pressure monitoring set must be determined prior to blood sampling to ensure that undiluted blood samples can be taken for blood analyses using the blood sampling system. To work out the configuration-dependent dead-space volume, draw liquid into a blood sampling container until blood can be seen at the intended sampling port. The withdrawn volume is then read off using the scale on the blood sampling container, and multiplied by a factor of 2. The result of this multiplication specifies the volume that must be removed prior to a withdrawal in order to ensure that undiluted blood samples can be taken for blood gas analysis. If the calculated volume cannot be withdrawn from this reservoir, the dead-space volume must be reduced for the application (for example, by using a catheter with a smaller internal volume) or more mixing volume must be removed prior to the blood sampling (for example, by removing mixed blood with an additional blood sampling container). A blood sample diluted with saline solution leads to invalid blood analysis results.

13.4.2

It should be ensured that the **2 ml reservoir** (NBSS & ABSS) of a pressure monitoring set is completely filled during blood sampling, since this is the only way to ensure that undiluted blood samples can be taken for a blood gas analysis. It is not necessary to completely fill a **6 ml reservoir** (NBSS<sup>XL</sup> & ABSS<sup>XL</sup>) if the configuration-dependent dead-space volume was determined using the method set out in the previous paragraph before the sampling, in order to remove at least twice the dead-space volume before the blood sampling using the scale on the 6 ml reservoir for the blood gas analysis. A blood sample diluted with saline solution leads to invalid blood analysis results. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set before blood sampling.

13.4.3

**If the control knob of a reservoir is turned too far**, it can damage the threaded rod and thus make the entire component permanently unusable. A predetermined breaking point on the threaded rod protects the housing from excessive force in order to maintain the leak-tightness of the system at all times. The control knob of a reservoir must only be opened and closed gently and until noticeable resistance is felt.


# Instructions for Use – English

13.4.4  The SWAN switch 4-way stopcock must not be punctured with **needles or cannulas**, as this may damage the valve and cause ingress/leakage of liquids, air and/or pathogens. Only use blood sampling containers or adapters with a Luer interface manufactured in accordance with the applicable standards and which is CE marked.

13.4.5  To ensure correct use of the SWAN Switch 4-way stopcock, only blood sampling **containers and adapters** manufactured in accordance with the applicable standards and which are CE marked may be connected at the **side port**. Lockable Luer lock systems are recommended for this application. When using Luer slip systems, it must be ensured that they are connected to the side port very carefully, twisting by 90°. Inserting Luer slip systems into the valve on the side port without twisting may damage the valve.


13.4.6  A **45° position of the SWAN switch 4-way stopcock** is not permitted because the flow channel of the sampling point is not closed and this can lead to the ingress/leakage of liquids, air and/or pathogens. For complication-free blood sampling, follow the application steps in 13.4.7 (NBSS) or 13.4.8 (NBSS<sup>XL</sup>).

## 13.4.7 Blood sampling with an NBSS

### 13.4.8 Blood sampling with an NBSS<sup>XL</sup>

### 13.5 Blood sampling with an ABSS or ABSS<sup>XL</sup>

Follow the instructions below for blood sampling with an ABSS or ABSS<sup>XL</sup>

13.5.1

For blood gas analyses, it is recommended that twice the **dead-space volume** (configuration-dependent liquid volume between the patient vessel and the sampling port) be removed prior to the blood sampling process. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set. The blood sampling systems (NBSS, NBSS<sup>XL</sup>, ABSS & ABSS<sup>XL</sup>) within the pressure monitoring set of the manufacturer CODAN pvb Critical Care are configured in such a way that, with a completely filled reservoir and using a catheter with an external diameter of  $\leq 1.3$  mm (18 G/4 Ch) and a length of  $\leq 20$  cm, at least twice the dead-space volume can be removed. If, however, it is clinically necessary and unavoidable to use a catheter with an external diameter greater than 1.3 mm and/or a length greater than 20 cm, the dead-space volume of the pressure monitoring set must be determined prior to blood sampling to ensure that undiluted blood samples can be taken for blood analyses using the blood sampling system. To work out the configuration-dependent dead-space volume, draw liquid into a blood sampling container until blood can be seen at the intended sampling port. The withdrawn volume is then read off using the scale on the blood sampling container, and multiplied by a factor of 2. The result of this multiplication specifies the volume that must be removed prior to a withdrawal in order to ensure that undiluted blood samples can be taken for blood gas analysis. If the calculated volume cannot be withdrawn from this reservoir, the dead-space volume must be reduced for the application (for example, by using a catheter with a smaller internal volume) or more mixing volume must be removed prior to the blood sampling (for example, by removing mixed blood with an additional blood sampling container). A blood sample diluted with saline solution leads to invalid blood analysis results.

13.5.2

It should be ensured that the **2 ml reservoir** (NBSS & ABSS) of a pressure monitoring set is completely filled during blood sampling, since this is the only way to ensure that undiluted blood samples can be taken for a blood gas analysis. It is not necessary to completely fill a **6 ml reservoir** (NBSS<sup>XL</sup> & ABSS<sup>XL</sup>) if the configuration-dependent dead-space volume was determined using the method set out in the previous paragraph before the sampling, in order to remove at least twice the dead-space volume before the blood sampling using the scale on the 6 ml reservoir for the blood gas analysis. A blood sample diluted with saline solution leads to invalid blood analysis results. For coagulation diagnostics, more mixing volume may need to be removed from the pressure monitoring set before blood sampling.

# Instructions for Use – English

13.5.3


Turning the control knob of the reservoir too far can damage the threaded rod and thus make the entire component permanently unusable. A predetermined breaking point on the threaded rod protects the housing from excessive force in order to maintain the leak-tightness of the system at all times. The control knob of a reservoir must only be opened and closed gently and until noticeable resistance is felt.

13.5.4


The sampling port of the ABSS, and ABSS<sup>XL</sup> must not be punctured with **needles or cannulas**, as this may damage the plastic membrane, causing ingress/leakage of liquids, air and/or pathogens. To puncture the membrane, only use the blood sampling adapter HG-01 (REF: 74.6934) together with syringes and/or equivalent blood sampling containers with a Luer connector.

## 13.5.5 Blood sampling with an ABSS

13.5.6 Blood sampling with an ABSS<sup>XL</sup>

**14 Disinfection instructions**

- 14.1  For disinfection, use **iodine-containing disinfectants** or standard **1-propanol, 2-propanol or ethanol-based disinfectants**. For reasons of user and patient safety, check with the manufacturer about the product-specific contact time of the disinfectant that must be complied with after application in order to adequately kill or inactivate pathogens. Disinfection without observing the required contact time may lead to the spread/transmission of pathogens and thus pose a danger to patients, users and third parties.
- 14.2  Following disinfection, make sure that the applied disinfectant has completely evaporated before using the components. If Luer interfaces are used without complete **evaporation of the disinfectant**, they may break, causing ingress/leakage of liquids, air and/or pathogens.
- 14.3 CODAN pvb Critical Care pressure monitoring sets are developed and manufactured in such a way that when the above-mentioned disinfectants are used and the product-specific contact and drying times and the maximum duration of use of the pressure monitoring set are complied with (**see section 15, page 64 – Duration of use**), any number of disinfections can be carried out without reducing the basic safety of the pressure monitoring set or damaging the affected component(s).
- 14.4  **Wetting the air separation membrane of the AeroGuard® filter system** with liquids (e.g. disinfectant or saline solution) before the component has been completely primed with saline solution can result in the AeroGuard® not allowing any more saline solution to pass through, thus making it impossible to prime and flush the pressure monitoring set. Before and during filling, make sure that no liquids can wet the filter system in order to maintain the patency of the pressure monitoring set.

# Instructions for Use – English

## 15 Duration of use

- 15.1  For reasons of hygiene, **the pressure monitoring set must be changed** at least every 96 hours and the cumulative duration of use of 30 days maximum should not be exceeded. Moreover, in the event of damage and/or any change in the performance of the individual components, e.g. due to environmental conditions, the pressure monitoring set must be replaced (**see section 19, page 66 – Replacing the pressure monitoring set/sterile accessories IBPM**).
- 15.2   
 **Reusing or reprocessing** the pressure monitoring set and sterile accessories is prohibited under all circumstances, since this can lead to lasting health damage due to infection or allergic/toxic reactions, accompanied by chronic long-term damage to tissues/organs. Reusing the blood sampling systems may result in clotting of the systems, which makes proper operation impossible. Reuse can also result in distortion of blood values. In addition, compliance with the recommended service life cannot be ensured if systems are reused. Moreover, material properties can be impaired by cleaning processes and as a result, the correct display of measured values, the leak-tightness of adhesively bonded and screw connections and the pressure resistance of the individual components, for example, can no longer be ensured. As a result, neither correct set-up nor proper use of the pressure monitoring set without air pockets and/or leakage of liquids is possible.
- 15.3   
 **Resterilisation** of the pressure monitoring set, accessories and/or individual components is prohibited. Otherwise the basic safety and essential performance characteristics of all components cannot be ensured, which will jeopardise patient safety.

## 16 Technical specifications and operating conditions

### 16.1 Electrical safety

- 16.1.1  → The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series have an **IP37 enclosure protection rating**. They are therefore protected against accessing dangerous parts with a tool and against the harmful penetration of solid foreign bodies with a diameter  $\geq 2.5$  mm as well as water in the event of temporary submersion (30 minutes in water 1 m deep).
- 16.1.2  → The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series have an **IK04 impact rating** and thus enable the measuring function and essential performance characteristics to be maintained after severe mechanical impact stress of up to 0.5 joule.
- 16.1.3  → The pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series are electrically isolated and have defibrillation-protected applied parts of the **CF (Cardiac Float) type**, and thus fulfil the strictest requirements for the insulation of accessible parts.

### 16.2 Conditions for continuous operation

- 16.2.1 For the correct continuous operation of the pressure monitoring set with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series, please take note of the technical specifications on page 2 of these Instructions for Use.
- 16.2.2  It is prohibited to use the pressure monitoring set and/or the accessories **outside the permitted conditions for continuous operation** (for example, in the event of transient operating conditions with irregular or severely fluctuating environmental parameters) because it may lead to loss of the essential performance characteristics and basic safety.

### 16.3 Electromagnetic compatibility (EMC)

The pressure monitoring sets with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series meet the requirements of the European standards EN 60601-1-2:2015 and EN 60601-2-34:2014, which are required as proof of electromagnetic compatibility (EMC), as well as the additional requirements in accordance with EN 60601-1-12:2015 + A1:2020 for the use of medical electrical equipment in an emergency setting. They are intended for operation within a defined and controlled electromagnetic environment. This environment and the tests required as proof of electromagnetic compatibility in accordance with EN 60601-1-2:2015, EN 1789:2020, EN 13718-1:2014+A1:2020, EUROCAE ED-14G:2015/RTCA/DO-160G:2010 + Cor.1:2014 (Section 8.5 Category S, Section 8.8 Category U, Section 20, Section 21 Category M) and EN 60601-2-34:2014 were successfully performed and are presented in the manufacturer's EMC guidance. The EMC guidance will be provided by the distributor and/or manufacturer upon request.

16.3.1

Temporary impairments due to electromagnetic interference are permissible, but normal operation must be ensured again after a maximum of 10 seconds. A reduction or loss of the essential performance characteristics due to **use outside the defined environmental conditions** and significantly greater electromagnetic interference can lead to an unacceptable risk for patients, users or third parties, for example:

- the loss of essential performance characteristics during use may lead to the patient experiencing circulatory crisis if there is a lack of monitoring;
- an electric shock may lead to the patient experiencing ventricular fibrillation and/or cardiovascular failure;
- the malfunction of other equipment may lead to cardiovascular failure or death of the patient due to the failure of vital equipment;
- if incorrect values are displayed (falsely high/falsely low), this may lead to severe systemic harm to the patient due to the incorrect or omitted administration of medication.

The responsible organisation must ensure that the pressure monitoring set with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series are only used within the environment described in the EMC guidance from the manufacturer CODAN pvb Critical Care.

16.3.2

Using the pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series **together with and/or immediately next to other stacked equipment** must be avoided because this may cause a reduction in the essential performance characteristics. However, if such an application is clinically necessary and unavoidable, all the equipment must be observed during use to make sure that it is functioning as intended.

16.3.3

The **use of electrical accessories from other companies** may result in increased electromagnetic interference emissions and/or decreased electromagnetic immunity of the pressure monitoring set, and lead to a reduction in the essential performance characteristics. Only use original CODAN pvb Critical Care accessories together with pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series (**see section 4.3, page 46 – Non-sterile accessories**).

16.3.4

**Portable radio equipment and other RF communications equipment** (for example, mobile telephones including their antenna cables, external antennas and/or other accessories) must not be used within 30 cm of the electrical accessories and pressure monitoring sets with pressure transducers of the **DPT-6000** or **DPT-9000/Xtrans®** series. Failure to observe this distance may lead to a reduction in the performance characteristics of the pressure monitoring set.

16.3.5

The pressure monitoring sets can be used during the **application of electrosurgical units** (300 W monopolar cutting, 100 W monopolar soft coagulation, 337 kHz). However this may lead to brief malfunctions and a reduction in the essential performance characteristics of the pressure monitoring set. All equipment must be observed during use to ensure that it is functioning as intended.

16.3.6

Using a pressure monitoring set **without wearing medical protective gloves and/or outside the environment described in the EMC guidance** from the manufacturer CODAN pvb Critical Care may lead to a reduction or loss of the essential performance characteristics due to harmful electrostatic discharge. The pressure monitoring set should only be used with medical protective gloves and within the environment described in the EMC guidance from the manufacturer CODAN pvb Critical Care in order to prevent harmful electrostatic discharge.

#### 16.4 Compatibility with electrosurgical units

The pressure monitoring sets can be used during the application of an electrosurgical unit. However this may lead to brief malfunctions and a reduction in the essential performance characteristics of the pressure monitoring set (**see section 16.3, page 64 – Electromagnetic compatibility**).

#### 16.5 Compatibility with magnetic fields (magnetic resonance imaging)

16.5.1

Non-clinical investigations have shown that the pressure monitoring sets **with pressure transducers from series DPT-6000** and **DPT-9000 (Xtrans®)** are MR conditional according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). A patient can be scanned together with a pressure monitoring set within the MR environment, under the following conditions:

- Static magnetic field of 1.5 Tesla, 3 Tesla and 7 Tesla with
- Maximum spatial gradient magnetic field of 24,600 G/cm (246 T/m)
- Maximum gradient field product of 1,033,000,000 G<sup>2</sup>/cm (1,033 T<sup>2</sup>/m)

It was possible to maintain the measuring accuracy of pressure transducers ( $\pm 4\%$  /  $\pm 4$  mmHg), as required by the European standard EN 60601-2-34, during MR-induced malfunction tests in which the devices were exposed to the static field B<sub>0</sub> at clinical level (200 mT), the time-varying magnetic field gradient (dB/dt) and the RF field (B<sub>1</sub>, E) of 1.5 Tesla, 3 Tesla and 7 Tesla MR systems. The devices have passed all of the malfunction tests according to the proprietary test procedures at CODAN pvb Critical Care GmbH.

All components made of plastic or other non-metallic materials in the pressure monitoring sets (e.g. pressure monitoring lines, one-way or multi-way stopcocks and reservoir systems) are MR safe and can be safely placed in the MRI tube. The MR image quality may be compromised by an incorrect proton signal if the plastic components are placed within or around the imaging area.

The DPT-6000 and DPT-9000 (Xtrans®) series pressure transducers and micro drip chamber filling systems must not be placed inside the MRI tube.

# Instructions for Use – English

- 16.5.2  → Non-clinical investigations have shown that the **systems for filling using a micro drip chamber** are **MR conditional** according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). A patient can be scanned together with a filling system with micro drip chamber within the MR environment, under the following conditions:
- Static magnetic field of 1.5 Tesla, 3 Tesla and 7 Tesla with
  - Maximum spatial gradient magnetic field of 24,600 G/cm (246 T/m)
  - Maximum gradient field product of 1,033,000,000 G<sup>2</sup>/cm (1,033 T<sup>2</sup>/m)
- The systems for filling using a micro drip chamber must not be placed inside the MRI tube.
- 16.5.3  → Non-clinical investigations have shown that the following **sterile accessories made of plastic or other non-metallic materials** are **MR safe** according to the established terminology from the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20):
- Systems for filling with a macro drip chamber, systems for filling with a spike, Heidelberger extensions, spiral lines, lines for filling, blood sampling adapters, pressure monitoring lines, arterial blood sampling systems (ABSS), arterial blood sampling systems XL (ABSS<sup>XL</sup>), needle-free blood sampling systems (NBSS), needle-free blood sampling systems XL (NBSS<sup>XL</sup>), reservoir systems, reservoir systems XL, sampling extensions, pressure monitoring set additions, CVP measuring limbs, one-way stopcocks, multi-way stopcocks, and sampling ports.
- The listed product groups can be safely placed in the tube of the MR system. The MR image quality may be compromised by an incorrect proton signal if the plastic components are placed within or around the imaging area.
- 16.5.4  → CODAN pvb Critical Care devices labelled with the symbol shown opposite are considered **MR unsafe** according to the established terminology of the International Electrotechnical Commission (IEC 62570) and the American Society for Testing and Materials (ASTM F2503-20). No medical device labelled with this symbol may be used within the MR environment.

## 17 Storage and transportation

- 17.1  → For the correct **storage and transportation** of the pressure monitoring sets with pressure transducers of the **DPT-6000** and **DPT-9000/Xtrans®** series, please observe the storage and transportation symbols on page 2 of these Instructions for Use and on the protective or sterile packaging of the device. Storage or transportation of the pressure monitoring set and/or accessories outside of the permitted environmental conditions is prohibited because it may lead to reduced useful life and loss of the essential performance characteristics and basic safety.
- 17.2  → Pressure monitoring sets and sterile accessories that are not immediately aseptically prepared, primed and used must not be transported or stored **if the sterile packaging has already been opened**. They must be disposed of because it is no longer possible to guarantee that the sterile condition has been maintained, and contamination with pathogens is possible (**see section 21, page 68 – Disposal**).

## 18 Residual risks

Based on the completed risk assessment in accordance with the EN ISO 14971:2019 + A11:2021 standard, the following residual risks arise with regard to the use of pressure monitoring sets and sterile accessories from the manufacturer CODAN pvb Critical Care: ventricular fibrillation, cardiovascular failure or crisis, drop in blood pressure, drop in cardiac output, drop in oxygen saturation, blood loss, severe systemic damage, chronic long-term damage to tissues and organs, tissue necrosis, thrombosis, (air) embolisms, spread of pathogens, local or systemic infections, vasospasms, sepsis, transient fever, local and allergic reactions, hypersensitivity reactions, pulmonary oedema, malignant tumors, genetic damage, hypervolaemia, malpuncture, vascular injuries, damage to visual organs, burns, or death of the patient.

All residual risks as well as the overall risk fall within the acceptable range after risk mitigation and cannot be reduced any further. The medical/clinical benefit of the pressure monitoring sets and sterile accessories far outweighs the remaining residual risks or the remaining overall risk. It is the responsibility of the attending physician to inform the patient about possible residual risks.

## 19 Replacing the pressure monitoring set/sterile accessories

During replacement, the safety of the patient must be ensured. Before a pressure monitoring set or accessories can be replaced, the new one(s) must be primed such that they are free of air bubbles. After the patient access port has been closed, the pressure monitoring set or accessories can be disconnected by hand (without tools) and the new one(s) connected. Please follow the instructions in sections 6–13 concerning safe provision, preparation, filling, initial operation and use. The replaced pressure monitoring set and/or accessories must be disposed of (**see section 21, page 68 – Disposal**).


## 20 Ending use of a pressure monitoring set

- 20.1 **General**
- 20.1.1  → To end use, follow the hygiene guidelines specified by the responsible organisation and **wear medical protective gloves** to prevent the spread/transmission of pathogens.
- 20.1.2  → When the use of a pressure monitoring set has to be terminated, for example because there is no longer a medical indication, or in the case of damage or a change in performance, or if the duration of use is exceeded, and the set has come into contact with **potentially infectious substances/materials** (e.g. blood, saliva and/or bodily excretions), it poses a risk and must be disposed of according to the regulations for contaminated/infectious medical waste (biohazard) to prevent the spread/transmission of pathogens (**see section 21, page 68 – Disposal**).



**20.2 Uncoupling a pressure monitoring set from the patient access port**

20.2.1 Disconnecting a pressure monitoring set from an **unsecured patient access port** can lead to uncontrolled bleeding, a severe drop in blood pressure and the need for a blood transfusion. Before uncoupling a pressure monitoring set, make sure that the patient access port is closed to prevent uncontrolled bleeding and to guarantee that it can be disconnected safely.


**20.3 Uncoupling a reservoir from non-sterile accessories**

**20.4 Uncoupling a pressure transducer (DPT-6003/6009 & DPT-9003/9009) from non-sterile accessories**

20.4.1 **Exposed contacts** of the electrical accessories must not be touched at the same time as the patient while the accessories are making an electrically conductive connection to a control monitor. Risk of electric shock! Make sure that an electrically conductive connection cannot be established between the patient and the exposed contacts of the electrical accessories.


20.4.2 **Spilling of liquids**, whatever the quantity, onto electrical contacts can lead to the loss of basic safety, including the essential performance characteristics. Make sure that no liquids can wet the electrical contacts of the accessories when ending use of the electrical accessories. In the event of liquid contact, clean and dry the electrical accessories and check their function.

**20.5 Uncoupling a pressure monitoring set with filling and irrigation system from a saline bag (DPT-6003, DPT-9003)**

20.5.1 Please note that the instructions in this section are only applicable for pressure monitoring sets with a filling and irrigation system and pressure transducer types DPT-6003 or DPT-9003. Pressure monitoring sets with pressure transducer types DPT-6009 or DPT-9009 should only be used with syringe pumps (see section 20.6, page 68 – Uncoupling a pressure monitoring set with filling line from a syringe (DPT-6003/6009 & DPT-9003/9009)).

# Instructions for Use – English


20.5.2

20.5.3  **The piercing device of a filling and irrigation system must not be disposed of unprotected because this could result in a needlestick injury with accompanying infection during or after disposal. Prior to disposal, push the piercing device into the opening of the roller clamp to prevent needlestick injuries and guarantee safe disposal.**

## 20.6 Uncoupling a pressure monitoring set with filling line from a syringe (DPT-6003/6009 & DPT-9003/9009)

20.6.1 Please note that the instructions in this section are only applicable for pressure monitoring sets with a filling line and pressure transducer types DPT-6003/6009 & DPT-9003/9009.

## 21 Disposal

21.1  For disposal, follow the **hygiene guidelines** specified by the responsible organisation and wear **medical protective gloves** to prevent the spread/transmission of pathogens.

21.2 Pressure monitoring sets, sterile accessories, consumables and packaging/transportation materials (e.g. protective & sterile packaging, mounting brackets, paper sleeves/labels, rubber bands and transportation strips) that must be disposed of **prior to use** due to damage or because they have passed their expiry date should be disposed of according to the regulations of the responsible organisation.

21.3  Pressure monitoring sets, sterile accessories, consumables and packaging/transportation materials (e.g. protective & sterile packaging, mounting brackets, paper sleeves/labels, rubber bands and transportation strips) that have to be disposed of **during or after use** because there is no longer a medical indication, or in the case of damage or a change in performance, or if the duration of use is exceeded, and they have come into contact with **potentially infectious substances/materials** (e.g. blood, saliva and/or bodily excretions), they pose a risk and must be disposed of according to the regulations for contaminated/infectious medical waste (biohazard) to prevent the spread/transmission of pathogens.

## 22 Customer service

Contact the manufacturer or your distributor for more information.

CODAN pvb Critical Care GmbH  
Römerstraße 18 · 85661 Forstinning, Deutschland  
Tel. +49 (0) 81 21 – 98 020  
codan@codanpvbcc.de · www.codancompanies.com

**23 Troubleshooting and fault rectification**

Phase	Possible complications	Possible causes	Fault rectification <b>A</b> / prevention <b>B</b>	IFU
Filling	Filling of the pressure monitoring set not possible or only to a limited extent	Insufficient rapid flushing	<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>A</b> Identify twisted or kinked tubes and correct if required <b>A</b> Check the stopcock positions and correct if required <b>A</b> Check the patency of the filling and irrigation system and correct if required	9 - 8 7
	Uncontrolled ingress/leakage of saline solution and/or air at Luer interfaces	Loose Luer interfaces	<b>A</b> Check the Luer interfaces and tighten by hand if required	6
Pressure monitoring	Signal transmission or zero calibration not possible	Flow channel not patent	<b>A</b> Identify twisted or kinked tubes and correct if required <b>A</b> Check the stopcock positions and correct if required	- 8
		Incorrect position or incipient occlusion of the patient access port	<b>A</b> Check the patency of the patient access port and replace if required	9
			<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>B</b> Carefully remove all blood residue from the pressure monitoring set after blood sampling	9 13
		Pressure transducer incorrectly connected to the electrical accessories	<b>A</b> Check the connection between the pressure transducer and electrical accessories and correct if required	9
	Electrical accessories damaged or defective	<b>B</b> Always handle, clean and store electrical accessories with the utmost care	-	
	Pressure curve displayed on the control monitor without pressure monitoring values	Incorrect assembly and/or zero calibration of the pressure monitoring set	<b>A</b> Check the positioning height of the pressure transducer and correct if required and carry out another zero calibration	9, 10
		Falsely high/low measured values (overshoot)	Flow channel not patent	<b>A</b> Identify twisted or kinked tubes and correct if required <b>A</b> Check the stopcock positions and correct if required
		Incorrect position or incipient occlusion of the patient access port	<b>A</b> Check the patency of the patient access port and replace if required	9
			<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>B</b> Carefully remove all blood residue from the pressure monitoring set after blood sampling	9 13
		Air bubbles in the tubing system	<b>A</b> Inspect the pressure monitoring set and, if required, remove air bubbles by means of rapid flushing or by gently tapping the affected components	8
		Incorrect assembly and/or zero calibration of the pressure monitoring set	<b>A</b> Check the positioning height of the pressure transducer and correct if required and carry out another zero calibration	9, 10
		Use of unapproved sterile accessories to supplement, extend or change the configuration and/or components of the pressure monitoring set	<b>B</b> Only use approved CODAN pvb Critical Care sterile accessories and contact the manufacturer or distributor for more information if required	4
	Blood forced into the tubing system	Insufficient flushing	<b>B</b> Use an S.E.R.O. damping stopcock	12
<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>A</b> Check the stopcock positions and correct if required <b>A</b> Check the patency of the filling system and correct if required			9 8 7.2	
Loose Luer interfaces		<b>A</b> Check the Luer interfaces and tighten by hand if required	6	
Blood sampling	Ingress/leakage of liquids and/or air at the flush system after opening/closing a reservoir	Flow channel not patent	<b>A</b> Identify twisted or kinked tubes and correct if required <b>A</b> Check the stopcock positions and correct if required	- 8
		Incorrect position or incipient occlusion of the patient access port	<b>A</b> Check the patency of the patient access port and replace if required	9
			<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>B</b> Carefully remove all blood residue from the pressure monitoring set after blood sampling	9 13
	Control knob turned too quickly	<b>B</b> Always operate the control knob of the reservoir carefully and without any noticeable resistance <b>B</b> Before operating the reservoir, turn the OFF tap of the 3-way zeroing stopcock in the direction of the patient (NBSS & ABSS)	13 -	
	Ingress/leakage of liquids and/or air at the sampling port after opening/closing a reservoir	Flow channel not patent	<b>A</b> Identify twisted or kinked tubes and correct if required <b>A</b> Check the stopcock positions and correct if required	- 8
		Incorrect position or incipient occlusion of the patient access port	<b>A</b> Check the patency of the patient access port and replace if required	9
			<b>A</b> Check the pressurisation of the saline bag or the flow rate of the syringe pump and correct if required <b>B</b> Carefully remove all blood residue from the pressure monitoring set after blood sampling	9 13
	Control knob turned too quickly	<b>B</b> Always operate the control knob of the reservoir carefully and without any noticeable resistance	13	
	Damaged sampling port after blood sampling	Use of unsuitable blood sampling containers	<b>B</b> Use of blood sampling containers and adapters manufactured in accordance with the applicable standards and which are CE marked	5, 13
		Use of cannulas/needles	<b>B</b> Use of blood sampling containers with Luer interfaces	13
Damaged Luer interface with stress cracks after blood sampling	Use of unapproved disinfectant	<b>B</b> Use of iodine-containing disinfectants or standard 1-propanol, 2-propanol or ethanol-based disinfectants	14	
	Inadequate evaporation of the disinfectant	<b>B</b> Before using the relevant components, allow the applied disinfectant to completely evaporate	14	

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CODAN Medical AG