

CODAN Oral syringes



Smooth
Precise
Reliable

Confusion excluded

CODAN Oral syringes have an extra wide oral cone that makes the syringe distinctive and prevents accidental connection to an administration system.

MADE IN DENMARK



CODAN REF	Description
62.1900	0.5 ml Oral Clear, centric 0.01 grad
62.1940	0.5 ml Oral Amber, centric 0.01 grad
62.1800	1 ml Oral Clear, centric 0.01 grad
62.1880	1 ml Oral Amber, centric 0.01 grad
62.3800	3 ml Oral Clear, centric 0.1 grad
62.3810	3 ml Oral Amber, centric 0.1 grad
62.4960	5 ml Oral Clear, centric 0.2 grad
62.4980	5 ml Oral Amber, centric 0.2 grad
62.6960	10 ml Oral Clear, eccentric 0.2 grad
62.6950	10 ml Oral Amber, eccentric 0.2 grad
62.7800	20 ml Oral Clear, eccentric 1.0 grad
62.7860	20 ml Oral Amber, eccentric 1.0 grad
62.8560	50/60 ml Oral Clear, eccentric 1.0 grad
62.8550	50/60 ml Oral Amber, eccentric 1.0 grad

The non-sterile CODAN Oral syringe program consist of reliable, precise, high-quality single use syringes, includes a range with oral tip dispensers with clear and amber barrel for UV light protection of sensitive medications. The CODAN Oral syringes also meets the highest standards for feeding premature infants when breast milk has to be administered via an oral syringe and provide support for sensitive handling. This applies also to all other areas in which CODAN Oral syringes are used.

The amber coloured CODAN Oral syringes provide optimal light protection, allowing safe oral administration of light-sensitive medications.

The CODAN characteristics apply to this range:

- Three-part syringe with thin silicone ring
- Easy and very even withdrawal and emptying
- Smooth plunger guidance
- Precise scale
- Reliable plunger stop
- Dual scale printing
- Environmentally friendly and PVC-free materials
- Latex-free
- Balanced portfolio with centric and eccentric connections



Manufacturer

CODAN Medical ApS
 Færgevej 4 · 4970 Rødby, Denmark
 Tel. +45 54 65 74 11
 codan@codanmedical.dk · www.codancompanies.com

The decisive connection