



CODAN

Luer & Wound and bladder-irrigation syringes



Smooth Precise Reliable

Confusion excluded

CODAN Luer syringes for simple and safe use are ideally supplemented with all common connectors and administration sets.

MADE IN DENMARK

CODAN Medical ApS · Færgevej 4 · 4970 Rødby, Denmark

Tel. +45 54 65 74 11 · codan@codanmedical.dk · www.codancompanies.com



CODAN REF	Description
62.1912	0.5 ml Luer U-100 Insulin, centric 0.02 grad
62.1612	1 ml Luer TBC, centric 0.01 grad
62.1640	1 ml Luer U-100 Insulin, centric 0.02 grad
62.2611	2 ml Luer, centric 0.1 grad
62.2659	2 ml Luer, centric 0.1 grad
62.8437	50/60 ml Luer, eccentric 1.0 grad
62.8436	50/60 ml Catheter Tip, centric 1.0 grad

The CODAN characteristics apply to this range:

- Three-part syringe with thin silicone ring
- Easy and very even withdrawal and emptying
- Smooth plunger guidance
- Precise scale
- Reliable plunger stop
- Dual scale printing
- Environmentally friendly and PVC-free materials
- Latex-free
- Balanced portfolio with centric and eccentric connections

CODAN Luer syringes are a natural part of the complete CODAN Product range, as they can be used in connection with a lot of other CODAN products and systems like:

- SWAN-LOCK®

The CODAN Luer syringe program as well as wound and bladder-irrigation syringes with catheter tip consist of reliable, precise, high-quality single use syringes, whereby a simple, quick and safe handling with maximum comfort to the users are safeguarded.

Manufacturer

CODAN Medical ApS
 Færgevej 4 · 4970 Rødby, Denmark
 Tel. +45 54 65 74 11
 codan@codanmedical.dk · www.codancompanies.com



The decisive connection