


CODAN Electrical Accessory



Reliability and flexibility in
IBPM signal transmission





CODAN pvb Critical Care GmbH

Electrical Accessory 2023

Reliability and flexibility in
invasive blood pressure monitoring
signal transmission

EA-2023-Q1-04

Copying or reproduction is prohibited!

Table of contents

Electrical accessory basics	8
Plug types overview	10
Products for dedicated applications	14
A	14
Abbott Laboratories	14
ADInstruments	14
American Optical	16
Agilent Technologies	17
Aries Medical	17
Arrow	18
Axcent Medical	18
B	19
Bard/Boston Scientific	19
BCI	19
BLT (Biolight)	20
Bionet	20
Bruker	21
C	21
Colin	21
Comen	23
Criticon	23
D	24
Datascope	24
Datex	25
Dräger	26
E	30
Econet Medical	30
Edan Instruments	31
Em-tec	32
F	32
Fazzini	32
Fukuda Denshi	33

G	34
GE (General Electric)	34
GS Elektromedizinische Geräte	35
Guardian	36
H	37
Hellige	37
HP (Hewlett-Packard)	38
HME	38
Hoyer	39
I	40
Invivo	40
Ivy Biomedical	40
J	41
Jostra	41
K	42
Kontron	42
L	42
Levitronix	42
Litton	43
M	43
Marquette	43
Mediana	44
Medical Data Electronics	45
Medos	46
Medrad	47
Medtronic	47
Mennen	48
Mindray	48
MIPM	50
MR-Equipment	51

Table of contents

N	51
Nihon Kohden	51
Nippon Colin	52
O	53
Ohmeda	53
P	54
Penlon	54
Philips	55
Physio-Control	57
Propaq	57
Protocol	58
Prucka	59
Pulsion	60
R	61
Raumedic	61
RDT (Remote Diagnostic Technologies)	61
Riegel Medical	62
Roche Medical	62
S	63
Saadat	63
Sarns	64
Schiller	65
Schwarzer	66
Siare Group	67
Siecov	68
Siemens	69
S&W (Simonsen & Weel)	71
Spacelabs	72
Stöckert	73
St. Jude Medical	74
T	75
Tektronix	75
Terumo	76

V	76
Vytech	76
W	77
Welch Allyn	77
Witt	77
Z	78
Zoll	78
Appendix for monitor models	79







Electrical accessory basics

Contactboards

Contactboards are the electrical docking stations for pressure transducer of the series Xtrans®. These Contactboards have no monitor plugs, thus they require separately available Monitor-Connectors in order to transmit pressure signals to the patient monitor. 1- to 2-fold Contactboards can be connected

with 1- to 2-fold Monitor-Connectors as well as 3- to 4-fold Contactboards with 3- to 4-fold Monitor-Connectors.

These boards have up to 4 Xtrans® transducer slots and can be purchased in different lengths (see table below).

Product type	REF	Historical no.	Length		GTIN	Picture
Contactboard, 1-fold - FX	75.1255	BDC901FX-M150	150 cm	1	4046765111625	
	75.1204	BDC901FX-M	320 cm	1	4046765108441	
Contactboard, 1-fold - Standard	75.1254	BDC901-M150	150 cm	1	4046765111618	
	75.1002	BDC901-M	320 cm	1	4046765087982	
	75.1205	BDC901-M500	500 cm	1	4046765106584	
Contactboard, 2-fold - Standard	75.1207	BDC902-M150	150 cm	1	4046765106607	
	75.1068	BDC902-M	320 cm	1	4046765025397	
	75.1240	BDC902-M500	500 cm	1	4046765109400	
Contactboard, 3-fold - Standard	75.1079	BDC903-M150	150 cm	1	4046765025724	
	75.1078	BDC903-M	320 cm	1	4046765025717	
	75.1261	BDC903-M500	500 cm	1	4046765427207	
Contactboard, 4-fold - Standard	75.1209	BDC904-M150	150 cm	1	4046765107727	
	75.1087	BDC904-M	320 cm	1	4046765026042	
	75.1252	BDC904-M500	500 cm	1	4046765111601	

Electrical accessory basics

Xtrans® Monitor-Connectors

Monitor-Connectors are the complements for Contactboards, allowing signal transmission between Xtrans® transducer and a patient monitor when combined. 1- to 2-fold Monitor-Connectors can be connected with 1- to 2-fold Contactboards as well as

3- to 4-fold Monitor-Connectors with 3- to 4-fold Contactboards. These connectors can be purchased with a variety of monitor interfaces (see 'Plug types overview' on the next pages) and have a standard length of 50 cm.



Monitor-Connector,
1-fold - Standard/FX



Monitor-Connector,
2-fold - Standard



Monitor-Connector,
3-fold - Standard



Monitor-Connector,
4-fold - Standard

Xtrans® Contactboards-GO

Contactboards-GO are an alternative solution to the combination of Contactboards and Monitor-Connectors, as they already have built-in monitor plugs. These boards can also be purchased with up to 4 Xtrans® pressure transducer slots. The standard lengths

of the boards range from 200 to 370 cm, depending on the application and the individual products. The monitor plugs assembled with the Contactboards-GO are shown in the category 'Plug types overview' on the next pages of this document.




Contactboard-GO - FX,
1-fold




Contactboard-GO,
1-fold



Contactboard-GO,
2-fold



Contactboard-GO,
3-fold



Contactboard-GO,
4-fold


DPT-6000 Interface cables

Interface cables allow signal transmission between DPT-6000 pressure transducer and a patient monitor. Usually, one interface cable is required per pressure transducer. However, for special patient monitors with double pin-out IBPM sockets, 2-fold interface cables can be purchased. The standard lengths

of the interface cables range from 200 to 350 cm, depending on the application and the individual products. The monitor plugs assembled with interface cables are shown in the category 'Plug types overview' on the next pages of this document.

















Interface cable,
1-fold



















Interface cable,
2-fold






Plug types overview

Pins	Colour	Material	Plug type	Picture	Usage
4	Red	Plastic	210		BLT (Biolight) page 20
4	Red	Plastic	176		Schiller page 66
4	Silver, black	Metal	110		Stöckert page 73
5	Silver, grey	Metal	142		Sarns page 64 Terumo page 76
5	Silver, blue	Metal	86		Medrad page 47
5	Silver, blue	Metal	168		Jostra page 41
5	Grey, blue	Plastic	152		Colin page 22 Schwarzer page 67
5	Grey, blue	Plastic	160		Em-tec page 32 HME page 38
5	Grey, blue	Plastic	173		Levitronix page 42
5	Grey, yellow	Plastic	186		Medos page 46
5	Black	Plastic	165		Schiller page 65
5	Grey	Plastic	223		Nihon Kohden page 51
5	Beige	Plastic	89		Dräger page 26
5	Grey	Plastic	195		St. Jude Medical page 74
5	White	Plastic	217		RDT (Remote Diagnostic Technologies) page 61
6	Grey, beige	Plastic	49		Invivo page 40 Mennen page 48 Philips page 56 Riegel Medical page 62



Pins	Colour	Material	Plug type	Picture	Usage
6	Grey, beige	Plastic	8		ADInstruments page 14 Penlon page 54 Dräger page 29 Propaq page 57 Edan Instruments page 31 Protocol page 58 Hoyer page 39 Saadat page 63 Ivy Biomedical page 40 Siare page 67 Mediana page 44 Spacelabs page 72 Medical Data Electr. .. page 45 Tektronix page 75 Mindray page 48 Zoll page 78
6	Grey	Plastic	204		
6	Grey, beige	Plastic	56		Aries Medical page 17 Medtronic page 47 Physio-Control page 57
6	Grey, beige	Plastic	127		BCI page 19 Dräger..... page 30 Bionet page 20 Econet Medical page 30 Criticon page 23 MR-Equipment page 51
6	Grey	Plastic	11		American Optical page 16 Colin page 21 Datascope page 24 Guardian page 36
6	Silver	Metal	187		Nippon Colin page 52 Roche Medical page 62
6	Grey, red	Plastic	162		MIPM page 50
7	Orange, beige	Plastic	230		Dräger page 27 Siemens page 70
7	Grey, red	Plastic	93		Abbott Laboratories page 14 St. Jude Medical page 74
8	Silver	Metal	192		ADInstruments page 15
8	Grey, yellow	Plastic	171		Raumedic page 61
8 (5)	Black	Plastic	147		Bard/Boston Scientific page 19
8 (5)	Black	Plastic	150		Prucka page 60 Vytech page 76
8 (5)	Black	Plastic	166		Witt page 77
8 (7)	Black	Plastic	67		Marquette page 43
10	Grey, red	Plastic	156		Spacelabs page 73



Plug types overview



Pins	Colour	Material	Plug type	Picture	Usage
10	Grey, red	Plastic	194		Siemens page 71
10	Orange, beige	Plastic	213		Dräger page 26 Siecov page 68 Siemens page 69
10	Grey	Plastic	172		GS Elektromedizinische Geräte page 35
10	White	Plastic	98		Dräger page 28
10 (6)	Grey, red	Plastic	100		Datex page 25 S&W (Simonsen & Weel) page 71
10 (6)	White, orange	Plastic	121		Hellige page 37 Litton page 43
11	Grey, red	Plastic	124		GE (General Electric) page 34 Ohmeda..... page 53 Hellige page 37 Prucka page 59 Marquette page 44
11	Grey, red	Plastic	209		GE (General Electric) page 35
11	Grey, red	Plastic	208		Datex page 25 GE (General Electric) page 34
11	Grey, blue	Plastic	94		Fukuda Denshi page 33
12	Grey, red	Plastic	48		Agilent Technologies page 17 Philips..... page 55 Axcent Medical page 18 Pulsion page 60 Comen page 23 Schwarzer page 66 HP (Hewlett-Packard) page 38 Welch Allyn page 77
12	Grey, red	Plastic	201		Philips page 56
12	Grey, red	Plastic	222		Mindray page 49 Penlon page 55
12	Grey, violet	Plastic	197		Mindray page 50
12	Grey, violet	Plastic	214		Fazzini page 32
12	Black, orange	Plastic	174		Fukuda Denshi page 33



Pins	Colour	Material	Plug type	Picture	Usage
12 (4)	Orange, beige	Plastic	112		Bruker page 21 Schiller page 65
12 (8)	Orange, beige	Plastic	73		Arrow page 18 Kontron page 42
14	Yellow, beige	Plastic	198		Nihon Kohden page 52
15	Silver	Metal	211		Dräger page 27 Siemens page 69
16	Orange, beige	Plastic	190		Dräger page 28 Siemens page 70

Products for dedicated applications



Abbott Laboratories					Monitor model(s): RadiAnalyzer Xpress, Wi-Box			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2142	BDP293RD1	50 cm	1	4046765071110	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1004	BDC901-293	320 cm	1	4046765087999	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	-	-	-	-	-	
Plug type: 93		Pins: 7, male		Material: Plastic		Additional information: -		



ADInstruments (1.1)					Monitor model(s): ML117			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



ADInstruments (1.2)					Monitor model(s): ML117			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



ADInstruments (2)					Monitor model(s): FE221 Bridge Amp, FE224 Quad Bridge Amp, FE228 Octal Bridge Amp, ML110 Bridge Amp, ML112 Quad Bridge Amp, ML221 Bridge Amp, ML224 Quad Bridge Amp, ML228 Octal Bridge Amp			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3220	BDPT-2019	350 cm	1	4046765088231	
		2	-	-	-	-	-	
Plug type: 192		Pins: 8, male		Material: Metal		Additional information: -		

Products for dedicated applications



American Optical (1.1)					Monitor model(s): 1006, 1008, 261006, 261008				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387		
		2	75.2024	BDP112DS2	50 cm	1	4046765026394		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451		
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333		
		2	-	-	-	-	-	-	
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187			



American Optical (1.2)					Monitor model(s): 1006, 1008, 261006, 261008				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681		
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667		
		2	-	-	-	-	-	-	
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11			



Agilent Technologies					Monitor model(s): M1006A, M1006B, M3000A Server			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		



Aries Medical					Monitor model(s): 700			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2038	BDP167PC1	50 cm	1	4046765026899	
		2	75.2039	BDP167PC2	50 cm	1	4046765026905	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1020	BDC901-167	320 cm	1	4046765024512	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3108	BDPT-1067	350 cm	1	4046765029777	
		2	-	-	-	-	-	
Plug type: 56		Pins: 6, male		Material: Plastic		Additional information: -		

Products for dedicated applications



Arrow					Monitor model(s): 7000 IABP, ACAT 1, ACAT 1 Plus, AutoCAT 2, AutoCAT 2 WAVE, K2000, KAAT, KAAT II, KAAT II Plus			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1021	BDC901-186	320 cm	1	4046765024543	
		2	75.1192	BDC902-186	370 cm	1	4046765105259	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3114	BDPT-1086	350 cm	1	4046765030056	
		2	-	-	-	-	-	
Plug type: 73		Pins: 8, male		Material: Plastic		Additional information: Unused pins not assembled		



Axcent Medical					Monitor model(s): CETUS-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		



Bard/Boston Scientific					Monitor model(s): CLEARSIGN Amplifier, CLEARSIGN II Amplifier, STAMP Amplifier, LabSystem 24 EP			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2190	BDPT558BA1	50 cm	1	4046765028602	
		2	75.2145	BDPT558BA2	50 cm	1	4046765028619	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1032	BDC901-558	320 cm	1	4046765067786	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3192	BDPT-5058	350 cm	1	4046765031268	
		2	-	-	-	-	-	
Plug type: 147		Pins: 5, male		Material: Plastic		Additional information: Unused pins not assembled		



BCI					Monitor model(s): 6200, 9100, Advisor			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	
Plug type: 127		Pins: 6, male		Material: Plastic		Additional information: -		

Products for dedicated applications



BLT (Biolight)					Monitor model(s): Anyview-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1248	BD901-209	320 cm	1	4046765111397	
		2	75.1249	BD902-209	370 cm	1	4046765111380	
		3	75.1250	BD903-209	370 cm	1	4046765111373	
		4	75.1251	BD904-209	370 cm	1	4046765111366	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3244	BDPT-2009	350 cm	1	4046765111090	
		2	-	-	-	-	-	
Plug type: 210		Pins: 4, male		Material: Plastic		Additional information: -		



Bionet					Monitor model(s): BM5, Compact-7			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	
Plug type: 127		Pins: 6, male		Material: Plastic		Additional information: -		



Bruker					Monitor model(s): DG2000P IBP MODULE, PHYSIOGARD, TM 910, SM785P, SM785PS, SM786P2, SM786P2S, SM786P2ST, SM786P2T, SM786PN, SM786PNS, SM786PNST, SM786PNT			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2066	BDP520BR1	50 cm	1	4046765027797	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3169	BDPT-5020	350 cm	1	4046765030711	
		2	-	-	-	-	-	
Plug type: 112		Pins: 4, male		Material: Plastic		Additional information: Unused pins not assembled		



Colin (1.1)					Monitor model(s): BP-408, CBM-3000, PILOT			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387	
		2	75.2024	BDP112DS2	50 cm	1	4046765026394	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333	
		2	-	-	-	-	-	
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187		

Products for dedicated applications




Colin (1.2)					Monitor model(s): BP-408, CBM-3000, PILOT			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667	
		2	-	-	-	-	-	
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11		




Colin (2)					Monitor model(s): BP-608EV			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2087	BDP563SR2	50 cm	1	4046765028718	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3150	BDPT-5063	350 cm	1	4046765031367	
		2	-	-	-	-	-	
Plug type: 152		Pins: 5, male		Material: Plastic		Additional information: -		




Comen					Monitor model(s): C-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		




Criticon					Monitor model(s): Dianamap Plus			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	
Plug type: 127		Pins: 6, male		Material: Plastic		Additional information: -		

Products for dedicated applications




Datascope (1.1)					Monitor model(s): Please see appendix for monitor models			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387	 
		2	75.2024	BDP112DS2	50 cm	1	4046765026394	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333	
		2	-	-	-	-	-	
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187		




Datascope (1.2)					Monitor model(s): Please see appendix for monitor models			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667	
		2	-	-	-	-	-	
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11		




Datex (1)					Monitor model(s): CARDIOCAP-series, CARDIOCAP II-series, AS/3 Light, S/5 Light, M-Modules (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2059	BDP506DX1	50 cm	1	4046765027551	 
		2	75.2060	BDP506DX2	50 cm	1	4046765027575	
		3	75.2061	BDP506DX3	50 cm	1	4046765027605	
		4	75.2120	BDP506DX4	50 cm	1	4046765027629	
	Contactboard-GO	1	75.1024	BDC901-5006	320 cm	1	4046765024574	
		2	75.1056	BDC902-506	370 cm	1	4046765025250	
		3	-	-	-	-	-	
		4	75.1103	BDC904-506	370 cm	1	4046765070427	
Contactboard-GO - FX	1	75.1037	BDC901FXD	320 cm	1	4046765072056		
DPT-6000	Interface cable	1	75.3260	BDPT-5006	320 cm	1	4046765432799	
		2	-	-	-	-	-	
Plug type: 100		Pins: 6, female		Material: Plastic		Additional information: Unused pins not assembled		




Datex (2)					Monitor model(s): S5/FM, F-FML, FM light, models from 01/05 via E-PRESTN module, I4 IMM MRI, AESTIVA, Avance			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 
		2	75.2172	BDP583DX2SM	50 cm	1	4046765076870	
		3	75.2173	BDP583DX3SM	50 cm	1	4046765076887	
		4	75.2174	BDP583DX4SM	50 cm	1	4046765076894	
	Contactboard-GO	1	75.1141	BDC901-583SM	320 cm	1	4046765076979	
		2	75.1143	BDC902-583SM	370 cm	1	4046765076993	
		3	75.1197	BDC903-583	370 cm	1	4046765105679	
		4	75.1198	BDC904-583	370 cm	1	4046765105686	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3187	BDPT-5083SM	350 cm	1	4046765072919	
		2	-	-	-	-	-	
Plug type: 208		Pins: 11, male		Material: Plastic		Additional information: Double pin-out		

Products for dedicated applications



Dräger (1)					Monitor model(s): Infinity-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2206	BDP174SI1-SM	50 cm	1	4046765108816	 
		2	75.2207	BDP174SI2-SM	50 cm	1	4046765108823	
		3	75.2208	BDP174SI3-SM	50 cm	1	4046765108830	
		4	75.2209	BDP174SI4-SM	50 cm	1	4046765108847	
	Contactboard-GO	1	75.1230	BDC901-1074-SM	320 cm	1	4046765108861	
		2	75.1231	BDC902-174-SM	370 cm	1	4046765108865	
		3	75.1246	BDC903-174-SM	370 cm	1	4046765108892	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1235	BDC901FXS-SM	320 cm	1	4046765108878		
DPT-6000	Interface cable	1	75.3236	BDPT-1074-SM	350 cm	1	4046765108809	
		2	-	-	-	-	-	
Plug type: 213		Pins: 10, male		Material: Plastic		Additional information: Monitor model(s) with integrated 10 pin block or separate 10 pin adapters needed for application (see appendix)		



Dräger (2)					Monitor model(s): Dual Hemo, Quad Hemo				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 	
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	75.2050	BDP253DR4	50 cm	1	4046765067625		
	Contactboard-GO	1	75.1149	BDC901-253A2	50 cm	1	4046765079543		
		2	75.1009	BDC902-253	35 cm	1	4046765088057		
		3	75.1122	BDC903-253	35 cm	1	4046765072315		
		4	75.1005	BDC904-253	35 cm	1	4046765088019		
Contactboard-GO - FX	1	75.1193	BDC901FX-253M2	200 cm	1	4046765105440			
DPT-6000	Interface cable	1	75.3119	BDPT-2053	25 cm	1	4046765030155		
		2	-	-	-	-	-		
Plug type: 89		Pins: 5, female		Material: Plastic		Additional information: -			



Dräger (3)					Monitor model(s): Hemo 2, Hemo 4			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 
		2	75.2219	BDP553SI2-M	35 cm	1	4046765111533	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1259	BDC901-553-M	35 cm	1	4046765111540	
		2	75.1260	BDC902-553-M	35 cm	1	4046765111557	
		3	-	-	-	1	-	
		4	75.1256	BDC904-553-M	35 cm	1	4046765111564	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	75.3143	BDPT-5053	50 cm	1	4046765031145	
Plug type: 211		Pins: 15, male		Material: Metal		Additional information: Double pin-out		



Dräger (4)					Monitor model(s): SC 5000, SC 6000P, SC 6002, SC 6002 XL, SC 6802 XL			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2256	BDP547SI1-SM	50 cm	1	4046765433949	 
		2	75.2257	BDP547SI2-SM	50 cm	1	4046765433963	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1311	BDC901-547SM	350 cm	1	4046765433987	
		2	75.1312	BDC902-547SM	370 cm	1	4046765434007	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3261	BDPT-5047SM	350 cm	1	4046765433925	
		2	-	-	-	-	-	
Plug type: 230		Pins: 7, male		Material: Plastic		Additional information: Successor of the type 137		

Products for dedicated applications


Dräger (5)					Monitor model(s): 57 00 102, 57 00 202, 57 00 602, Dialog 2000, PB 8800 Parameter-Box				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2056	BDP504DR1	50 cm	1	4046765027483		
		2	75.2057	BDP504DR2	50 cm	1	4046765027490		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	-	
		2	-	-	-	-	-	-	
Plug type: 98		Pins: 10, male		Material: Plastic		Additional information: -			

Dräger (6)					Monitor model(s): Infinity Delta, Infinity Delta XL, Infinity Gamma X XL, Infinity Kappa, Infinity Kappa XLT, Infinity Vista XL				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	75.2245	BDP596DR2	50 cm	1	4046765429850		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	75.1001	BDC902-596	350 cm	1	4046765087975		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	-	
		2	75.3161	BDPT5096C2	350 cm	1	4046765068394		
Plug type: 190		Pins: 16, male		Material: Plastic		Additional information: Double pin-out			


Dräger (7.1)					Monitor model(s): Vitalert 2000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		


Dräger (7.2)					Monitor model(s): Vitalert 2000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		

Products for dedicated applications


Dräger (8)					Monitor model(s): Vista 120, Vista 120S			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	



Plug type: 127 **Pins:** 6, male **Material:** Plastic **Additional information:** -





Econet Medical					Monitor model(s): COMPACT 5, COMPACT 7, COMPACT 9			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	


Plug type: 127 **Pins:** 6, male **Material:** Plastic **Additional information:** -




Edan Instruments (1.1)					Monitor model(s): iM50, iM60, iM70, iM8, iM8 VET, iM80, iM8A, iM8B VET, M50, M80			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		


Edan Instruments (1.2)					Monitor model(s): iM50, iM60, iM70, iM8, iM8 VET, iM80, iM8A, iM8B VET, M50, M80			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
	2	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		

Products for dedicated applications


Em-tec					Monitor model(s): SonoTT FlowComputer			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3226	BDPT-5071-220	220 cm	1	4046765103323	
		2	-	-	-	-	-	



Plug type: 160 **Pins:** 5, male **Material:** Plastic **Additional information:** -





Fazzini					Monitor model(s): PMS-series			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3250	BDPT-FAZZINI	350 cm	1	4046765428648	
		2	-	-	-	-	-	




Plug type: 214 **Pins:** 12, male **Material:** Plastic **Additional information:** -









Fukuda Denshi (1)					Monitor model(s): HM-800, HS-8000, HS-8312, DS-8007, DS-8100, DS-8200, DS-8400, DS-8500			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2052	BDP297FD1	50 cm	1	4046765067861	
		2	75.2053	BDP297FD2	50 cm	1	4046765067878	
		3	75.2054	BDP297FD3	50 cm	1	4046765068745	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3124	BDPT-2097	350 cm	1	4046765067854	
		2	-	-	-	-	-	
Plug type: 94		Pins: 11, male		Material: Plastic		Additional information: -		




Fukuda Denshi (2)					Monitor model(s): HM-800, DS-5100E, DS-7100, DS-7200, HB300, HB301, HB-60, HS-500, HS-710, HS-710E, HS-720, HS-720E, HU-71, HU-72, HU-73			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2101	BDP597FD1	50 cm	1	4046765029210	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3162	BDPT-5097	350 cm	1	4046765031695	
		2	-	-	-	-	-	
Plug type: 174		Pins: 12, male		Material: Plastic		Additional information: For HS-710 to HS-720E branch cable with REF: CJ-7546 (Fukuda Denshi product) needed for application		

Products for dedicated applications



GE (General Electric) (1)					Monitor model(s): Carescape B450, B650, B850, AESPIRE, Aisys Carestation, SOLAR 8000i			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2148	BDP533MQ1SM	50 cm	1	4046765071813	 
		2	75.2149	BDP533MQ2SM	50 cm	1	4046765071820	
		3	75.2150	BDP533MQ3SM	50 cm	1	4046765071837	
		4	75.2151	BDP533MQ4SM	50 cm	1	4046765071844	
	Contactboard-GO	1	75.1127	BDC901-533SM	320 cm	1	4046765072926	
		2	75.1130	BDC902-533SM	370 cm	1	4046765073107	
		3	75.1131	BDC903-533SM	370 cm	1	4046765073084	
		4	75.1180	BDC904-533SM	370 cm	1	4046765102838	
Contactboard-GO - FX	1	75.1132	BDC901FXM	320 cm	1	4046765073091		
DPT-6000	Interface cable	1	75.3186	BDPT-5033SM	350 cm	1	4046765071851	
		2	-	-	-	-	-	
Plug type: 124		Pins: 11, male		Material: Plastic		Additional information: Single pin-out; Adapter needed for double pin-out sockets		



GE (General Electric) (2)					Monitor model(s): DASH iMM, i4, B40, E-PRESTN			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 
		2	75.2172	BDP583DX2SM	50 cm	1	4046765076870	
		3	75.2173	BDP583DX3SM	50 cm	1	4046765076887	
		4	75.2174	BDP583DX4SM	50 cm	1	4046765076894	
	Contactboard-GO	1	75.1141	BDC901-583SM	320 cm	1	4046765076979	
		2	75.1143	BDC902-583SM	370 cm	1	4046765076993	
		3	75.1197	BDC903-583	370 cm	1	4046765105679	
		4	75.1198	BDC904-583	370 cm	1	4046765105686	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3187	BDPT-5083SM	350 cm	1	4046765072919	
		2	-	-	-	-	-	
Plug type: 208		Pins: 11, male		Material: Plastic		Additional information: Double pin-out		




GE (General Electric) (3)					Monitor model(s): PDM (Patient Data Module)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2175	BDP283GE1SM	50 cm	1	4046765076900	 
		2	75.2176	BDP283GE2SM	50 cm	1	4046765076917	
		3	75.2177	BDP283GE3SM	50 cm	1	4046765076924	
		4	75.2178	BDP283GE4SM	50 cm	1	4046765076931	
	Contactboard-GO	1	75.1186	BDC901-283	320 cm	1	4046765104542	
		2	75.1187	BDC902-283	370 cm	1	4046765104559	
		3	75.1181	BDC903-283	370 cm	1	4046765102821	
		4	75.1188	BDC904-283	370 cm	1	4046765104566	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3229	BDPT-2083	350 cm	1	4046765104535	
		2	-	-	-	-	-	
Plug type: 209		Pins: 11, male		Material: Plastic		Additional information: Double pin-out and configuration for GE PDM (Patient Data Module)		




GS Elektromedizinische Geräte					Monitor model(s): Corpuls 3			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2100	BDP589CP1	50 cm	1	4046765029180	 
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	75.1063	BDC902-589	320 cm	1	4046765068714	
		3	-	-	-	-	-	
		4	75.1263	BDC904-589-370	370 cm	1	4046765111632	
Contactboard-GO - FX	1	75.1270	BDC901FXCP	320 cm	1	4046765428242		
DPT-6000	Interface cable	1	75.3158	BDPT-5089	350 cm	1	4046765031633	
		2	75.3249	BDPT-5089C2	350 cm	1	4046765428235	
Plug type: 172		Pins: 10, male		Material: Plastic		Additional information: Double pin-out; Adapter with ODU-socket needed		

Products for dedicated applications




Guardian (1.1)					Monitor model(s): -				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387		
		2	75.2024	BDP112DS2	50 cm	1	4046765026394		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451		
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333		
		2	-	-	-	-	-		
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187			




Guardian (1.2)					Monitor model(s): -				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681		
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667		
		2	-	-	-	-	-		
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11			



Hellige (1)					Monitor model(s): DASH 2000 - 5000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2148	BDP533MQ1SM	50 cm	1	4046765071813	 
		2	75.2149	BDP533MQ2SM	50 cm	1	4046765071820	
		3	75.2150	BDP533MQ3SM	50 cm	1	4046765071837	
		4	75.2151	BDP533MQ4SM	50 cm	1	4046765071844	
	Contactboard-GO	1	75.1127	BDC901-533SM	320 cm	1	4046765072926	
		2	75.1130	BDC902-533SM	370 cm	1	4046765073107	
		3	75.1131	BDC903-533SM	370 cm	1	4046765073084	
		4	75.1180	BDC904-533SM	370 cm	1	4046765102838	
Contactboard-GO - FX	1	75.1132	BDC901FXM	320 cm	1	4046765073091		
DPT-6000	Interface cable	1	75.3186	BDPT-5033SM	350 cm	1	4046765071851	
		2	-	-	-	-	-	
Plug type: 124		Pins: 11, male		Material: Plastic		Additional information: -		



Hellige (2)					Monitor model(s): MA88, MS-series, SERVOMED-series, SM-series, SMC 108, SMS-series, SMV-series, 206 010 xx, 227 064 xx, 236 015 xx, 236 016 xx, 236 037 xx, 236 060 xx				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 	
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3135	BDPT-5029	350 cm	1	4046765030773		
		2	-	-	-	-	-		
Plug type: 121		Pins: 6, male		Material: Plastic		Additional information: Unused pins not assembled			

Products for dedicated applications



HP (Hewlett-Packard)					Monitor model(s): 78205 (A-D) Opt. 006, 78341A, 78342A, 78345A, 78346A, 78353A, 78353B, 78354A, 78532A, 78534A, 78551AX, 78551D, 78553A, 78801A, 78801B, 78833A, 78834, 8805 (A-D), M1006A, M1006A, M1006B, M1280, M3000A			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	 
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		



HME					Monitor model(s): LifePulse-series				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	 	
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-			
DPT-6000	Interface cable	1	75.3226	BDPT-5071-220	220 cm	1	4046765103323		
		2	-	-	-	-	-		
Plug type: 160		Pins: 5, male		Material: Plastic		Additional information: -			



Hoyer (1.1)					Monitor model(s): -			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



Hoyer (1.2)					Monitor model(s): -			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		

Products for dedicated applications



Invivo					Monitor model(s): 3000, 3202, EXPRESSION MR 200 & 400, M12 3550, Magnitude 3150M, Millennium 3500, Omni-Trak 3150 MRI, Omni-Trak TVS, Precess 3160				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2036	BDP157ME1	50 cm	1	4046765026844		
		2	75.2037	BDP157ME2	50 cm	1	4046765026868		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1123	BDC901-1057	320 cm	1	4046765024413		
		2	75.1191	BDC902-157	370 cm	1	4046765105044		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3105	BDPT-1057	350 cm	1	4046765029739		
		2	-	-	-	-	-	-	
Plug type: 49		Pins: 6, male		Material: Plastic		Additional information: -			



Ivy Biomedical (1.1)					Monitor model(s): 401, 402, 403, 403D, 404, 405C, 405D, 405T, 700, 701, 704, 706, VITAL-GUARD 450C				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431		
		2	75.2026	BDP114SL2	50 cm	1	4046765026462		
		3	75.2027	BDP114SL3	50 cm	1	4046765026493		
		4	75.2164	BDP114SL4	50 cm	1	4046765026516		
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468		
		2	75.1212	BDC902-114	370 cm	1	4046765106089		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998	-		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340		
		2	-	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204			



Ivy Biomedical (1.2)					Monitor model(s): 401, 402, 403, 403D, 404, 405C, 405D, 405T, 700, 701, 704, 706, VITAL-GUARD 450C			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Jostra					Monitor model(s): DPM 20-440, DPM 20-441			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3171	BDPT-5086	350 cm	1	4046765031626	
		2	-	-	-	-	-	
Plug type: 168		Pins: 5, male		Material: Metal		Additional information: -		

Products for dedicated applications




Kontron					Monitor model(s): 6700-series, 7000-series, Fetalogic 4000, K2000, KAAT, KAAT II, Minimon-series, Supermon-series, KOLORMON-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1021	BDC901-186	320 cm	1	4046765024543	
		2	75.1192	BDC902-186	370 cm	1	4046765105259	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3114	BDPT-1086	350 cm	1	4046765030056	
		2	-	-	-	-	-	
Plug type: 73		Pins: 8, male		Material: Plastic		Additional information: Unused pins not assembled		




Levitronix					Monitor model(s): CentriMag, B-Pump			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	75.1173	BDC902-590	370 cm	1	4046765025328	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3175	BDPT-5090	350 cm	1	4046765031640	
		2	-	-	-	-	-	
Plug type: 173		Pins: 5, male		Material: Plastic		Additional information: -		



Litton					Monitor model(s): 236-060-XX, FM660				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3135	BDPT-5029	350 cm	1	4046765030773		
		2	-	-	-	-	-	-	
Plug type: 121		Pins: 6, male		Material: Plastic		Additional information: Unused pins not assembled			



Marquette (1)					Monitor model(s): 7000- & 7010-series				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3113	BDPT-1079	350 cm	1	4046765030018		
		2	-	-	-	-	-	-	
Plug type: 67		Pins: 7, male		Material: Plastic		Additional information: Unused pins not assembled			

Products for dedicated applications



Marquette (2)					Monitor model(s): MonoEagle, Eagle 2XXX-4XXX, CathLab, TRAM 600 Tramscope, Tramrac, Solar 7000-9000, Mac LAB Tram RAC			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2148	BDP533MQ1SM	50 cm	1	4046765071813	 
		2	75.2149	BDP533MQ2SM	50 cm	1	4046765071820	
		3	75.2150	BDP533MQ3SM	50 cm	1	4046765071837	
		4	75.2151	BDP533MQ4SM	50 cm	1	4046765071844	
	Contactboard-GO	1	75.1127	BDC901-533SM	320 cm	1	4046765072926	
		2	75.1130	BDC902-533SM	370 cm	1	4046765073107	
		3	75.1131	BDC903-533SM	370 cm	1	4046765073084	
		4	75.1180	BDC904-533SM	370 cm	1	4046765102838	
Contactboard-GO - FX	1	75.1132	BDC901FXM	320 cm	1	4046765073091		
DPT-6000	Interface cable	1	75.3186	BDPT-5033SM	350 cm	1	4046765071851	
		2	-	-	-	-	-	
Plug type: 124		Pins: 11, male		Material: Plastic		Additional information: -		



Mediana (1.1)					Monitor model(s): D500, M50, YM6000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	 
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



Mediana (1.2)					Monitor model(s): D500, M50, YM6000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Medical Data Electronics (1.1)					Monitor model(s): Escort 100, Escort 100T, Escort 101, Escort 101T, Escort 102, Escort 102T, Escort 300, Escort 300A, Escort PRISM SE, MPM			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		

Products for dedicated applications



Medical Data Electronics (1.2)					Monitor model(s): Escort 100, Escort 100T, Escort 101, Escort 101T, Escort 102, Escort 102T, Escort 300, Escort 300A, Escort PRISM SE, MPM				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069		
		2	75.2216	-	50 cm	1	4046765111052		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632			
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083		
		2	-	-	-	-	-		
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space			



Medos					Monitor model(s): Delta Stream				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	75.1134	BDC903-203-90	90 cm	1	4046765074845		
		4	75.1195	BDC904-203-90	90 cm	1	4046765105563		
Contactboard-GO - FX	1	-	-	-	-	-			
DPT-6000	Interface cable	1	75.3168	MEDP20-015	25 cm	1	4046765053123		
		2	-	-	-	-	-		
Plug type: 186		Pins: 5, male		Material: Plastic		Additional information: -			



Medrad					Monitor model(s): Veris 8600 MRI				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	75.1194	BDC902-584	370 cm	1	4046765105600		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3157	BDPT-5084	350 cm	1	4046765031619		
		2	-	-	-	-	-		
Plug type: 86		Pins: 5, male		Material: Metal		Additional information: -			



Medtronic					Monitor model(s): LIFEPAK 12, LIFEPAK 15				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2038	BDP167PC1	50 cm	1	4046765026899		
		2	75.2039	BDP167PC2	50 cm	1	4046765026905		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	75.1020	BDC901-167	320 cm	1	4046765024512		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3108	BDPT-1067	350 cm	1	4046765029777		
		2	-	-	-	-	-		
Plug type: 56		Pins: 6, male		Material: Plastic		Additional information: -			

Products for dedicated applications


Mennen					Monitor model(s): 2110, CFE, Envoy-series, HORIZON-series, MERCURY, MR-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2036	BDP157ME1	50 cm	1	4046765026844	
		2	75.2037	BDP157ME2	50 cm	1	4046765026868	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1123	BDC901-1057	320 cm	1	4046765024413	
		2	75.1191	BDC902-157	370 cm	1	4046765105044	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3105	BDPT-1057	350 cm	1	4046765029739	
		2	-	-	-	-	-	
Plug type: 49		Pins: 6, male		Material: Plastic		Additional information: For CFE and Envoy-series 8-6 pin adapter needed for application		

Mindray (1.1)					Monitor model(s): PM-series, MEC-1000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		


Mindray (1.2)					Monitor model(s): PM-series, MEC-1000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		


Mindray (2)					Monitor model(s): PM6000, BeneView T5/T6/T8/T9			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2247	BDP156HP1-RSM	50 cm	1	4046765432300	
		2	75.2246	BDP156HP2-RSM	50 cm	1	4046765432317	
		3	75.2248	BDP156HP3-RSM	50 cm	1	4046765432324	
		4	75.2249	BDP156HP4-RSM	50 cm	1	4046765432331	
	Contactboard-GO	1	75.1281	BDC901-156RSM	350 cm	1	4046765432287	
		2	75.1282	BDC902-156RSM	370 cm	1	4046765432294	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3255	BDPT-1056RSM	350 cm	1	4046765432270	
		2	-	-	-	-	-	
Plug type: 222		Pins: 12, male		Material: Plastic		Additional information: Successor of the type 180		

Products for dedicated applications


Mindray (3)					Monitor model(s): IM2204, BeneView T1			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2189	BDP215MY2	50 cm	1	4046765080013	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1176	BDC901-215	350 cm	1	4046765101916	
		2	75.1177	BDC902-215	350 cm	1	4046765101978	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	-	-	-	-	-	



Plug type: 197
Pins: 12, male
Material: Plastic
Additional information: Double pin-out





MIPM					Monitor model(s): Tesla M3, Tesla DUO			Picture
Product type	Signals	REF	Historical no.	Length		GTIN		
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1190	BDC901-576	320 cm	1	4046765104528	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3227	BDPT-5076	350 cm	1	4046765104276	
		2	-	-	-	-	-	



Plug type: 162
Pins: 6, male
Material: Plastic
Additional information: -







MR-Equipment					Monitor model(s): 9500			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	75.2238	BDP537DR3	50 cm	1	4046765429522	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1271	BDC901-537	320 cm	1	4046765428259	
		2	-	-	-	-	-	
		3	75.1273	BDC903-537	370 cm	1	4046765429515	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3138	BDPT-5037	350 cm	1	4046765030964	
		2	-	-	-	-	-	
Plug type: 127		Pins: 6, male		Material: Plastic		Additional information: -		



Nihon Kohden (1)					Monitor model(s): AP-800PA, Life Scope-series, AA-series, AY-series, BSM-series, SVM-series, Vismo-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2250	BDP178NK1SM	50 cm	1	4046765432393	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1289	BDC901-178SM	350 cm	1	4046765432355	
		2	75.1290	BDC902-178SM	370 cm	1	4046765432362	
		3	75.1291	BDC903-178SM	370 cm	1	4046765432379	
		4	75.1292	BDC904-178SM	370 cm	1	4046765432386	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3256	BDPT-1078SM	350 cm	1	4046765432348	
		2	-	-	-	-	-	
Plug type: 223		Pins: 5, male		Material: Plastic		Additional information: Successor of the type 66. For AA-series, AY-series, BSM-series, Vismo-series REF: JP-910P (Nihon Kohden product) needed for application. For SVM-series REF: JP-751P (Nihon Kohden product) needed for application.		

Products for dedicated applications



Nihon Kohden (2)					Monitor model(s): Life Scope-series, AA-series, AY-series, BSM-series, Vismo-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2234	BDP559NK1-SM	50 cm	1	4046765428266	
		2	75.2233	BDP559NK2-SM	50 cm	1	4046765428273	
		3	75.2232	BDP559NK3-SM	50 cm	1	4046765428280	
		4	75.2235	BDP559NK4-SM	50 cm	1	4046765428297	
	Contactboard-GO	1	75.1214	BDC901-559NK	320 cm	1	4046765106218	
		2	75.1215	BDC902-559NK	370 cm	1	4046765106225	
		3	75.1216	BDC903-559NK	370 cm	1	4046765106201	
		4	75.1217	BDC904-559NK	370 cm	1	4046765106232	
Contactboard-GO - FX	1	75.1272	BDC901FXNK-SM	320 cm	1	4046765428303		
DPT-6000	Interface cable	1	75.3233	BDPT-5059NK	350 cm	1	4046765106249	
		2	-	-	-	-	-	
Plug type: 198		Pins: 14, male		Material: Plastic		Additional information: -		



Nippon Colin (1.1)					Monitor model(s): 308ETI, CBM300D			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387	
		2	75.2024	BDP112DS2	50 cm	1	4046765026394	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333	
		2	-	-	-	-	-	
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187		




Nippon Colin (1.2)					Monitor model(s): 308ETI, CBM300D			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667	
		2	-	-	-	-	-	
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11		




Ohmeda					Monitor model(s): FM-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2148	BDP533MQ1SM	50 cm	1	4046765071813	
		2	75.2149	BDP533MQ2SM	50 cm	1	4046765071820	
		3	75.2150	BDP533MQ3SM	50 cm	1	4046765071837	
		4	75.2151	BDP533MQ4SM	50 cm	1	4046765071844	
	Contactboard-GO	1	75.1127	BDC901-533SM	320 cm	1	4046765072926	
		2	75.1130	BDC902-533SM	370 cm	1	4046765073107	
		3	75.1131	BDC903-533SM	370 cm	1	4046765073084	
		4	75.1180	BDC904-533SM	370 cm	1	4046765102838	
Contactboard-GO - FX	1	75.1132	BDC901FXM	320 cm	1	4046765073091		
DPT-6000	Interface cable	1	75.3186	BDPT-5033SM	350 cm	1	4046765071851	
		2	-	-	-	-	-	
Plug type: 124		Pins: 11, male		Material: Plastic		Additional information: -		

Products for dedicated applications



Penlon (1.1)					Monitor model(s): SP-series, PM-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
	Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998	
2		-	-	-	-	-		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



Penlon (1.2)					Monitor model(s): SP-series, PM-series			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
	Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632	
2		-	-	-	-	-		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Penlon (2)					Monitor model(s): PM6000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2247	BDP156HP1-RSM	50 cm	1	4046765432300	 
		2	75.2246	BDP156HP2-RSM	50 cm	1	4046765432317	
		3	75.2248	BDP156HP3-RSM	50 cm	1	4046765432324	
		4	75.2249	BDP156HP4-RSM	50 cm	1	4046765432331	
	Contactboard-GO	1	75.1281	BDC901-156RSM	350 cm	1	4046765432287	
		2	75.1282	BDC902-156RSM	370 cm	1	4046765432294	
		3	-	-	-	-	-	
	Contactboard-GO - FX	1	-	-	-	-	-	
2		-	-	-	-	-		
DPT-6000	Interface cable	1	75.3255	BDPT-1056RSM	350 cm	1	4046765432270	
		2	-	-	-	-	-	
Plug type: 222		Pins: 12, male		Material: Plastic		Additional information: Successor of the type 180		



Philips (1)					Monitor model(s): HeartStart MRx, IntelliVue-series, MP-series, MX-series, M-modules, SureSigns-series, Capnography, Xper Flex Cardio Physiomonitring System (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	 
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
	Contactboard-GO - FX	1	75.1147	BDC904-156	370 cm	1	4046765078980	
2		75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: For multi pressure monitoring with IntelliVue X3 Y-Adapter REF: 989803199741 (Philips product) needed for application		

Products for dedicated applications



Philips (2)					Monitor model(s): IntelliVue X3			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2237	BDP156X32	50 cm	1	4046765429492	
		3	-	-	-	-	-	
		4	75.2236	BDP156X34	50 cm	1	4046765429485	
	Contactboard-GO	1	-	-	-	-	-	
		2	75.1226	BDC902-1056X3	320 cm	1	4046765107703	
		3	-	-	-	-	-	
		4	75.1227	BDC904-1056X3	370 cm	1	4046765107710	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	75.3251	BDPT-1056X3	350 cm	1	4046765429508	
Plug type: 201		Pins: 12, male		Material: Plastic		Additional information: Double pin-out; No Y-Adapter REF: 989803199741 (Philips product) needed for application		



Philips (3)					Monitor model(s): Expression MR200, Expression MR400			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2036	BDP157ME1	50 cm	1	4046765026844	
		2	75.2037	BDP157ME2	50 cm	1	4046765026868	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1123	BDC901-1057	320 cm	1	4046765024413	
		2	75.1191	BDC902-157	370 cm	1	4046765105044	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3105	BDPT-1057	350 cm	1	4046765029739	
		2	-	-	-	-	-	
Plug type: 49		Pins: 6, male		Material: Plastic		Additional information: -		



Physio-Control					Monitor model(s): LIFEPAK 12, LIFEPAK 15, VSM-1, VSM-4A, VSM-5			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2038	BDP167PC1	50 cm	1	4046765026899	
		2	75.2039	BDP167PC2	50 cm	1	4046765026905	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1020	BDC901-167	320 cm	1	4046765024512	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3108	BDPT-1067	350 cm	1	4046765029777	
		2	-	-	-	-	-	
Plug type: 56		Pins: 6, male		Material: Plastic		Additional information: -		



Propaq (1.1)					Monitor model(s): IMediserive, Encore 204EL, Encore 206EL, 104 EL, 106 EL, 104 LCD, 106 LCD, 244 CS, 246 CS			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		

Products for dedicated applications



Propaq (1.2)					Monitor model(s): IMediserve, Encore 204EL, Encore 206EL, 104 EL, 106 EL, 104 LCD, 106 LCD, 244 CS, 246 CS			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Protocol (1.1)					Monitor model(s): Encore 204EL, Encore 206EL, 104 EL, 106 EL, 104 LCD, 106 LCD, 244 CS, 246 CS			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



Protocol (1.2)					Monitor model(s): Encore 204EL, Encore 206EL, 104 EL, 106 EL, 104 LCD, 106 LCD, 244 CS, 246 CS			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Prucka (1)					Monitor model(s): MacLab 7200, GE-Innova			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2148	BDP533MQ1SM	50 cm	1	4046765071813	
		2	75.2149	BDP533MQ2SM	50 cm	1	4046765071820	
		3	75.2150	BDP533MQ3SM	50 cm	1	4046765071837	
		4	75.2151	BDP533MQ4SM	50 cm	1	4046765071844	
	Contactboard-GO	1	75.1127	BDC901-533SM	320 cm	1	4046765072926	
		2	75.1130	BDC902-533SM	370 cm	1	4046765073107	
		3	75.1131	BDC903-533SM	370 cm	1	4046765073084	
		4	75.1180	BDC904-533SM	370 cm	1	4046765102838	
Contactboard-GO - FX	1	75.1132	BDC901FXM	320 cm	1	4046765073091		
DPT-6000	Interface cable	1	75.3186	BDPT-5033SM	350 cm	1	4046765071851	
		2	-	-	-	-	-	
Plug type: 124		Pins: 11, male		Material: Plastic		Additional information: -		

Products for dedicated applications



Prucka (2)					Monitor model(s): CardioLab II, CardioLab II plus, 4000				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2084	BDP561PA1	50 cm	1	4046765028671		
		2	75.2085	BDP561PA2	50 cm	1	4046765028688		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	75.1182	BDC902-561	370 cm	1	4046765103347		
		3	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
	2	-	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-		
		2	-	-	-	-	-		
Plug type: 150		Pins: 5, male		Material: Plastic		Additional information: Unused pins not assembled			



Pulsion					Monitor model(s): PiCCO plus, PiCCO2, PulsioFlex			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
Contactboard-GO - FX	1	75.1147	BDC904-156	370 cm	1	4046765078980		
	2	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		



Raumedic					Monitor model(s): MPR 1 Datalogger, MPR2 logO Datalogger			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1136	BDC901-587-2M	200 cm	1	4046765024765	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3231	BDPT-5087-2M	200 cm	1	4046765105013	
		2	-	-	-	-	-	
Plug type: 171		Pins: 8, male		Material: Plastic		Additional information: -		



RDT (Remote Diagnostic Technologies)					Monitor model(s): Tempus Pro			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2239	BDP247TP1	50 cm	1	4046765429713	
		2	75.2240	BDP247TP2	50 cm	1	4046765429706	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1275	BDC901-247	200 cm	1	4046765429690	
		2	75.1276	BDC902-247	200 cm	1	4046765429676	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1274	BDC901FXTP	200 cm	1	4046765429683		
DPT-6000	Interface cable	1	75.3253	BDPT-2047	200 cm	1	4046765429669	
		2	-	-	-	-	-	
Plug type: 217		Pins: 5, male		Material: Plastic		Additional information: Y-Adapter REF: 01-2108 (RDT product) needed for application		

Products for dedicated applications



Riegel Medical					Monitor model(s): MultiCare 409, MultiCare 429			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2036	BDP157ME1	50 cm	1	4046765026844	
		2	75.2037	BDP157ME2	50 cm	1	4046765026868	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1123	BDC901-1057	320 cm	1	4046765024413	
		2	75.1191	BDC902-157	370 cm	1	4046765105044	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3105	BDPT-1057	350 cm	1	4046765029739	
		2	-	-	-	-	-	
Plug type: 49		Pins: 6, male		Material: Plastic		Additional information: -		



Roche Medical (1.1)					Monitor model(s): 108, 120, 126, 127, 128, 220, 222, 223, 224, 3010, 5120, 7320, 7324			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2023	BDP112DS1	50 cm	1	4046765026387	
		2	75.2024	BDP112DS2	50 cm	1	4046765026394	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1016	BDC901-112-1M	100 cm	1	4046765024451	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3100	BDPT-1012	350 cm	1	4046765029333	
		2	-	-	-	-	-	
Plug type: 11		Pins: 6, female		Material: Plastic		Additional information: Equivalent to plug type 187		



Roche Medical (1.2)					Monitor model(s): 108, 120, 126, 127, 128, 220, 222, 223, 224, 3010, 5120, 7320, 7324			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1098	BDC901-112M	320 cm	1	4046765069681	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3172	BDPT-1012M	350 cm	1	4046765069667	
		2	-	-	-	-	-	
Plug type: 187		Pins: 6, female		Material: Metal		Additional information: Equivalent to plug type 11		



Saadat (1.1)					Monitor model(s): Aria, Alborz B5, Alborz B7, NOVIN S1600, NOVIN S1800			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		

Products for dedicated applications



Saadat (1.2)					Monitor model(s): Aria, Alborz B5, Alborz B7, NOVIN S1600, NOVIN S1800				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069		
		2	75.2216	-	50 cm	1	4046765111052		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632			
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083		
		2	-	-	-	-	-		
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space			



Sarns					Monitor model(s): 9000				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3142	BDPT-5052	350 cm	1	4046765031138		
		2	-	-	-	-	-	-	
Plug type: 142		Pins: 5, male		Material: Metal		Additional information: -			



Schiller (1)					Monitor model(s): MAGLIFE C-Plus, MAGSCREEN, Defigard DG 2002, PHYSIOGARD TM 910				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2066	BDP520BR1	50 cm	1	4046765027797		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3169	BDPT-5020	350 cm	1	4046765030711		
		2	-	-	-	-	-		
Plug type: 112		Pins: 4, male		Material: Plastic		Additional information: Unused pins not assembled			



Schiller (2)					Monitor model(s): ARGUS LCM plus D, ARGUS LCM plus E				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3154	BDPT-5081	350 cm	1	4046765031565		
		2	-	-	-	-	-		
Plug type: 165		Pins: 5, male		Material: Plastic		Additional information: -			

Products for dedicated applications



Schiller (3)					Monitor model(s): ARGUS LCX, ARGUS PB-1000, ARGUS PB-2000, ARGUS PB-2200, ARGUS PRO LifeCare				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2103	BDP599SC1	50 cm	1	4046765029265		
		2	75.2104	BDP599SC2	50 cm	1	4046765029272		
		3	75.2227	BDP599SC3	50 cm	1	4046765427993		
		4	75.2231	BDP599SC4	50 cm	1	4046765428006		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3164	BDPT-5099	350 cm	1	4046765031725		
		2	-	-	-	-	-	-	
Plug type: 176		Pins: 4, male		Material: Plastic		Additional information: -			



Schwarzer (1)					Monitor model(s): -			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021	-	
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	-
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		



Schwarzer (2)					Monitor model(s): Cardis 5000			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2087	BDP563SR2	50 cm	1	4046765028718	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3150	BDPT-5063	350 cm	1	4046765031367	
		2	-	-	-	-	-	
Plug type: 152		Pins: 5, male		Material: Plastic		Additional information: -		



Siare Group (1.1)					Monitor model(s): Neptune 12, Neptune 15			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		

Products for dedicated applications



Siare Group (1.2)					Monitor model(s): Neptune 12, Neptune 15			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		



Siecov					Monitor model(s): 804/700			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2206	BDP174SI1-SM	50 cm	1	4046765108816	
		2	75.2207	BDP174SI2-SM	50 cm	1	4046765108823	
		3	75.2208	BDP174SI3-SM	50 cm	1	4046765108830	
		4	75.2209	BDP174SI4-SM	50 cm	1	4046765108847	
	Contactboard-GO	1	75.1230	BDC901-1074-SM	320 cm	1	4046765108861	
		2	75.1231	BDC902-174-SM	370 cm	1	4046765108865	
		3	75.1246	BDC903-174-SM	370 cm	1	4046765108892	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1235	BDC901FXS-SM	320 cm	1	4046765108878		
DPT-6000	Interface cable	1	75.3236	BDPT-1074-SM	350 cm	1	4046765108809	
		2	-	-	-	-	-	
Plug type: 213		Pins: 10, male		Material: Plastic		Additional information: -		



Siemens (1)					Monitor model(s): 342R, Cartridge, Hemo2 Pod, Hemo4 Pod, HemoMed Pod, SC-series, Sirecust-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2206	BDP174SI1-SM	50 cm	1	4046765108816	
		2	75.2207	BDP174SI2-SM	50 cm	1	4046765108823	
		3	75.2208	BDP174SI3-SM	50 cm	1	4046765108830	
		4	75.2209	BDP174SI4-SM	50 cm	1	4046765108847	
	Contactboard-GO	1	75.1230	BDC901-1074-SM	320 cm	1	4046765108861	
		2	75.1231	BDC902-174-SM	370 cm	1	4046765108865	
		3	75.1246	BDC903-174-SM	370 cm	1	4046765108892	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1235	BDC901FXS-SM	320 cm	1	4046765108878		
DPT-6000	Interface cable	1	75.3236	BDPT-1074-SM	350 cm	1	4046765108809	
		2	-	-	-	-	-	
Plug type: 213		Pins: 10, male		Material: Plastic		Additional information: Monitor models with integrated 10 pin block or separate 10 pin adapters needed for application (see appendix)		



Siemens (2)					Monitor model(s): Hemo 2, Hemo 4			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2219	BDP553SI2-M	35 cm	1	4046765111533	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1259	BDC901-553-M	35 cm	1	4046765111540	
		2	75.1260	BDC902-553-M	35 cm	1	4046765111557	
		3	-	-	-	1	-	
		4	75.1256	BDC904-553-M	35 cm	1	4046765111564	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	75.3143	BDPT-5053	50 cm	1	4046765031145	
Plug type: 211		Pins: 15, male		Material: Metal		Additional information: Double pin-out		

Products for dedicated applications



Siemens (3)					Monitor model(s): SC 5000, SC 6000P, SC 6002, SC 6002 XL, SC 6802 XL			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2256	BDP547SI1-SM	50 cm	1	4046765433949	
		2	75.2257	BDP547SI2-SM	50 cm	1	4046765433963	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1311	BDC901-547SM	350 cm	1	4046765433987	
		2	75.1312	BDC902-547SM	370 cm	1	4046765434007	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3261	BDPT-5047SM	350 cm	1	4046765433925	
		2	-	-	-	-	-	
Plug type: 230		Pins: 7, male		Material: Plastic		Additional information: -		



Siemens (4)					Monitor model(s): Infinity Delta, Infinity Delta XL, Infinity Gamma X XL, Infinity Kappa, Infinity Kappa XLT, Infinity Vista XL			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2245	BDP596DR2	50 cm	1	4046765429850	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	75.1001	BDC902-596	350 cm	1	4046765087975	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	75.3161	BDPT5096C2	350 cm	1	4046765068394	
Plug type: 190		Pins: 16, male		Material: Plastic		Additional information: Double pin-out		



Siemens (5)					Monitor model(s): Sensis HiSIB, Vibe			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2243	BDP221DR2	50 cm	1	4046765429799	
		3	-	-	-	-	-	
		4	75.2244	BDP221DR4	50 cm	1	4046765429805	
	Contactboard-GO	1	-	-	-	-	-	
		2	75.1174	BDC902-221-100	100 cm	1	4046765101633	
		3	-	-	-	-	-	
		4	75.1175	BDC904-221-100	100 cm	1	4046765101664	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	75.3221	BDPT2021C2-100	100 cm	1	4046765101794	
Plug type: 194		Pins: 10, male		Material: Plastic		Additional information: Double pin-out		



S&W (Simonsen & Weel)					Monitor model(s): 8041, Diascope			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2059	BDP506DX1	50 cm	1	4046765027551	
		2	75.2060	BDP506DX2	50 cm	1	4046765027575	
		3	75.2061	BDP506DX3	50 cm	1	4046765027605	
		4	75.2120	BDP506DX4	50 cm	1	4046765027629	
	Contactboard-GO	1	75.1024	BDC901-5006	320 cm	1	4046765024574	
		2	75.1056	BDC902-506	370 cm	1	4046765025250	
		3	-	-	-	-	-	
		4	75.1103	BDC904-506	370 cm	1	4046765070427	
Contactboard-GO - FX	1	75.1037	BDC901FXD	320 cm	1	4046765072056		
DPT-6000	Interface cable	1	75.3260	BDPT-5006	320 cm	1	4046765432799	
		2	-	-	-	-	-	
Plug type: 100		Pins: 6, female		Material: Plastic		Additional information: Unused pins not assembled		

Products for dedicated applications



Spacelabs (1.1)					Monitor model(s): Alpha 7, 3700, 4XX-series, 5XX-series, 7XX-series, 9XXXX-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204. Monitor models with integrated 6 pin block or separate 6 pin adapters needed for application (see appendix)		



Spacelabs (1.2)					Monitor model(s): Alpha 7, 3700, 4XX-series, 5XX-series, 7XX-series, 9XXXX-series (for further information see appendix)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space. Monitor models with integrated 6 pin block or separate 6 pin adapters needed for application (see appendix)		



Spacelabs (1.3)					Monitor model(s): Xprezzon, 90496-B, 90496-C, 91496-B Ultraview, 91496-C Ultraview			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2088	BDP567SL2	50 cm	1	4046765028787	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	-	-	-	-	-	
Plug type: 156		Pins: 10, male		Material: Plastic		Additional information: Double pin-out		



Stöckert					Monitor model(s): S5 System, SC System, 22-01-00, 22-05-00, 22-06-00, 22-07-00, 22-10-20, 22-20-20, 27-60-00, 27-70-20, 27-80-20			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2062	BDP518ST1	50 cm	1	4046765027742	
		2	75.2063	BDP518ST2	50 cm	1	4046765027759	
		3	75.2064	BDP518ST3	50 cm	1	4046765027766	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1184	BDC901-518-100	100 cm	1	4046765103941	
		2	-	-	-	-	-	
		3	75.1189	BDC903-518-1M8	180 cm	1	4046765104573	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3133	BDPT-5018	350 cm	1	4046765030681	
		2	-	-	-	-	-	
Plug type: 110		Pins: 4, male		Material: Metal		Additional information: -		

Products for dedicated applications



St. Jude Medical (1)					Monitor model(s): WorkMate Claris			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1287	BDC901-226SM	350 cm	1	4046765432447	
		2	75.1288	BDC902-226SM	370 cm	1	4046765432454	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3257	BDPT-2026SM	350 cm	1	4046765432430	
		2	-	-	-	-	-	
Plug type: 195		Pins: 5, male		Material: Plastic		Additional information: -		



St. Jude Medical (2)					Monitor model(s): RadiAnalyzer Xpress, Wi-Box			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2142	BDP293RD1	50 cm	1	4046765071110	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	75.1004	BDC901-293	320 cm	1	4046765087999	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-	
		2	-	-	-	-	-	
Plug type: 93		Pins: 7, male		Material: Plastic		Additional information: -		



Tektronix (1.1)					Monitor model(s): 413A, 414, 414-21			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		



Tektronix (1.2)					Monitor model(s): 413A, 414, 414-21			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		

Products for dedicated applications



Terumo					Monitor model(s): 8000, 9000, System 1				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-		
		2	-	-	-	-	-		
		3	-	-	-	-	-		
		4	-	-	-	-	-		
	Contactboard-GO	1	-	-	-	-	-		-
		2	-	-	-	-	-		-
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3142	BDPT-5052	350 cm	1	4046765031138		
		2	-	-	-	-	-		
Plug type: 142		Pins: 5, male		Material: Metal		Additional information: -			



Vytech					Monitor model(s): MostCare, MostCare Up				
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture	
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2084	BDP561PA1	50 cm	1	4046765028671		
		2	75.2085	BDP561PA2	50 cm	1	4046765028688		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
	Contactboard-GO	1	-	-	-	-	-		-
		2	75.1182	BDC902-561	370 cm	1	4046765103347		
		3	-	-	-	-	-		-
		4	-	-	-	-	-		-
Contactboard-GO - FX	1	-	-	-	-	-	-		
DPT-6000	Interface cable	1	-	-	-	-	-		
		2	-	-	-	-	-		
Plug type: 150		Pins: 5, male		Material: Plastic		Additional information: Unused pins not assembled			

Welch Allyn					Monitor model(s): Encore 204EL, Encore 206EL, 244 CS, 246 CS			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2116	BDP156PH1	50 cm	1	4046765069339	
		2	75.2117	BDP156PH2	50 cm	1	4046765069322	
		3	75.2205	BDP156PH3	50 cm	1	4046765108052	
		4	75.2119	BDP156PH4	50 cm	1	4046765069308	
	Contactboard-GO	1	75.1013	BDC901-1056	320 cm	1	4046765024406	
		2	75.1052	BDC902-156	370 cm	1	4046765025151	
		3	75.1074	BDC903-156	370 cm	1	4046765025663	
		4	75.1147	BDC904-156	370 cm	1	4046765078980	
Contactboard-GO - FX	1	75.1040	BDC901FXP	320 cm	1	4046765025021		
DPT-6000	Interface cable	1	75.3103	BDPT-1056	350 cm	1	4046765029524	
		2	-	-	-	-	-	
Plug type: 48		Pins: 12, male		Material: Plastic		Additional information: -		

Witt					Monitor model(s): Calysto, Philips series 4, Physiomonitoring 5, Xper Cardio (Old Module)			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	-	-	-	-	-	
		2	75.2111	BDP582WT2	50 cm	1	4046765028992	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	-	-	-	-	-		
DPT-6000	Interface cable	1	75.3237	BDPT-5082	350 cm	1	4046765109394	
		2	-	-	-	-	-	
Plug type: 166		Pins: 5, male		Material: Plastic		Additional information: Unused pins not assembled		

Products for dedicated applications

Zoll (1.1)					Monitor model(s): M-series, MD-series, X-Series, CCT			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2025	BDP114SL1	50 cm	1	4046765026431	
		2	75.2026	BDP114SL2	50 cm	1	4046765026462	
		3	75.2027	BDP114SL3	50 cm	1	4046765026493	
		4	75.2164	BDP114SL4	50 cm	1	4046765026516	
	Contactboard-GO	1	75.1017	BDC901-114	320 cm	1	4046765026468	
		2	75.1212	BDC902-114	370 cm	1	4046765106089	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1101	BDC901FXSL	320 cm	1	4046765071998		
DPT-6000	Interface cable	1	75.3247	BDPT-1014-SM	350 cm	1	4046765029340	
		2	-	-	-	-	-	
Plug type: 8		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 204		

Zoll (1.2)					Monitor model(s): M-series, MD-series, X-Series, CCT			
	Product type	Signals	REF	Historical no.	Length		GTIN	Picture
Xtrans® (DPT-9000)	Monitor-Connector	1	75.2217	-	50 cm	1	4046765111069	
		2	75.2216	-	50 cm	1	4046765111052	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
	Contactboard-GO	1	-	-	-	-	-	
		2	-	-	-	-	-	
		3	-	-	-	-	-	
		4	-	-	-	-	-	
Contactboard-GO - FX	1	75.1229	BDC901FXPP90-2M	200 cm	1	4046765108632		
DPT-6000	Interface cable	1	75.3243	-	200 cm	1	4046765111083	
		2	-	-	-	-	-	
Plug type: 204		Pins: 6, male		Material: Plastic		Additional information: Equivalent to plug type 8, angled for limited space		

Appendix for monitor models

Plug type	Usage	Monitor Model(s)
11 / 187	Datascope (1.1 + 1.2) .. page 24	2000, 2000A, 2001, 2001A, 2002, 2002A, 2200I, 3000, 870, CS100 IABP, CS300 IABP, P2, P3, Passport, Passport 2, Spectrum, Spectrum 83, Spectrum 90, System 90T, System 95, System 96, System 97, System 98, System 98XT IABP, Trio
100	Datex (1) page 25	CARDIOCAP CC-104, CARDIOCAP CCI-104, CARDIOCAP CMO-104, CARDIOCAP II CG-2CS, CARDIOCAP II CG-2G, CARDIOCAP II CG-2GS, CARDIOCAP II CH-1R, CARDIOCAP II CH-1RS, CARDIOCAP II CH-2, CARDIOCAP II CH-2S, M-ESTP, M-ETP, M-P, AS/3 Light, M-COP, M-COPSV, M-ESTP, M-ESTPR, M-ETP, M-ETPR, M-NESTPR, M-NETPR, M-P, M-PP, M-PT, S/5 Light, M-PRESTN, M-PRETN
213	Dräger (1) page 26	Infinity Delta [w/10-16 Pin Y-Adapter], Infinity Delta XL [w/10-16 Pin Y-Adapter], Infinity Gamma [w/10-7 Pin Y-Adapter], Infinity Gamma X XL [w/10-16 Pin Y-Adapter], Infinity Gamma XL [w/10-7 Pin Y-Adapter], Infinity Hemo2 Pod [w/10 Pin Block], Infinity Hemo2 Pod [w/7 Pin Block & 10-7 Pin Adapter], Infinity Hemo4 Pod [w/10 Pin Block], Infinity Hemo4 Pod [w/7 Pin Block & 10-7 Pin Adapter], Infinity HemoMed Pod [w/10 Pin Block], Infinity HemoMed Pod [w/7 Pin Block & 10-7 Pin Adapter], Infinity Kappa [w/10-16 Pin Y-Adapter], Infinity Kappa XLT [w/10-16 Pin Y-Adapter], Infinity Vista [w/10-7 Pin Single or Y-Adapter], Infinity Vista XL [w/10-16 Pin Y-Adapter]
49	Mennen page 48	2110, CFE [w/8-6 Pin Adapter], Envoy 2BP [w/8-6 Pin Adapter], Envoy 2BP/CO [w/8-6 Pin Adapter], HORIZON 1100, HORIZON 2000, HORIZON 9000 CathLab, HORIZON 9000WS CathLab, HORIZON XL Bedside, MERCURY, MR-1240, MR-1250, MR-1260, MR-1261, MR-1270, MR-1282, MR-1340, MR-1350, MR-1360, MR-1381, MR-9200 CathLab
223	Nihon Kohden (1) page 51	AP-800PA, Life Scope L 2102a, AA-672P SMART INPUT UNIT [W/JP-910P ADAPTER], AA-674P SMART INPUT UNIT [W/JP-910P ADAPTER], AY-661P STD INPUT UNIT [W/JP-910P ADAPTER], AY-663P STD INPUT UNIT [W/JP-910P ADAPTER], BSM-2301(A, J, or K) [W/JP-910P ADAPTER], BSM-2303(A, J, or K) [W/JP-910P ADAPTER], BSM-2304(A, J, or K) [W/JP-910P ADAPTER], BSM-2351(A, J, or K) [W/JP-910P ADAPTER], BSM-2354(A, J, or K) [W/JP-910P ADAPTER], BSM-4101(A, J, or K) [W/JP-910P ADAPTER], BSM-4102(A, J, or K) [W/JP-910P ADAPTER], BSM-4103(A, J, or K) [W/JP-910P ADAPTER], BSM-4104(A, J, or K) [W/JP-910P ADAPTER], BSM-4111(A, J, or K) [W/JP-910P ADAPTER], BSM-4112(A, J, or K) [W/JP-910P ADAPTER], BSM-4113(A, J, or K) [W/JP-910P ADAPTER], BSM-4114(A, J, or K) [W/JP-910P ADAPTER], BSM-5105(A, J, or K) [W/JP-910P ADAPTER], BSM-5106(A, J, or K) [W/JP-910P ADAPTER], BSM-5135(A, J, or K) [W/JP-910P ADAPTER], BSM-5136(A, J, or K) [W/JP-910P ADAPTER], BSM-6501 Life Scope TR [W/JP-910P ADAPTER], BSM-6701 Life Scope TR [W/JP-910P ADAPTER], BSM-9101(A, J, or K) [W/JP-910P ADAPTER], SVM-7603 [W/JP-751P ADAPTER], SVM-7623 [W/JP-751P ADAPTER], Vismo PVM-2703 [W/JP-910P ADAPTER]
198	Nihon Kohden (2) page 52	AA-672P SMART INPUT UNIT, AA-674P SMART INPUT UNIT, AY-661P STD INPUT UNIT, AY-663P STD INPUT UNIT, BSM-2301(A, J, or K), BSM-2303(A, J, or K), BSM-2304(A, J, or K), BSM-2351(A, J, or K), BSM-2354(A, J, or K), BSM-4101(A, J, or K), BSM-4102(A, J, or K), BSM-4103(A, J, or K), BSM-4104(A, J, or K), BSM-4111(A, J, or K), BSM-4112(A, J, or K), BSM-4113(A, J, or K), BSM-4114(A, J, or K), BSM-5105(A, J, or K), BSM-5106(A, J, or K), BSM-5135(A, J, or K), BSM-5136(A, J, or K), BSM-6501 Life Scope TR, BSM-6701 Life Scope TR, BSM-9101(A, J, or K), Vismo PVM-2703
48	Philips (1) page 55	HeartStart MRx, IntelliVue MP2 (M8102A), IntelliVue MP5 (M8105A), MP20 (M8001A), MP20 JUNIOR (M8001A), MP30 (M8002A), MP40 (M8003A), MP50 (M8004A), MP60 (M8005A), MP70 (M8007A), MP90 (M8010A), MX800 (M865240), IntelliVue X2 (M3002A), M1006B, M3000A, M3001A, M3012A, M3014A (C05, C07, C10), M3015A (C06), M3015B (C08), Capnography, SureSigns VM6, SureSigns VM8, Xper Flex Cardio Physiomonitring System
213	Siemens (1) page 69	342R, Cartridge, Hemo2 Pod [w/10 Pin Block], Hemo2 Pod [w/7 Pin Block & 10-7 Pin Adapter], Hemo4 Pod [w/10 Pin Block], Hemo4 Pod [w/7 Pin Block & 10-7 Pin Adapter], HemoMed Pod [w/10 Pin Block], HemoMed Pod [w/7 Pin Block & 10-7 Pin Adapter], HemoMed Pod [w/10 Pin Block], HemoMed Pod [w/7 Pin Block & 10-7 Pin Adapter], SC 5000 [w/10-7 Pin Adapter], SC 6000P [w/10-7 Pin Adapter], SC 6002 [w/10-7 Pin Y-Adapter], SC 6002XL [w/10-7 Pin Y-Adapter], SC 6802XL [w/10-7 Pin Y-Adapter], SC 7000 [w/10-16 Pin Y-Adapter], SC 8000 [w/10-16 Pin Y-Adapter], SC 9000 [w/10-16 Pin Y-Adapter], Sirecust 620, Sirecust 630, Sirecust 730, Sirecust 732, Sirecust 402-P, Sirecust 403-2P, Sirecust 403-P
8 / 204	Spacelabs (1.1 + 1.2) .. page 72	Alpha 7, 3700, 413A, 414, 414-21, 511, 512, 512D, 514, 514T, 521, 522, 524, 702S, 703S, 90402, 90403, 90425, 90428, 90432, 90470, 90470-04, 90470-05, 90470-06, 90470-10 (w/7-6 Pin Y-Adapter), 90496-B (w/10-6 Pin Y-Adapter), 90496-C (w/10-6 Pin Y-Adapter), 90602A, 90603, 90622A, 90623A, 91496-B Ultraview (w/10-6 Pin Y-Adapter), 91496-C Ultraview (w/10-6 Pin Y-Adapter)

CODAN Worldwide



CODAN is known internationally as a manufacturer and supplier of medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

CODAN Companies

CODAN Medizinische Geräte GmbH · Deutschland
CODAN pvb Critical Care GmbH · Deutschland
CODAN pvb Medical GmbH · Deutschland
CODAN 11, S.A. · Portugal
CODAN US Corporation · California · USA
CODAN Inc. · California · USA
CODAN NORGE AS · Norge
CODAN TRIPLUS AB · Sverige
CODAN Limited · Great Britain
CODAN FRANCE Sarl · France
CODAN Medical AG · Schweiz
CODAN ARGUS AG · Schweiz
CODAN BV · Nederland
CODAN s.r.l. · Italia
CODAN Medical GmbH · Österreich
CODAN Medical ApS · Danmark
CODAN DEHA ApS · Danmark
CODAN MEDITECH s.r.o. · Česká republika

CODAN Product range

- Infusion sets
- Transfusion sets
- Extension lines and manifold connectors
- Infusion and transfusion accessories
- Infusion filters and filter systems
- Neonatology/Paediatric products
- Withdrawal, preparation and administration systems
- CODAN CYTO®
- Chemoprotect® products
- Single use syringes
- Invasive blood pressure monitoring systems
- Infusion pumps
- Other CODAN Products

Compliance of the established quality management systems with the provisions of EN ISO 13485 and the Council Directive 93/42/EEC has been certified by the relevant, competent notified bodies:

TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH
85661 Forstinning, Germany

CODAN US Corporation · Santa Ana, CA 92704, USA

CODAN ARGUS AG · 6340 Baar, Switzerland

Presafe Denmark A/S

CODAN Medical ApS · 4970 Rødby, Denmark

Manufacturer

CODAN pvb Critical Care GmbH
Römerstraße 18 · 85661 Forstinning, Germany
Tel. +49 (0) 81 21 · 98 02 0
codan@codanpvbcc.de · www.codancompanies.com



The decisive connection